

DENTAL

Catalogue



RAIZ & SAJJAD
SURGICAL PVT. LTD.

Manufactures and Exports of Fine Quality Surgical & Dental Instruments





Factory view

Raiz & Sajjad Surgical (PVT) Ltd

We feel immense pleasure to present our latest catalogues in which we have featured a most comprehensive range of all quality Dental Instruments required for day-to-day use for SURGEONS, OPERATING ROOMS, NURSES CENTRAL SUPPLY, CLINICS, FLOOR USE, EMERGENCY ROOMS, PROFESSIONAL AND OTHER MEMBERS OF THE DENTAL CARE TEAM. Without their loyalty, support and steady flow of innovative ideas, we would not be marking this momentous occasion.

We have **ISO 9001:2000, ISO: 13485:2003, cGMP Approved and CE Marked**. We wear this badge with pride. We have implemented not to our product, but rather to a Company's documented quality policies, processes and procedures.

All instruments are manufactured according to the recognized standard patterns currently in use in Hospitals throughout the world and are being produced under the strict control of experts to avoid any flaw and to refrain any complaint at any stage. If there is any way we can be of service to you, please contact us without any hesitation.

We express our thanks to our entire customers for the confidence they have placed in us in the past. We further assure them, that we will Continue our best efforts to gain their confidence in the future.

Quality is never an accident; It's always the results of an intelligent effort.

Briefly Corporate History:

➤ **Introduction & Start up year:**

To introduce the firm briefly, M/s. **Raiz & Sajjad Surgical (PVT) Ltd.** was established in 1969 by who was well known as a specialized skilled worker and maker of Dental Instruments.

➤ **Raiz & Sajjad Surgical (PVT) Ltd. Story:**

He used to supply his hand made products within the entire territory of sub-continent prior when the India and the Pakistan were combined.

➤ **Dynamic Expansion:**

Then the son of the founder, entered into the business under the supervision of his great father IN 1991. He sorted out the markets through concrete and restless efforts and he succeeded. He also started export to Europe and gradually developed and expanded his business around the Globe. This company is totally involved in Dental field, covering operation room, High Grade Quality to Floor Grade etc.

➤ **Business Review:**

Since then the Raiz & Sajjad Surgical (PVT) Ltd. are exporting the major lots of quality products to USA, UK, Germany, UAE and South America etc., and their export Performance is raising up yearly and now looking ahead. By the grace of God Almighty, This firm is recognized as a leading Manufacturing / Exporting Unit of Sialkot, Pakistan.

Raiz & Sajjad Surgical (PVT) Ltd. is a unique leader in producing the best quality standard instruments and is well familiar with quality controlling process by using most modern technology of International quality control systems/programme. Our all instruments are carefully passivated, electrolytically polished and ultrasonically cleaned with best chemicals available. Therefore, we **GUARANTEE** of our instruments against any eventual manufacturing defects and material, to either Replaced or credit.

Customer Satisfaction is our priority, business, recognition & basic aim.



Company Profile & Introductory Letter

Introductory Letter

Subject: Supply of Dental Instruments

Dear Sir/Madam,

We introduce ourselves as pioneer Manufacturers cum Exporters of Dental Instruments in super quality here in Pakistan beg to our most excellent service to cater against your valuable requirements.

Our organization, backed by the history of over four decades fully equipped with the most modern machinery, hardening, tempering, plating plants having adequate testing facilities by the highly skilled and most experienced craftsmen / technicians, is producing unrevealed quality carrying an unconditional GUARANTEE for:-

- Use of best steel of proper specification in the production, flawless workmanship and superb finish.
- Electrolytically passivated and ultrasonically cleaned with the best-imported chemicals to provide maximum resistance to corrosion.
- Exact match of samples approved by the buyers, as well as, uniform standard of fabrication and finish for any quantities.
- Immediate settlement of any item found to be of substandard quality replaced to the customers.

CATALOGUES:

Illustrated Dental Catalogue is being provided for your future reference.

DELIVERY:

All orders are being supplied within the stipulated period of time. However, the orders of special and non-existing items can also be executed upon receipt of sketch, drawings or specifications from the customer.

SAMPLES:

In order to just take an idea about our quality, please ask for necessary samples of your interest against which we shall be very glad to submit without any further delay.

PRICE:

Now-a-day, due to great competition, this business is being conducted on quotation basis, therefore, you are requested to kindly provide us a list of your desire items against which we shall appreciate to quote / offer our most competitive prices for favour of your kind approval.

Looking forward to your early kind encouraging reply.

With Cordial Greetings,

Our Export Market

We are manufacturing our quality products and exporting more than 57 countries of the World including major markets such as USA, UK, EU, South America, Middle East, Gulf and Asian countries.

Our Strategy

- ❖ Fair dealing, Straight Forwardness, Honesty, Sincerity, Truth
- ❖ The Material used in the production against order is exactly as per approved samples.
- ❖ Normally, in Cutting instrument like Scissors, Knives etc; we use (420 Type) Quality stainless Steel.
- ❖ For Non-Cutting instruments like Forceps, Tweezers etc; we use (410 Type) Quality Stainless Steel.
- ❖ Moreover, we can use also the material as per our customer requirements.
- ❖ We always cordially welcome our valued customer's opinions, suggestion and recommendation.
- ❖ We give replacement against defective goods, if the malfunction occurred at our end.

The best you advise us, the better we will serve you

Our Goals & Objectives:

All of our goals and objectives are reviewed on a periodic basis to ensure that they are used as a means of continuous improvement for the organization.

This is our family Business, so we know it better

QUALITY POLICY

Raiz & Sajjad Surgical (PVT) Ltd. is committed to Produce Quality Goods as per Customer Requirements International Standard through Continuous Improvement.


Objectives;

- To Improve customers and users expectations by reducing the complaints to below 4%.
- To meet legal and regulatory requirements at international level.
- To keep company rejection below 5%.
- To Improve supplier's performance up-to 90 100%.
- To Improve workers performance up-to 90 100%.

The quality policy is the main guiding principles that provide the focus of Quality to the whole organization and improve the effectiveness of the quality management system.

Quality Management System

Process Model



Quality Begins With People Who Care

All Instruments are guaranteed to be free from defects in workmanship and material. Any Instrument, which proves defective in workmanship or material, will either be issued creditor replaced at our discretion, without charge.

All instruments shown in this catalogue are made of stainless steel, unless otherwise indicated.
Illustration signify general design only


RSS is registered trademark of
Raiz & Sajjad Surgical (PVT) Ltd.
Sialkot – Pakistan



(cGMP)
Approved



Process Flow Chart





Instruments Care:

To Our Customers:

Concerning your question on why staining, pitting or rusting can develop on stainless steel instruments; there are some basic points for the care and handling of Dental instruments which will extend the life of the instrument, eliminate staining or corrosion and lower the cost in use.

These are as follow:

1. Instruments should be rinsed as soon as possible after use to remove blood, serum and saline. Do not use a sponge basin as any residue from normal saline solution can cause corrosion to develop over a period of time.
2. Instruments should be thoroughly dried after cleaning and before wrapping or storing, as moisture left in the box lock or serrations can cause rust or corrosion to develop.
3. Do not leave instruments in a cold sterilization soak solution for extended periods as these solutions tend to be corrosive (after cold sterilization, rinse the instruments in distilled or demineralized water and dry thoroughly).
4. If using cold sterilization technique, change the solution according to directions, as prolonged use will cause the solution become corrosive.
5. Do not use caustic cleaners or any other cleaner except those specified for cleaning Dental instruments
6. Do not clean instruments with steel wool or scouring powder, or use high PH detergents. If cleaning by hand, use specific instruments cleaner or a neutral detergent soap and a soft brush (rinse with demineralized or distilled water).
7. Check your autoclave. Iron, sodium, calcium, magnesium or copper in your water can cause spotting, staining or corrosion to occur. In hard water areas, clean the lime deposits from steam lines periodically or use distilled or demineralized water for your autoclave. Install a steam filter if possible. It will remove most of dirt, rust and pipe scale.
8. Make sure your muslin wraps are rinsed sufficiently (6 to 8 times) by the laundry to remove all chlorine bleach as any residue can cause staining or corrosion during autoclaving.
9. Do not leave instruments in the following solutions for extended periods as corrosion can result:-
(a) Aluminum, Barium, Calcium, Ferrous or Stannous Chloride solutions.
(b) Phenol, Lysol or Iodine. (c) Dakin's Solution or Zephiran. (d) Any acid, Mercury or Potassium solution.
10. Do not sterilize stainless steel instruments with chrome plated instruments if possible, as any break in the surface of the chrome plate will allow an electrolytic action to develop between the dissimilar metals, causing pitting and/or a rusty appearance to develop.
11. Preheat instruments in autoclave before allowing steam to enter, if possible, to prevent an excessive amount of condensation forming on the surface of the instruments.
12. Don't open the autoclave too quickly after the sterilization cycle as excessive condensation will develop on the surface of the instruments (crack the clave to allow drying before opening fully). If you follow the above rules and use instrument cleaners and lubricant milk such as our Ami-Clean or Ami-Tec products, your instruments should remain rust, spot and stain free. I hope this will answer any question you might have on this topic.

When proper care is taken these stainless steel instruments can last for many years.

Sterilization Method:

Sterilization in Autoclave:

1. After having carefully washed dried and oiled instruments place instruments that can be opened in the open position; put these instruments in sterilization packets and place in appropriate trays.
2. Be careful not to over load autoclave. When an autoclave is too full it's possible that proper sterilization has not taken place.
3. After the sterilization cycle has taken place, follow manufacture's suggestions for drying cycle.

Cold Sterilization:

When using a solution for sterilizing instruments careful attention must made that the solution kills specific bacteria and virus, many solutions have no effect on bacteria or virus when in the presence of blood; also many solution are very corrosive on the instruments, most solutions recommend a soak time of 10 hrs. Because of the long soak time the chemical action can damage the instrument. Reduce when possible soak time; where indicated use a corrosion inhibitor.

Important Information:

To our Customer:

- 1 To provide "Excellent Service and Economically Price" is the motto of our business.
- 2 All instruments are made of Stainless Steel otherwise specified.
- 3 All types of "Mirror-Finish" "Satin/Dull-Finish" "Sand Blast-Finish" instruments are available in stock but choice of polishes will be supplied on demand without any extra charges.
- 4 All Disposable/Floor Grade Instruments are stamped with words "STAINLESS" and "PAKISTAN". Buyer's own brand may also be engraved on demand. "Die-Sunk" procedure generally used.
- 5 The etching stamps are only used for "High Quality" Instruments and keep in stock as per our customer's demand. These are also available in "Mirror-Finish" "Satin/Dull-Finish" "Sand Blast-Finish" and are processed through boil test and passivated, which is a warranty against oxidation.
- 6 Each and every piece is individually packed in airtight polythene bag in order to avoid by contamination but can be packed as per customer's demand.
- 7 All general items are stored in big quantities in finishes three "Mirror-Finish" "Satin/Dull-Finish" "Sand Blast-Finish" in order to fulfill the urgent requirements of our valuable customer's immediately and price would be same for both finishes.
- 8 We offer Freight prepaid and collect prices as per Customer's requirements.
- 9 Compromise on running inquiry is our one of the greatest characteristics of our dealing.
- 10 While ordering, specify quantity, catalogue number and description completely.
State if order is to be dispatched by UPS, FedEx, DHL or other means of transportation.



A Traditional of Excellence & Production Capabilities

- 11 Easy replacement in case of any damage.
- 12 Advise carrier promptly of any damage or shortage.
- 13 Merchandise will not be accepted for exchange or credit unless or until prior permission for its return has been granted.
- 14 Samples are provided at free of cost.
- 15 Special items can also be supplied if sketch, drawing or samples are provided to us or refer any reference from any reputed catalogue for our guidance.

Specification of material used production:

For	Cutting Instruments	Stainless Steel AISI/ASTM (Rock Well Hardness	(50° to 55°)	Hot drop forged
For	Extracting Forceps & Retractor	Stainless Steel AISI/ASTM (Rock Well Hardness	410/420 (45° to 50°)	Hot drop forged
For	Scalar, Gauges and SCALPEL / HANDLES	Stainless Steel AISI/ASTM (Rock Well Hardness	(50° to 55°)	Hot drop forged
For	Impression Trays	Stainless Steel Sheet (Non-Magnetic)		(18% CHROMIUM & 8% NICKEL)
For	Root Elevators and Pliers	Stainless Steel Sheet (Magnetic & Non-Magnetic)	300/410 (45° to 50°)	Hot drop forged

Raiz & Sajjad Surgical (PVT) Ltd. is an employee owned company of Craftsmen and has got a good reputation for quality and reliability. Our significant manufacturing steps are being described here-under briefly:

- 1. MATERIAL USED:** Always a best quality stainless steel material is purchased from steel manufacturers. Two kinds of series i.e. 300, 410 and 420 of stainless steel is used for production of Dental Instruments.
- 2. MATERIAL TESTING:** The purchased raw material is thoroughly tested under PQA. authority. A rigid quality control policy is being adopted before releasing to production. This test is taken to check the standard of steel as per desired specifications and is considered the most important factor of our business.
- 3. FORGING & HAMMERING:** After testing, the stainless steel strips/bars are divided into different sizes as per our requirements and our highly experienced craftsmen are specialized in forging (Hot/Cold Drop) to perform all standard press operations in forming and deep drawing of all type of instruments. We anneal the forging and prepare for the milling process.
- 4. MACHINE & TOOL WORK:** Our highly skilled technicians have a great knowledge about machine milling process. We form accurately male/female portion of box-joint forceps, as well as, serrations/ratches in uniformity. We are equipped to mill, drill, bore turn and grind all type of metals. We have in house facility of custom tool and die production, maintenance and repair.
- 5. FILLING:** We are specialized in handmade fabrication. First, we assemble male/female portion after setting/adjusting of Extracting Forceps and Scissors, tips, serrations, ratches in uniformity and file precisely to make it semi-finish under strictly quality control specifications. We are now moving from basic hand manufacturing techniques to mechanization to production run.
- 6. HEAT TREATMENT:** All instruments are being tempered maintaining "Rock-Well Hardness" for Extracting Forceps & Retractor 45° to 50° degree and Scalar, Gauges and Scalpel 50° to 55° degree. A test is taken to check the hardness after tempering of each lot to maintain the uniformity and, before proceeding further, we clean all instruments by various chemicals. In addition, we will buy specialized Q/C equipment in near future to meet special customer Q/C. Specifications.
- 7. RIVETING & WELDING:** We perform a wide variety of assembly operations from prototype quantities to production run. We rivet precisely each instrument such as, Extracting Forceps & Retractor, etc. But, for some kind of instruments, we use spot welding.
- 8. GRINDING & POLISHING:** For precise grinding and excellent finish, we specialize in applying polishing wheels to belts to improve the overall finish in high uniformity utilizing best of polishing luster. We use electro-plating plants to polish serrations/ratches and joints of stainless steel instruments. For satin/matt finish instruments, we use glass heads to achieve uniformity. All instruments are being polished under rigid quality control specifications.
- 9. PASSIVATION:** This method actually provides the instruments with its corrosion resistance properties. Passivation, which follows the final polishing step, is a chemical bath, which creates an oxidized layer on the surface of instruments. We are now equipped with most modern process.



CONTENTS

RSD1

EXTRACTING FORCEPS

Page # 01 - 74

RSD2

ROOT ELEVATORS

Page # 75 - 123

RSD3

**DIAGNOSTIC, PROBES, SCALERS, FILLING
& LABORATORY INSTRUMENTS**

Page # 124 - 242

RSD4

**MAXILLOFACIAL SURGERY
INSTRUMENTS**

Page # 243 - 328

RSD5

**ORTHODONTIC INSTRUMENTS
& PLIERS**

Page # 329 - 354

RSD6

DENTAL IMPLANT INSTRUMENTS

Page # 355 - 372

RSD7

IMPRESSION TRAY


Page # 373 - 388

RSD6

HOLLOWWARE

Page # 389 - 397

Handle Selection For Extracting Forceps



English Pattern




English Pattern



English Pattern,
Anatomical Handle



American Pattern



Deep Grip Forceps
With Hole Pattern



Hole Pattern



A



B



C



D



E



F



Extraction Forceps Function

Dental forceps are pliers-like instruments that are used to grasp teeth.

A dentist will usually have a number of different ones on hand, each having a design that's tailored to:

The specific shape of the tooth it's intended to remove (such as incisors vs. molars having two roots vs. molars having three roots).

The location of the tooth in the mouth (such as front vs. back, or left lower molar vs. right upper molar.)

How they're used.

A dentist will grasp a tooth with their forceps and then firmly and deliberately rock it back and forth as much as it will.

Because the bone that surrounds the tooth is compressible, the socket will gradually start to expand.

As it does, the range of the dentist's side-to-side motions will increase.

In addition to this rocking motion, the dentist will also rotate the tooth back and forth.

This twisting action helps to rip and tear the tooth away from the ligament that binds it in place.

At some point, the socket will be enlarged enough, and the ligament torn away enough, that the tooth can be easily removed.

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)




RSD1-001
1
upper laterals
and canines



RSD1-002
2
upper laterals
and canines



RSD1-003
3
upper laterals
and canines



RSD1-004
4
lower laterals
and canines



RSD1-005
7
upper premolars



RSD1-006
8
lower premolars



RSD1-007
13
lower premolars

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)


Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)




RSD1-008
13S
lower premolars,
for children



RSD1-009
17
upper molars,
right



RSD1-010
18
upper molars,
left



RSD1-011
18A
upper molars



RSD1-012
18B
upper molars



RSD1-013
19
upper third molars

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-014
20
lower third molar



RSD1-015
21
lower molar



RSD1-016
22
lower molar



RSD1-017
22S
lower molar,
for children



RSD1-018
22G
lower molar



RSD1-019
23
lower molar, right




RSD1-020
24
lower molar, right

EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-028
33
lower roots



RSD1-029
33A
lower roots



RSD1-030
33B
lower roots




RSD1-031
33L
lower roots



RSD1-032
33S
lower roots,
for children



RSD1-033
Fig. 33C
Lower roots



RSD1-034
Fig. 34
Upper incisors
and premolars
deep-gripping

RSD1-035
Fig. 35
Upper incisors
and premolars
deep-gripping



RSD1-036
Fig. 35N
Upper incisors
and premolars
deep-gripping



RSD1-037
Fig. 33M
Lower roots
with long beak



RSD1-038
Fig. 36
Lower premolars
deep-gripping



RSD1-039
Fig. 36N
Lower premolars
deep-gripping



RSD1-040
37
upper incisors,
for children

Quality Symbol

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-041
38
lower incisors and
canines, for children



RSD1-042
39
upper molars
for children



RSD1-043
39 R
upper molars right
for children



RSD1-046
39 L
upper molars left
for children



RSD1-044
39 A
upper molars
for children



RSD1-045
40
lower molars
for children

EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-047
Fig. 41
upper roots



RSD1-048
Fig. 44
upper roots



RSD1-049
Fig. 45
lower roots



RSD1-050
Fig. 46
Lower root



RSD1-051
Fig. 46N
Lower premolars
deep-gripping



RSD1-055
50
upper roots



RSD1-052
Fig. 48
lower roots



RSD1-053
Fig. 49
upper roots



RSD1-054
Fig. 49L
Very thin lower roots

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)


استخراج الملقط (نمط اللغة الإنجليزية)




RSD1-056
51
upper roots



RSD1-057
51 A
upper roots



RSD1-058
51 L
upper roots




RSD1-059
Fig. 51 C
Upper roots



RSD1-060
51 S
upper roots,
for children



RSD1-061
52
upper roots



RSD1-062
52 A
upper teeth
and resection

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)


استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-063
54
for separating
upper molars




RSD1-064
55
for separating
upper molars



RSD1-065
56
for separating
lower molars



RSD1-066
Fig. 59
Lower root



RSD1-067
65 R
upper molars,
right



RSD1-068
65 L
upper molars,
left



RSD1-069
67
upper third molars



RSD1-070
67 A
upper third molars

EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-071
67 1/2 R
 upper third molars,
 right



RSD1-072
67 1/2 L
 upper third molars,
 left



RSD1-073
 Fig. 67N
 Upper third molars



RSD1-074
68
 lower roots



RSD1-075
68 A
 lower roots



RSD1-076
73
 lower molars



RSD1-077
73 R
 lower molars,
 right



RSD1-078
73 L
 lower molars,
 left

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-079
73 A
lower molars



RSD1-080
73 S
lower molars



RSD1-081
73 B
lower molars
for children



RSD1-082
74
lower roots



RSD1-083
74 D
lower roots



RSD1-084
73 K
lower molars
for children

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-085
74M
lower roots



RSD1-086
74N
lower roots



RSD1-087
74XN
lower roots and
crowded incisors



RSD1-088
75
lower premolars



RSD1-089
76
upper roots



RSD1-090
76N
upper roots

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)


Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Dabiers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-098
86A
lower molars



RSD1-099
86B
lower molars




RSD1-100
86C
lower molars



RSD1-101
87
lower molars



RSD1-102
89
upper molars,
right



RSD1-103
90
upper molars,
left


EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)




RSD1-104
92
upper second molars,
right



RSD1-105
93
upper second molars,
left



RSD1-106
94
upper molars,
right



RSD1-107
91
for separating
lower molars



RSD1-108
99
lower molars



RSD1-109
Fig. 97
Very thin
upper roots



RSD1-110
95
upper molars,
left

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)




RSD1-111
99 1/4
lower molars



RSD1-112
99 1/2
lower molars



RSD1-113
Fig. 233
Lower root
(parallel beaks)



RSD1-114
Fig. 235
Upper root
(parallel beaks)



RSD1-115
100
upper roots



RSD1-116
101
upper premolars
and roots

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-117
104
upper roots



RSD1-118
102
upper molars,
right



RSD1-119
103
upper molars,
left



RSD1-120
105
lower incisors



RSD1-121
106
lower incisors



RSD1-122
107
upper canines

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)


Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-123
112
upper roots



RSD1-124
113
upper roots



RSD1-126
136
upper premolars
and roots



RSD1-127
137
lower incisors
and roots



RSD1-125
114
lower roots



RSD1-128
143
lower roots

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)


Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)


استخراج الملقط (نمط اللغة الإنجليزية)



**RSD1-134
STIKSA**
lower molars
and wisdoms



**RSD1-135
STIKSA**
lower molars and
wisdoms roots



**RSD1-136
HERTZ**
lower molars
and wisdoms roots



**RSD1-137
HERTZ**
lower incisors

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-138

Fig. 158

Children's Upper Teeth Molars Universal



RSD1-139

Fig. 159


Children's upper Molars Universal



RSD1-140

Fig. 160

Children's Upper Lower Universal (Hawk's Bill)



RSD1-141

Fig. 161

Children's Lower Centrals & Roots



RSD1-142

Fig. 162


Children's Upper Centrals



RSD1-143

Fig. 163

Children's Upper Centrals



RSD1-144

Fig. 37

Children's Upper Incisors & Canines

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)


Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-145
Fig. 123

Children's Lower Teeth Universal




RSD1-146
Fig. 138

Children's Upper Universal & Roots




RSD1-147
Fig. 157

Children's Upper Morals Universal




RSD1-148
Fig. 39

Children's Upper Morals Universal



RSD1-149
Fig. 39R

Children's Upper Morals Right



RSD1-150
Fig. 39L

Children's Upper Morals Left

WISDOM EXTRACTING FORCEPS

Щипцы экстракционные для зубов мудрости (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні для зубів мудрості (англійська форма)

Forceps de Extracción (Forme Inglesa)

الحكمة استخراج الملقط



RSD1-151
Routurier
Fig. 22 1/2 L
Lower Molars and
Wisdoms, Left



RSD1-152
Routurier
Fig. 22 1/2 R
Lower Molars and
Wisdoms, Right



RSD1-153
Routurier
Fig. 33 1/2 L
Lower Molar and
Wisdom Roots, Left



RSD1-154
Routurier
Fig. 33 1/2 R
Lower Molar and
Wisdom Roots, Right



RSD1-155
Fig. 186L
Lower Molars



RSD1-156
Fig. 186R
Lower Molars

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)


Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)


استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-157
Busch
lower roots



RSD1-158
Busch
lower roots



RSD1-159
Witzel
lower roots



RSD1-160
Fig. 146
Upper and lower roots
deep-gripping



RSD1-161
Fig. 333S
Very thin lower
roots



RSD1-162
Fig. 151
upper roots



RSD1-163
Fig. 351S
Very thin
upper roots

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)


Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)


استخراج الملقط (نمط اللغة الإنجليزية)




RSD1-164
Fig. 121
Upper Third
Molars



RSD1-165
Fig. 136
Upper Premolars
and Roots



RSD1-166
Fig. 137
Lower Incisors
and Roots



RSD1-167
Fig. 167
Universal Patter
For Lower Teeth
Molars & Premolars




RSD1-168
Fig. 168
Upper Canines



RSD1-169
Fig. 169
Lower Canines &
Premolars



RSD1-170
Fig. 233
Lower Root
(Parallel Beaks)



RSD1-171
Fig. 235
Upper Root
(Parallel Beaks)

EXTRACTING FORCEPS (For Children)

Щипцы экстракционные детские
Daviers - forme anglaise
Pinzeper per estrazione per bambini

Щипці екстракційні дитячі
Fórceps de extracción - forma inglesa
ستخراج الملقط (للأطفال)

Pediatric Tooth Extracting Forceps



RSD1-172
English Pattern
Klein Fig. 4
Pedodontic Upper root

RSD1-175
English Pattern
Klein Fig. 1
Pedodontic, Upper Incisors



RSD1-176
English Pattern Klein Fig. 5
Pedodontic Lower Incisors.



RSD1-173
English Pattern
Klein Fig. 2
Pedodontic Upper Premolars

Baby Tooth Extraction Forceps




RSD1-179
Set of 7 Pieces



RSD1-177
English Pattern
Klein Fig. 6
Pedodontic Lower Molars



RSD1-174
English Pattern
Klein Fig. 3
Pedodontic Upper Molars Universal



RSD1-179-A
Tooth Forceps Children Set For 7Pcs




RSD1-178
English Pattern
Klein Fig. 7
Pedodontic Lower Molars,

	Fig.
01 Tooth Forceps	1
02 Tooth Forceps	2
03 Tooth Forceps	3
04 Tooth Forceps	4
05 Tooth Forceps	5
06 Tooth Forceps	6
07 Tooth Forceps	7

EXTRACTING FORCEPS (For Children)

Щипцы экстракционные детские
Daviers - forme anglaise
Pinzette per estrazione per bambini

Щипці екстракційні дитячі
Fórceps de extracción - forma inglesa
ستخراج الملقط (للأطفال)



RSD1-187
Baby Extracting Forceps
Set of 07 Pcs
Fig 1
Fig 2
Fig 3
Fig 4
Fig 5
Fig 6
Fig 7


EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-188
Fig. 151 SK
 Lower Incisors,
 Premolars, Roots



RSD1-189
Fig. 222 SK
 Lower Third Molars



RSD1-190
Fig. 560
 Upper Incisors

Deep Gripping Root Forceps



RSD1-191
Fig. 146
 Upper Roots



RSD1-192
Fig. 149
 Lower Roots



RSD1-193
Fig. 153
 Upper Roots

ROOTS SPLINTER FORCEPS

Щипцы для удаления корневых осколков
Pincers Fragments
Wurzelsplitterzangen


Щипці для видалення корневих осколків
Pinzas Para Fragmentos
چنور منشقة الملقط



RSD1-194
Splinter Forceps
Straight
18cm



RSD1-195
Splinter Forceps
18cm
45°



RSD1-197
Stieglitz
Root Fragments



RSD1-198
Stieglitz
Root Fragments



RSD1-199
Stieglitz Splinter Forceps
14cm 90°



RSD1-200
Stieglitz Splinter Forceps
14cm Straight

EXTRACTING FORCEPS (For Children)

Щипцы экстракционные детские

Daviers - forme anglaise

Pinzeper per estrazione per bambini


Щипці екстракційні дитячі

Fórceps de extracción - forma inglesa


ستخراج الملقط (للأطفال)




RSD1-201
Koehler
Luxation Forceps
Right & Left Molars



RSD1-202
Upper Canines
& Premolars



RSD1-203
Lower Molars
Either Side



RSD1-204
Lower Third Molars



RSD1-205
Albrecht
16cm



RSD1-206
Albrecht
16cm

deciduous teeth

EXTRACTING FORCEPS (For Children)

Щипцы экстракционные детские

Daviers - forme anglaise

Pinzeper per estrazione per bambini

Щипці екстракційні дитячі

Fórceps de extracción - forma inglesa

ستخراج الملقط (للأطفال)



RSD1-207
Fig.10 SK
Upper Molars Sup, 10SK



RSD1-208
Fig.17 SK
Universal
Lower Molars



RSD1-209
Fig.23 SK
Universal
Upper Molars



RSD1-210
Fig. 69 SK
Tomas
Upper and Lower Roots



RSD1-211
Fig. 99 SK
Upper Incisors,
Canines, Premolars



RSD1-212
Fig. 150 SK
Upper Incisors,
Premolars, Roots



RSD1-213
Fig. 151 SK
Lower Incisors,
Premolars, Roots



RSD1-214
Fig. 222 SK
Lower Third Molars


EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)

Deep Gripping Root Forceps


RSD1-215
Upper Centrals
Fig. 34A



RSD1-216
Upper Premolar
Fig. 35A



RSD1-217
Upper Roots
Fig. 331



RSD1-218
Lower Premolar
Fig. 336



RSD1-219
Lower Premolar
Fig. 36A



EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)

Deep Gripping Root Forceps



RSD1-220

Upper Third Molars


Fig. 67N



RSD1-221

Lower Third Molars

Fig. 79N



RSD1-222

Upper Premolars

Fig. 7N



RSD1-223

Lower Premolars

Fig. 36N

EXTRACTING FORCEPS (English Pattern Hole Handle)


Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)

(Deep Gripe Extracting Forceps)



- No Fracture
- No Tissue Damage
- Extraction With Easy



Thin Tip



RSD1-224
Fig. 34
Upper Anterior

Wide Tip



RSD1-225
Fig. 34A
Upper Anterior

Thin Tip



RSD1-226
Fig. 35
Upper Universal

Wide Tip



RSD1-227
Fig. 35A
Upper Universal

Thin Tip




RSD1-228
Fig. 321
Lower Universal

A-TRAUMATIC EXTRACTION FORCEPS
Improved tapered beaks for complete root contact.
Cross serrations for a never slip grasp.
Easily enter the sub gingival.

EXTRACTING FORCEPS (English Pattern Hole Handle)

Щипцы экстракционные (английская форма)
Daviers (Forme Anglaise)
Zahnzangen (Englische Form)


Щипці екстракційні (англійська форма)
Forceps de Extracción (Forme Inglesa)
استخراج الملقط (نمط اللغة الإنجليزية)



Deep Root Pick Forceps



RSD1-229
Fig. 321A
Lower Root Molars



RSD1-230
Fig. 36
Lower Premolars
Deep-Gripping



Wide Tip



RSD1-231
Fig. 36A
Lower Premolars
Deep-Gripping



RSD1-232
Fig. 67
Upper Thrid Molars




RSD1-233
Fig. 67A
Upper Thrid Molars

EXTRACTING FORCEPS (English Pattern Hole Handle)

Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)

DEEP ROOT PICK FORCEPS



RSD1-234
Fig. 79
 Lower Third Molars



RSD1-235
Fig. 79A
 Apex Forceps Lower Universal
 Deep Root Pick



RSD1-236
Fig. 79 T
 Lower Third Molars



RSD1-237
Fig. 44
 Apex Forceps Upper Root
 Deep Root Pick




RSD1-238
Fig. 45
 Apex Forceps Lower Root
 Deep Root Pick

EXTRACTING FORCEPS (English Pattern Hole Handle)

Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)


Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



DEEP ROOT PICK FORCEPS



RSD1-239
 Fig. 51
 Apex Forceps
 Upper Root
 Deep Root Pick



RSD1-240
 Fig. 46XL
 Apex Forceps Lower Root
 Deep Root Pick



RSD1-241
 Upper Root
 Fig. 44



RSD1-242
 Lower Root
 Fig. 45



RSD1-243
 Upper Root
 Fig. 51

EXTRACTING FORCEPS (American Pattern)

 Щипцы экстракционные (американская форма)
 Davigers - forme américaine
 Fórceps de extracción - forma americana

 Щипці екстракційні (американська форма)
 Pinzere per estrazione - forma americana
 استخراج المنقط (النمط الأمريكي)

Carbide Tipped Root Tip Forceps



Anatomical Handle



RSD1-244
T.C.
Carbide Serrated Beak
Lower Root Tips & Molars



RSD1-245
Root Tip-Lower
Carbide Serrated Beak
Lower Root Tips & Molars



RSD1-246
Root Tip-Upper
Carbide Serrated Beak
Lower Root Tips & Molars

German Stainless Steel Extraction Forceps



RSD1-247
150AS
Serrated Beak
Upper Anteriors,
Split Beaks



RSD1-248
151AS
Serrated Beak
Lower Anteriors,
Split Beaks



RSD1-249
23
1st & 2nd Lower Molars
Universal, Cow Horn

Wrap Around Forcep



RSD1-250
Wrap around extraction forcep
right angle "Wrap Around" Design
Provides Superior Access
Serrated Beak

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)
 Daviers - forme américaine
 Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)
 Pinzere per estrazione - forma americana
 استخراج المنقط (النمط الأمريكي)

Carbide Tip Extraction Forceps - Better Grip & Tactile Feel - Less Trauma



RSD1-252
 151 T.C.
 Carbide Serrated Beak
 Lower Incisors, Cuspids,
 Bicuspsids & Roots, Universal

RSD1-253
 1STD T.C.
 Carbide Serrated Beak
 Upper Incisors, Cuspids,

RSD1-254
 17 T.C.
 Carbide Serrated Beak
 1st & 2nd Lower Molars,
 Universal



RSD1-255
 53L T.C.
 Carbide Serrated Beak
 1st & 2nd Upper Molars,
 Left

RSD1-256
 53R T.C.
 Carbide Serrated Beak
 1st & 2nd Upper Molars,
 Right



RSD1-257
 65 T.C.
 Carbide Serrated Beak
 Upper Incisors, Roots &
 Fragments



RSD1-251
 150 T.C.
 Carbide Serrated Beak
 Upper Incisors, Cuspids,
 Bicuspsids & Roots,
 Universal

RSD1-258
 69 T.C.
 Carbide Serrated Beak
 Upper & Lower Fragments
 And Roots

RSD1-259
 150S T.C. Pedro forcep
 Carbide Serrated Beak
 Pedro, Upper Primary Teeth
 and Roots, Universal

RSD1-260
 151S T.C. Pedro forcep
 Carbide Serrated Beak
 Pedro, Lower Primary Teeth
 and Roots, Universal

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-261
45° Angled Root Forcep
"Fish-Hook" Serrations
Upper Atracumatic Extraction



RSD1-262
90° Angled Root Forcep
"Fish-Hook" Serrations
Upper Atracumatic Extraction



RSD1-263
Fig. 01
witzel
Radici inferiori
Lower roots



RSD1-264
Fig. 02
witzel
Radici superiori
Upper roots



RSD1-265
Fig. 01
Straight for pins
and fragments




RSD1-266
Fig. 00 - Lindo Levien
Straight for roots and
fragments - TC coating




RSD1-267
Fig. 01 - Lindo Levien
Curved for roots and
fragments - TC coating



RSD1-268
Fig. 02
Curved for pins
and fragments



Witzel
Upper Roots, Universal
RSD1-269
17cm



Witzel
Lower Roots, Universal
RSD1-270
16.5cm



RSD1-271
Witzel 221
Universal Pattern For
Upper Teeth



RSD1-272
Upper Roots, Universal
400 (300)
16.5cm



RSD1-273
Lower Roots, Universal
401 (301)
16cm



RSD1-274
Witzel 223
Universal Pattern For
Lower Teeth




RSD1-274-A
Root Fragment Forceps Set For 3 Pcs

- 01 Root Fragment Forceps
- 02 Root Fragment Forceps
- 03 Root Fragment Forceps



RSD1-275
With Diamond
13cm



RSD1-276
With Diamond
13cm



RSD1-277
Ralk
Root Fragments



RSD1-278
13cm



RSD1-279
13cm



RSD1-280
12.5cm

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzere per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-281
W1
Winter
upper centrals
and canines




RSD1-282
1
upper incisors
and roots



RSD1-283
2
upper and
lower roots



RSD1-284
4
upper and
lower premolars



RSD1-285
5
Physick
lower third molars



RSD1-286
6
separating forceps
spatulate

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana


Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-287
7
upper roots




RSD1-288
8
upper and lower
third molars



RSD1-289
10S
upper roots



RSD1-290
10 H
upper molars



RSD1-291
13
upper incisors



RSD1-292
14
lower incisors


EXTRACTING FORCEPS (American Pattern)

 Щипцы экстракционные (американская форма)
 Davigers - forme américaine
 Fórceps de extracción - forma americana


 Щипці екстракційні (американська форма)
 Pinzeper per estrazione - forma americana
 استخراج المنطق (النمط الأمريكي)



RSD1-293
15
Harris
 lower molars



RSD1-295
16
 lower molars




RSD1-296
17
 lower molars



RSD1-294
16S
 lower molars
 for children



RSD1-297
18 L
Harris
 upper molars,
 left



RSD1-298
18R
Harris
 upper molars,
 right

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)


Daviers - forme américaine

Férceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Férceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzere per estrazione - forma americana


ستخراج المنقط (النمط الأمريكي)



RSD1-305
30S
upper teeth,
for children




RSD1-306
32 Parmly
upper premolars
and molars



RSD1-307
32A Parmly
upper premolars
and roots



RSD1-308
30S
lower anteriors,
and roots



RSD1-309
35
upper roots



RSD1-310
36
lower anteriors,
and roots

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)


Daviers - forme américaine

Férceps de extracción - forma americana


Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)




RSD1-311
37
lower anteriors
and roots



RSD1-313
53 R
upper molars,
right



RSD1-314
53 L
upper molars,
left



RSD1-312
38
upper anteriors
and roots



RSD1-315
62
incisors, premolars,
deciduous teeth



RSD1-316
65
upper incisors
and roots

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana


ستخراج المنقط (النمط الأمريكي)



RSD1-317
69 Tomas
upper and
lower roots



RSD1-318
85 Tomas
lower canines,
premolars and molars



RSD1-319
88R Nevius
upper molars,
right



RSD1-320
85A Tomas
lower molars,
premolars and molars



RSD1-321
88L Nevius
upper molars,
left

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)


Daviers - forme américaine

Férceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-322
99A
Kells
upper incisors,
canines, premolars

RSD1-324
99C
Kells
upper incisors,
canines, premolars

RSD1-325
101
Hull
Deciduous teeth

RSD1-326
103
lower incisors,
premolars

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



**RSD1-327
#150
Cryer**
upper incisors,
premolars,
roots



**RSD1-328
#150A
Cryer**
upper incisors,
premolars,
roots



**RSD1-329
#150S
Cryer**
upper teeth,
for children



**RSD1-330
Fig. 150X**
Upper premolars,
incisors and roots



**RSD1-331
Fig. 150XAS**
Upper premolars,
incisors and
roots -
also for children



**RSD1-332
#151A
Cryer**
lower incisors,
premolars,
roots



**RSD1-333
#151 S
Cryer**
lower incisors,
premolars,
roots



**RSD1-334
#151
Cryer**
lower teeth,
for children

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)


Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-341
222
lower third molars



RSD1-342
286
upper incisors
and roots



RSD1-343
287 Swing
lower molars



RSD1-344
288
Lower canines
and premolars

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana


Щипці екстракційні (американська форма)

Pinzoper per estrazione - forma americana


ستخراج المنقط (النمط الأمريكي)



RSD1-345
Cohen
universal forceps



RSD1-346
10 S Nevius
upper canines



RSD1-347
1 A Henahan
upper incisors



RSD1-348
10 H Nevius
upper canines



RSD1-349
4 Standard
upper third molars

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Férceps de extracción - forma americana


Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-350
FH Woodward
Lower molars



RSD1-351
FS Woodward
Lower molars



RSD1-352
89
Upper molars,
right



RSD1-353
90
Upper molars,
left

EXTRACTING FORCEPS (American Pattern)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-354
23 S
Lower molars




RSD1-355
17 SK
Lower molars



RSD1-356
150 SK
Upper anteriors
and roots



RSD1-357
151 SK
lower anteriors
and roots



RSD1-358
23 SK
Lower molars



RSD1-359
Fig. 49S
Extracting Forceps
For Lower Molars
Spacial Design



RSD1-360
Fig. 51E
Extracting Forceps
For Upper Molars
Spacial Design



RSD1-361
Fig. 51ES
Extracting Forceps
For Upper Molars
Spacial Design



EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)

Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج المنقط (مع تركيب مقابض)



RSD1-362
1
upper laterals
and canines



RSD1-363
2
upper laterals
and canines




RSD1-364
3
upper incisors
and canines



RSD1-365
4
lower incisors
and canines



RSD1-366
7
upper premolars



RSD1-367
13
lower premolars



RSD1-368
13 S
lower premolars,
for children

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)

Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج المنقط (مع تركيب مقابض)



RSD1-369
17
upper molars,
right



RSD1-370
18
upper molars,
left



RSD1-371
18 A
upper molars



RSD1-372
18 B
upper molars



RSD1-373
19
upper third molars



RSD1-374
22
lower molars



RSD1-375
22 S
lower molars,
for children

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)

Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج الملقط (مع تركيب مقابض)



RSD1-376
29
upper roots




RSD1-377
29 S
upper roots,
for children



RSD1-378
30
upper roots



RSD1-379
30 S
upper roots,
for children



RSD1-380
33
lower roots



RSD1-381
33 A
lower roots



RSD1-382
33 B
lower roots

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)

Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج الملقط (مع تركيب مقابض)



RSD1-383
33 L
lower roots



RSD1-384
33 S
lower roots,
for children



RSD1-385
37
upper incisors,
for children



RSD1-386
39
upper Molars,
for children



RSD1-387
39 A
upper molars,
for children



RSD1-388
39 R
upper molars right,
for children



RSD1-389
39 L
upper molars left,
for children

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)


Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج الملقط (مع تركيب مقابض)



RSD1-390
41
upper roots




RSD1-391
44
upper roots



RSD1-392
45
lower roots



RSD1-393
48
lower roots



RSD1-394
51
upper roots



RSD1-395
51 A
upper roots



RSD1-396
51 S
upper roots,

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)

Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج المنقط (مع تركيب مقابض)



RSD1-397
51 L
upper roots



RSD1-398
52
upper roots



RSD1-399
54
for separating
upper molars



RSD1-400
65 R
upper molars,
right



RSD1-401
65 L
upper molars, left



RSD1-402
67
upper third molars



RSD1-403
67 A
upper third molars

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)


Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج المنقط (مع تركيب مقابض)



RSD1-404
67 1/2 R
upper third molars,
right



RSD1-405
67 1/2 L
lower third molars,
left.



RSD1-406
73
lower molars



RSD1-407
74
lower roots



RSD1-408
74 M
lower roots



RSD1-409
74 N
lower roots



RSD1-410
74 XN
lower roots and crowded
incisors

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)

Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج المنقط (مع تركيب مقابض)



RSD1-411
75
lower premolars



RSD1-412
76
upper roots




RSD1-413
76N
upper roots



RSD1-414
81
lower premolars
and roots



RSD1-415
79A
lower third molars



RSD1-416
79
lower third molars

EXTRACTING FORCEPS (Anatomical Handle)

Щипцы экстракционные (с эргономичной формой ручки)

Daviers (avec manche d'ajustage)

Zahnzangen (mit Paßformgriff)

Щипці екстракційні (з ергономічною формою ручки)

Forceps de Extracción (con mango adaptado a la mano)

استخراج المنقط (مع تركيب مقابض)



RSD1-417
86 A
lower molars



RSD1-418
86 B
lower molars



RSD1-419
86 C
lower molars



RSD1-420
99
lower molars



RSD1-421
99 1/4
lower molars



RSD1-422
136
upper premolars
and roots



RSD1-423
143
lower molars

EXTRACTING FORCEPS (American Pattern Hole Handle)

 Щипцы экстракционные (американская форма)
 Daviers - forme américaine
 Fórceps de extracción - forma americana

 Щипці екстракційні (американська форма)
 Pinzere per estrazione - forma americana
 استخراج المنقط (النمط الأمريكي)



RSD1-424
 Fig. 1
 Upper Anterior



RSD1-425
 Lower Molar
 Fig. 23 Cow Horn



RSD1-426
 Fig. 88L
 Upper Molar Left



RSD1-427
 Fig. 88R
 Upper Molar Right



RSD1-428
 Fig. 150
 Upper Universal

EXTRACTING FORCEPS (American Pattern Hole Handle)

Щипцы экстракционные (американская форма)

Daviers - forme américaine

Fórceps de extracción - forma americana

Щипці екстракційні (американська форма)

Pinzeper per estrazione - forma americana

ستخراج المنقط (النمط الأمريكي)



RSD1-429
Fig. 150S
Upper Universal



RSD1-430
Fig. 151
Lower Universal



RSD1-431
Fig. 151S
Lower Universal



RSD1-432
Fig. 106
Lower Incisors



RSD1-433
Fig. 174
Extra Narrow Lower Anterior
Apical Retention Forceps

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)


Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-434
Meissner
Extracting Forceps



RSD1-435
Meissner
Extracting Forceps



RSD1-436
Meissner
Extracting Forceps



RSD1-437
Comellas Forceps
Anatomical Handle

EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-438

Extracting Forceps

Set of 10 Pieces

Extracting Forceps Fig. 2

Extracting Forceps Fig. 7

Extracting Forceps Fig. 13

Extracting Forceps Fig. 17

Extracting Forceps Fig. 18

Extracting Forceps Fig. 22

Extracting Forceps Fig. 33

Extracting Forceps Fig. 51A

Extracting Forceps Fig. 67A

Extracting Forceps Fig. 79



RSD1-439

Extracting Forceps

Set of 10 Pieces

Extracting Forceps Fig. 2

Extracting Forceps Fig. 7

Extracting Forceps Fig. 13

Extracting Forceps Fig. 17

Extracting Forceps Fig. 18

Extracting Forceps Fig. 22

Extracting Forceps Fig. 33

Extracting Forceps Fig. 51A

Extracting Forceps Fig. 67A

Extracting Forceps Fig. 79



EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)

Daviers (Forme Anglaise)

Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)

Forceps de Extracción (Forme Inglesa)

استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-440

Extracting Forceps (English Pattern)

Set of 10 Pieces

Extracting Forceps Fig. 1
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33A
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79



RSD1-441

Extracting Forceps (Anatomical Handle)

Set of 10 Pieces

Extracting Forceps Fig. 1
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33A
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79



RSD1-442

Extracting Forceps (Hole Pattern Handle)

Set of 10 Pieces

Extracting Forceps Fig. 1
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33A
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79

EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-443

Extracting Forceps (English Pattern)

Set of 13 Pieces

Extracting Forceps Fig. 2
 Extracting Forceps Fig. 7
 Extracting Forceps Fig. 13
 Extracting Forceps Fig. 17
 Extracting Forceps Fig. 18
 Extracting Forceps Fig. 22
 Extracting Forceps Fig. 33
 Extracting Forceps Fig. 51A
 Extracting Forceps Fig. 67A
 Extracting Forceps Fig. 79
 Extracting Forceps Fig. 86C
 Extracting Forceps Fig. 89
 Extracting Forceps Fig. 90



RSD1-444

Extracting Forceps (English Pattern)

Set of 10 Pieces

Extracting Forceps Fig. 13S
 Extracting Forceps Fig. 22S
 Extracting Forceps Fig. 30S
 Extracting Forceps Fig. 33S
 Extracting Forceps Fig. 37
 Extracting Forceps Fig. 38
 Extracting Forceps Fig. 39
 Extracting Forceps Fig. 39L
 Extracting Forceps Fig. 39R
 Extracting Forceps Fig. 51S

EXTRACTING FORCEPS (English Pattern)

 Щипцы экстракционные (английская форма)
 Dabiers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



Fig. 1
upper
laterals
and canines



RSD1-445
 Set of 10 Pieces
 Fig. 1
 Fig. 7
 Fig. 13
 Fig. 17
 Fig. 18
 Fig. 22
 Fig. 51A
 Fig. 67
 Fig. 74N
 Fig. 79

Fig. 7
upper premolars



Fig. 13
lower premolars



Fig. 17
upper molars, right



Fig. 18
upper molars, left



Fig. 22
lower roots



Fig. 51A
upper roots



Fig. 67
upper third molars



Fig. 74N
lower molar



Fig. 79
lower third molars



EXTRACTING FORCEPS

 Щипцы экстракционные (английская форма)
 Dangers (Forme Anglaise)
 Zahnzangen (Englische Form)

 Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-446

Extracting Forceps English adult Adult Set of 10 Pcs

	Fig.
01 Tooth Forceps	1
02 Tooth Forceps	7
03 Tooth Forceps	13
04 Tooth Forceps	17
05 Tooth Forceps	18
06 Tooth Forceps	22
07 Tooth Forceps	33
08 Tooth Forceps	51
09 Tooth Forceps	67
10 Tooth Forceps	79



RSD1-447

Tooth Forceps Set of 8 Pcs

01 Tooth Forceps
02 Tooth Forceps
03 Tooth Forceps
04 Tooth Forceps
05 Tooth Forceps
06 Tooth Forceps
07 Tooth Forceps
08 Tooth Forceps



RSD1-448

Extracting Forceps English adult Adult Set of 10 Pcs

	Fig.
01 Tooth Forceps	1
02 Tooth Forceps	7
03 Tooth Forceps	13
04 Tooth Forceps	17
05 Tooth Forceps	18
06 Tooth Forceps	22
07 Tooth Forceps	33
08 Tooth Forceps	51
09 Tooth Forceps	67
10 Tooth Forceps	79



RSD1-449

Tooth Forceps Set For 4 Pcs

01 Tooth Forceps
02 Tooth Forceps
03 Tooth Forceps
04 Tooth Forceps

EXTRACTING FORCEPS

Щипцы экстракционные (английская форма)
 Daviers (Forme Anglaise)
 Zahnzangen (Englische Form)

Щипці екстракційні (англійська форма)
 Forceps de Extracción (Forme Inglesa)
 استخراج الملقط (نمط اللغة الإنجليزية)



RSD1-450
Set of 37 Pcs

- | | | |
|---------------------------------|-----------------|---|
| 01 - Extracting Forceps | Fig.1 | 1 |
| 02 - Extracting Forceps | Fig.7 | 1 |
| 03 - Extracting Forceps | Fig.13 | 1 |
| 04 - Extracting Forceps | Fig.17 | 1 |
| 05 - Extracting Forceps | Fig.18 | 1 |
| 06 - Extracting Forceps | Fig.18A | 1 |
| 07 - Extracting Forceps | Fig.19 | 1 |
| 08 - Extracting Forceps | Fig.20 | 1 |
| 09 - Extracting Forceps | Fig.22 | 1 |
| 10 - Extracting Forceps | Fig.29 | 1 |
| 11 - Extracting Forceps | Fig.33 | 1 |
| 12 - Extracting Forceps | Fig.51 | 1 |
| 13 - Extracting Forceps | Fig.51A | 1 |
| 14 - Extracting Forceps | Fig.67A | 1 |
| 15 - Extracting Forceps | Fig.79 | 1 |
| 16 - Extracting Forceps | Fig.136 | 1 |
| 17 - Extracting Forceps | Fig.137 | 1 |
| 18 - Extracting Forceps | Fig.53 L | 1 |
| 19 - Extracting Forceps | Fig.53 R | 1 |
| 20 - Extracting Forceps | Fig.65 | 1 |
| 21 - Root Elevators | Fig.1B | 1 |
| 22 - Root Elevators | Fig.3B | 1 |
| 23 - Root Elevators | Fig.5B | 1 |
| 24 - Root Elevators | Fig.4B SA Apexo | 1 |
| 25 - Root Elevators | Fig.5B SA Apexo | 1 |
| 26 - Cotton Tweezers | | 1 |
| 27 - Mouth Mirrors With Handles | Fig.5 | 1 |
| 28 - Explorers | | 1 |
| 29 - Cement Spatulas | | 1 |
| 30 - Filling Instruments | | 1 |
| 31 - Filling Instruments | | 1 |
| 32 - Filling Instruments | | 1 |
| 33 - Filling Instruments | | 1 |
| 34 - Gum Scissors | | 1 |
| 35 - Excavators | | 1 |
| 36 - Bone Curettes | | 1 |
| 37 - Dental Aspirated Syringes | | 1 |



RSD1-451 Available in Gold Coated Set of 25 Pcs
RSD1-452 Available in Polish Set of 25 Pcs

- | | | |
|---------------------------------|-----------------|---|
| 01 - Extracting Forceps | Fig.13 | 1 |
| 02 - Extracting Forceps | Fig.17 | 1 |
| 03 - Extracting Forceps | Fig.18 | 1 |
| 04 - Extracting Forceps | Fig.19 | 1 |
| 05 - Extracting Forceps | Fig.20 | 1 |
| 06 - Extracting Forceps | Fig.22 | 1 |
| 07 - Extracting Forceps | Fig.29 | 1 |
| 08 - Extracting Forceps | Fig.33 | 1 |
| 09 - Extracting Forceps | Fig.51 | 1 |
| 10 - Extracting Forceps | Fig.51A | 1 |
| 11 - Extracting Forceps | Fig.67A | 1 |
| 12 - Extracting Forceps | Fig.79 | 1 |
| 13 - Root Elevators | Fig.1B | 1 |
| 14 - Root Elevators | Fig.3B | 1 |
| 15 - Root Elevators | Fig.4B SA Apexo | 1 |
| 16 - Root Elevators | Fig.5B SA Apexo | 1 |
| 17 - Cotton Tweezers | | 1 |
| 18 - Mouth Mirrors With Handles | | 1 |
| 19 - Explorers | | 1 |
| 20 - Cement Spatulas | | 1 |
| 21 - Filling Instruments | | 1 |
| 22 - Filling Instruments | | 1 |
| 23 - Filling Instruments | | 1 |
| 24 - Excavators | | 1 |
| 25 - Dental Aspirated Syringes | | 1 |

ROOT ELEVATOR

FUNCTION OF ROOT ELEVATORS IN SURGERY

To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

HANDLES FOR ROOT ELEVATORS POINTS A, B, C, D, E, F, G, H

Thumb Handle



Shortblast Handle



Half Gold Handle



Round With Edge Handle



Thom with Flat Handle



Finger Thom Handle



Anatomical Handle



Eye Style Handle





FUNCTION OF ROOT ELEVATORS IN SURGERY

To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

LG Dental Root Elevators To look at one of these instruments you might mistake it for a narrow screwdriver.

That's because just like one an elevator has a handle and then a specially designed "blade" or tip portion.

How they're used.

The tip of an elevator is designed to be wedged into the ligament space between a tooth and its surrounding bone (as shown in our animation).

As the elevator is forced into and twisted around in this space, the tooth is rocked around and pressed against the walls of its socket.

This helps to both expand the shape of the socket and separate the tooth from its ligament.

As this work is continued, the tooth gradually becomes more and more mobile.

In some cases, the dentist may be able to completely remove the tooth just using their elevator.

If not, they will choose a point where they feel they have accomplished as much as they can with it and will then switch to using extraction forceps.



Area: distal



Area: mesial




ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المساعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Atraumatic Extraction Elevators & Socket Scoops
Socket Scoops For Cleaning Out Extraction Sites



RSD2-017
Fig. 4B SA
Schmeckeblie Apexo
Concave tip



RSD2-018
Fig. 5B SA
Schmeckeblie
Apexo
Concave tip



RSD2-019
Fig. 4 SA
Schmeckeblie
Apexo
Flat tip



RSD2-020
Fig. 5 SA
Schmeckeblie
Apexo
Flat tip



RSD2-021
Misch Socket Scoop
4mm Angled Left



RSD2-022
Misch Socket Scoop
4mm Angled Right



RSD2-023
Potts Elevator-
Left



RSD2-024
Potts Elevator-
Right



RSD2-025
Potts "T-Handle"
Elevator-Left



RSD2-026
Potts "T-Handle"
Elevator-Right

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Apexo Root Elevators



Apexo
Root Elevator



FUNCTION: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Luxator Root Elevators



Luxator
Root Elevator



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Bein Root Elevators



Seldin Root Elevators




ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المساعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Flohr Root Elevators




Coupland Root Elevators



Apical Root Elevators



Apical Root Elevators



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Laxutor Root Elevators

Long Type



Laxutor
Root Elevator

RSD2-073
Laxutor
1.5mm
Str.



RSD2-074
Laxutor
2mm
Str.



RSD2-075
Laxutor
3mm
Str.



RSD2-076
Laxutor
5mm
Str.



RSD2-077
Laxutor
3mm
Cvd.



FUNCTION:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-078
Laxutor
5mm
Cvd.



RSD2-079
Laxutor
3mm
Upper
Angle



RSD2-080
Laxutor
3mm
Double
Angle



RSD2-081
Laxutor
3mm
Str.



RSD2-082
Laxutor
5mm
Str.




ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Cryer Root Elevators



Miller Root Elevators



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Lindo Levian Root Elevators



Lindo Levian
Root Elevator



Friedmann Roots Elevator




Friedmann
Roots Elevator



Heidbrink Root Elevators



Heidbrink
Root Elevator



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-122
301
Apexo

For displacement of
deeply seated roots



RSD2-123
303
Apexo



RSD2-124
301W
Seldin



RSD2-125
304W
Seldin



RSD2-126
1
Heidbrink
Straight
Upper and lower anterior teeth and roots



RSD2-127
2
Heidbrink
Straight



upper root elevator



upper root elevator



RSD2-128
80



RSD2-129
81



RSD2-130
81A



RSD2-131
81E



RSD2-132
40
Heidbrink
Straight



RSD2-133
41
Heidbrink
Straight

Upper anterior roots and small teeth For displacement of anterior roots

Upper and lower anterior teeth and roots

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر



									
RSD2-134 Luxator 4mm straight	RSD2-135 Luxator 3.5mm straight	RSD2-136 Apexo 3mm Straight	RSD2-137 Apexo 2.5mm Straight	RSD2-138 Apexo 2mm Straight	RSD2-139 Luxator 4mm curved	RSD2-140 Luxator 3.5mm curved	RSD2-141 Luxator 3mm curved	RSD2-142 Luxator 2.5mm curved	RSD2-143 Luxator 2mm curved

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-144
Straight
3mm

						
RSD2-145 Curved 3mm	RSD2-146 Straight 5mm	RSD2-147 Curved 5mm	RSD2-148 Straight 2mm	RSD2-149 Curved 2mm	RSD2-150 Straight 4mm	RSD2-151 Curved 4mm

LUXATING ELEVATORS



Upper root elevator with helicoidal tip



Upper root elevator with helicoidal tip

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-152
Fig. 190



RSD2-153
Fig. 191



RSD2-154
27
27 Cryer



RSD2-155
28
28 Cryer



RSD2-156
27 S
27 S Cryer



RSD2-157
28 S
28 S Cryer



RSD2-158
710/29
(344)



RSD2-159
710/30
(345)



RSD2-160
Fig. 1
Hyllin 4,5mm



RSD2-161
Fig. 2
Hyllin 3mm



RSD2-162
Fig. 3
Hyllin 3mm



RSD2-163
Fig. 4
Hyllin 3mm



RSD2-164
Hyllin 153mm



RSD2-165
Read




RSD2-166
Read



RSD2-167
London Hospital



RSD2-168
Coleman



RSD2-169



Hyllin lower roots elevator

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Cryer root or fragment elevator

CRYER



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Seldin root or fragment elevator

SELDIN



RSD2-186
34S

Thinner tip and more narrow than the 34 shank



RSD2-187
1L

RSD2-188
1R

RSD2-189
4L

RSD2-190
4R

RSD2-191
51
Curtis

RSD2-192
52
Curtis

RSD2-193
11
Stout
3rd Molars

RSD2-194
11A
Stout
3rd Molars

Upper and lower molar roots

RSD2-195
34
Upper and lower
anterior roots

RSD2-196
2S
Lower 3rd
Molars

RSD2-197
3S
Upper and lower
anterior roots

RSD2-198
8
Crane Pick

RSD2-199
44C
Bradley
For displacement
of roots and 3rd
Molars

RSD2-200
A
Cogswell

RSD2-201
B
Cogswell

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-202

4
Schmeckebier
Apexo

Schmeckebier and Miller designs used for displacement of broken or deep-seated root tips.



RSD2-203

5
Schmeckebier
Apexo



RSD2-204

71
Miller
Apexo



RSD2-205

72
Miller
Apexo



RSD2-206

6
Potts

For impacted upper 3rd Molars without bony obstruction.



RSD2-207

7
Potts



RSD2-208

73
Miller
Apexo



RSD2-209

74
Miller
Apexo



RSD2-210

19

Shorter leverage for lower 1st and 2nd Molar roots. Also useful for premolar and molar deciduous teeth.



RSD2-211

20



lower roots elevator

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-212
38
Serrated
For removal of root fragments

RSD2-213
39
Serrated

RSD2-214
190

RSD2-215
191

RSD2-216
15
Woodward
Lower Molar roots

RSD2-217
16
Woodward
Lower Molar roots

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-218
23



ENGLISH PATTERNS



Upper root elevator with helicoidal tip



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

ERGO DYNAMIC ELEVATORS

Ergo Dynamic Elevators offer an alternative to standard size elevator handles.



RSD2-236
301
Apexo



RSD2-237
301
Apexo
Ergo Dynamic



RSD2-238
#77R Elevator
Serrated
Ergo Dynamic



RSD2-239
#B Cogswell
Elevator
Ergo Dynamic



RSD2-240
#34S
Seldin Elevator
Straight
Ergo Dynamic



RSD2-241
Luxating Elevator
Straight
3mm
Ergo Dynamic



RSD2-242
Luxating Elevator
Curved
3mm
Ergo Dynamic



RSD2-243
Luxating Elevator
Straight
5mm
Ergo Dynamic



RSD2-244
Luxating Elevator
Curved
5mm
Ergo Dynamic

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المساعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Laxutor Root Elevators Set of 9 Pcs



RSD2-245
2mm Str.



RSD2-246
3mm Str.



RSD2-247
5mm Str.



RSD2-248
3mm Cvd.



RSD2-249
5mm Cvd.



RSD2-250
3mm Angle



RSD2-251
3mm Str.



RSD2-252
3mm Str.



RSD2-253
5mm Str.



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Mesial & Distal Root Elevators



Root Elevators



RSD2-254
Fig. 33
Mesial
2 mm

RSD2-255
Fig. 34
Mesial
3 mm

RSD2-256
Fig. 35
Mesial
4 mm

RSD2-257
Fig. 36
Distal
2 mm

RSD2-258
Fig. 37
Distal
3 mm

RSD2-259
Fig. 38
Distal
4 mm

Apexo Root Elevators



RSD2-260
Fig. 301
Apexo
2.5mm



RSD2-261
Fig. 302
Apexo
2.5mm



RSD2-262
Fig. 303
Apexo
2.5mm



RSD2-263
Fig. 301W
Seldin
3mm



RSD2-264
Fig. 304W
Seldin
3mm

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المساعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Wisdom Roots Elevator



RSD2-265
Wisdom
for 18 and 38



RSD2-266
Wisdom
for 28 and 48

Roots Elevator



RSD2-267
Fig. 60



RSD2-268
Fig. 61



RSD2-269
Fig. 62



RSD2-270
Fig. 63



Luxier Roots Elevator



RSD2-271
Luxier
3mm, Straight



RSD2-272
Luxier
3mm, Curved



RSD2-273
Luxier
5mm, Straight



RSD2-274
Luxier
3mm, Angular


ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Hylin Roots Elevator



Molars Roots Elevator

Schlemmer Roots Elevator



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Medan Roots Elevator



Medan Roots Elevator



Roots Elevator



ROOT ELEVATOR


Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.



Root Elevators



RSD2-323
Flohr 3mm

FUNCTION: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Kopp Roots Elevator



RSD2-329
Fig. 1
Kopp 3.6mm
Upper Roots



RSD2-330
Fig. 2
Kopp 3.6mm
Upper Roots



RSD2-331
Fig. 3
Kopp 4.2mm
Upper Roots



RSD2-332
Fig. 4
Kopp 4.2mm
Upper Roots



RSD2-333
Fig. 1
Kopp 3.6mm
Lower Roots



RSD2-334
Fig. 2
Kopp 3.6mm
Lower Roots



RSD2-335
Fig. 3
Kopp 4.2mm
Lower Roots



RSD2-336
Fig. 4
Kopp 4.2mm
Lower Roots



Bein
RSD2-337



Bein
RSD2-338



Bein
RSD2-339



Bein
RSD2-340



Gartner
RSD2-341



Gartner
RSD2-342



Gartner
RSD2-343

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Kopp Roots Elevator



RSD2-344
Kopp
Right



RSD2-345
Kopp
Left

Heber Roots Elevator



RSD2-346
Heber
Bayonet



RSD2-347
Heber
Bayonet

Bein Roots Elevator



RSD2-348
Bein
3mm



RSD2-349
Bein
4mm

Bein Roots Elevator



RSD2-350
Bein
2mm, curved
serrated



RSD2-351
2mm
straight
serrated



RSD2-352
3mm
straight



RSD2-353
4mm
Curved



RSD2-354
Bein
2mm, curved
serrated



RSD2-355
2mm
straight
serrated



RSD2-356
3mm
straight



RSD2-357
4mm
Curved



RSD2-358
Bein
2mm, curved
serrated



RSD2-359
2mm
straight
serrated



RSD2-360
3mm
straight



RSD2-361
4mm
Curved

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Atraumatic Extraction Bayonet Root Elevators



RSD2-362
3mm Luxating
Elevator Curved



RSD2-363
3mm Luxating
Elevator Straight



RSD2-364
5mm Luxating
Elevator Curved



RSD2-365
5mm Luxating
Elevator Straight



RSD2-366
2mm Grooved
Luxating Elevator



RSD2-367
3mm Grooved
Luxating Elevator



RSD2-368
4mm Grooved
Luxating Elevator



RSD2-369
5mm Grooved
Luxating Elevator

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Luxating Elevator With Periotome Tip



RSD2-370
3mm Luxating
Elevator / Periotome



RSD2-371
4mm Luxating
Elevator / Periotome



RSD2-372
5mm Luxating
Elevator / Periotome

Atraumatic Extraction Bayonet Root Elevators



RSD2-373
Bayonet Elevator
Straight



RSD2-374
Bayonet Elevator
Offset "In" Small



RSD2-375
Bayonet Elevator
Offset "Out" Small



RSD2-376
Bernard Bayonet
Elevator Offset Left



RSD2-377
Bernard Bayonet
Elevator Offset Right

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

W. Barry Roots Elevator



W. Barry
Fig. 1
(125mm)



RSD2-378
W. Barry
Fig. 1



RSD2-379
W. Barry
Fig. 2



RSD2-380
W. Barry
Fig. 3



RSD2-381
W. Barry
Fig. 4

Winter Roots Elevator



Winter
Fig. 1L
(125mm)



RSD2-382
Winter
Fig. 1L



RSD2-383
Winter
Fig. 1R



RSD2-384
Winter
Fig. 11L



RSD2-385
Winter
Fig. 11R



RSD2-386
Winter
Fig. 12L



RSD2-387
Winter
Fig. 12R

Quality Symbol

Pott Roots Elevator



RSD2-388
Winter
Fig. 1L



RSD2-389
Winter
Fig. 1R



RSD2-390
Lecluse
bayonet



RSD2-391
Pott
Left



RSD2-392
Pott
Right

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

Heidbrink Roots Elevator

Syndesmotomes Roots Elevator

Heidebrink Root tip Picks

8mm

RSD2-397
Fig. 2
Left

RSD2-398
Fig. 3
Right

RSD2-393

RSD2-394

RSD2-395

RSD2-396
Fig. 1

RSD2-399
Straight

RSD2-400
Left

RSD2-401
Right

Periotome Roots Elevator

Syndesmotomes

RSD2-402
Periotome
2.1mm
Straight

RSD2-403
Periotome
2.6mm
Straight

RSD2-404
Periotome
2.1mm
Angular

RSD2-405
Periotome
2mm
Curved

RSD2-406
Periotome
2.0mm

RSD2-407
Periotome
2.5mm

RSD2-408
Periotome
3.0mm

RSD2-409
Periotome
4.0mm

ROOT ELEVATOR


Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To loosen small root fragments from bony socket

Heidbrink Roots Elevator

Warwick-James



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.




Heidbrink universal mini-elevator for fragments

HOWARD

Slightly heavier design than Heidbrink picks.

HEIDBRINK

APICAL



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

FUNCTION: To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر



Lecluse elevator for lower third molar dislocation

CROSS-BAR HANDLE

RSD2-448
Fig. 2
Chompret

RSD2-449
Fig. 3
Chompret

RSD2-450
Fig. 4
Chompret



Chompret syndesmotome, sickle pattern



RSD2-451
11LX



RSD2-452
11RX

RSD2-447
Fig. 1
Chompret



RSD2-453
12LX



RSD2-454
12RX



RSD2-455
14LX



RSD2-456
14RX



RSD2-457
1LX



RSD2-458
1RX

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



RSD2-459
775/3
Winter
11 L



RSD2-460
775/4
Winter
11 R



RSD2-461
775/5
Winter
12 L



RSD2-462
775/6
Winter
12 R



RSD2-463
775/7
Winter
13 L



RSD2-464
775/8
Winter
13 R



RSD2-465
winter
Left



RSD2-466
winter
right



RSD2-467
775/1
Winter
1 L



RSD2-468
775/2
Winter
1 R



RSD2-469
Winter
1L




RSD2-470
Winter
1R

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

LUXACUT



Luxacut Luxation Instrument

Over recent years, there has been an increased demand for reduced trauma during extraction procedures. The use of elevators has become widespread for the loosening of teeth before the application of forceps. Luxation instruments tips are more finely tapered than elevators tips. Their characteristics are a combination of a Periosteal and an elevator. The sharp tip enables more efficient cutting of the periodontal ligament, whilst the tapering shaft assists with controlled expansion of the bony socket. The extraction operation may be performed with the minimum amount of tissue damage. Healing is therefore quicker and the operation is less distressing for the patient and surgical team alike.


The Luxacut Kit comprises of seven Luxation instruments, a straight, curved, inverted curved & contra blade version in Three different widths. Also available in dual edge blades

- RSD2-471-L2S Straight Blade- Green-2 mm- Apical and/or Interproximal
- RSD2-473-L3C Curved Blade- Black- 3 mm General use/Starter
- RSD2-472-L3S Straight Blade- Grey- 3 mm Interproximal
- RSD2-475-L3CA Contra Angle- Blue- 3 mm Lingual and/or distal-Molars
- RSD2-474-L3IC Inverted Curved- Light Blue- 3 mm- General, Lingual and/or Distal
- RSD2-477-L5C Curved Blade- Brown- 5 mm- Large molar roots- General use
- RSD2-476-L5S Straight Blade- Yellow- 5 mm- Interproximal
- RSD2-478-DE3 dual edge/straight grey 3 mm/1.5 mm
- RSD2-479-DE5 dual edge/straight brown 5 mm Deep fractured/Decayed Roots

Luxacut Kit



RSD2-480
LWC7



RSD2-481
LGWC7

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

RSD2-482

LUXATING ELEVATORS SET OF 10 PCS



RSD2-493

LUXATING ELEVATORS SET OF 10 PCS - BLUE TITANIUM COATED TIPS



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

RSD2-504

LUXATING ELEVATORS SET OF 10 PCS - GOLD COATED TIPS



RSD2-515

LUXATING ELEVATORS SET OF 08 PCS - BLUE TITANIUM COATED TIPS



ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

RSD2-524

PERIOTOME LUXATING ELEVATORS SET OF 04 PCS - GOLD COATED TIPS



RSD2-525
PERIOTOME GOLD TIP - L3S



RSD2-526
PERIOTOME GOLD TIP - De3



RSD2-527
PERIOTOME GOLD TIP - L2S



RSD2-528
PERIOTOME GOLD TIP - L3CA



RSD2-529 (Double Curved 1.5mm)
RSD2-530 (Double Curved 2mm)
RSD2-531 (Double Curved 3mm)
RSD2-532 (Double Curved 4mm)
RSD2-533 (Double Curved 5mm)



RSD2-534
(Inverted Curved 3mm)



RSD2-535
(Dual Edge 3mm / 1.5mm)



RSD2-536
(Dual Edge 5mm / 3mm)

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر



RSD2-537


Luxators (White Handle) Set of 10 Pcs

01	1.5mm Straight	06	1.5mm Curved
02	2mm Straight	07	2mm Curved
03	3mm Straight	08	3mm Curved
04	4mm Straight	09	4mm Curved
05	5mm Straight	10	5mm Curved

RSD2-538

Luxators Set of 8 Pcs

01	Luxator
02	Luxator
03	Luxator
04	Luxator
05	Luxator
06	Luxator
07	Luxator
08	Luxator



RSD2-539

Luxators Set of 7 Pcs

01	Luxator
02	Luxator
03	Luxator
04	Luxator
05	Luxator
06	Luxator
07	Luxator

RSD2-540

Luxators Set of 8 Pcs

01	Luxator
02	Luxator
03	Luxator
04	Luxator
05	Luxator
06	Luxator
07	Luxator
08	Luxator



RSD2-541

Luxators Set of 8 Pcs

01	Luxator
02	Luxator
03	Luxator
04	Luxator
05	Luxator
06	Luxator
07	Luxator
08	Luxator

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر



RSD2-542

Root Elevators Set of 6 Pcs

- 01 Bein 3mm
- 02 Bein 5mm
- 03 Cryer 31
- 04 Cryer 32
- 05 Apical 302
- 06 Apical 303



RSD2-543

Root Elevators Set of 6 Pcs

- 01 Bein 3mm
- 02 Bein 5mm
- 03 Cryer 31
- 04 Cryer 32
- 05 Flohr Fig. 2
- 06 Flohr Fig. 3

RSD2-544
Elevators Set of 7 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator

RSD2-545
Elevators Set of 8 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator
- 08 Elevator

RSD2-546
Root Elevators Set of 10 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator
- 08 Elevator
- 09 Elevator
- 10 Elevator

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

RSD2-547
Elevators Set of 10 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator
- 08 Elevator
- 09 Elevator
- 10 Elevator

RSD2-548
Elevators Set of 13 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator
- 08 Elevator
- 09 Elevator
- 10 Elevator
- 11 Elevator
- 12 Elevator
- 13 Elevator

RSD2-549
Elevators Set of 8 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator
- 08 Elevator

RSD2-550
Elevators Set of 8 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator
- 08 Elevator

RSD2-551
Elevator Set of 8 Pcs



- 01 Elevator
- 02 Elevator
- 03 Elevator
- 04 Elevator
- 05 Elevator
- 06 Elevator
- 07 Elevator
- 08 Elevator

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

PDL ELEVATORS

Distal, 2.5mm
Micro Edge

Mesial, 2.5mm
Micro Edge

Distal, 4mm
Micro Edge

Mesial, 4mm
Micro Edge

2.5mm
Micro Edge, curved

2.5mm
Micro Edge, Straight

3.4mm
Straight

RSD2-552
PDL Elevators
Set of 7 Pcs

Distal, 4mm
Micro Edge

Mesial, 4mm
Micro Edge

Distal, 2.5mm
Micro Edge

Mesial, 2.5mm
Micro Edge

2.5mm
Micro Edge, curved

2.5mm
Micro Edge, Straight

3.4mm
Straight

RSD2-553
PDL Elevators
Set of 7 Pcs

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر




Area: distal



Area: mesial



RSD2-554 PDL Elevators Set of 10 Pcs



RSD2-555
Luxa-Tool, straight, 3.0 mm,

RSD2-556
Luxa-Tool, curved, 3.0 mm,

RSD2-557
Proximotor, curved, 2.5 mm,
distal elevator,

RSD2-558
Proximotor, curved, 2.5 mm,
mesial elevator,

RSD2-559
Proximotor, curved, 4.5 mm,
distal elevator,

RSD2-560
Proximotor, curved, 4.5 mm,
mesial elevator,

RSD2-561
Syndesmotome, curved, 4.0 mm,
for distal access

RSD2-562
Syndesmotome, curved, 4.0 mm,
for mesial access

RSD2-563
Syndesmotome, curved left,
4.0 mm, for lingual/buccal luxation,

RSD2-564
Syndesmotome, curved right,
4.0 mm, for lingual/buccal luxation,

ROOT ELEVATOR

Элеватор корневой экстракционный
Elévateurs à racines
Wurzelheber

Елеватор корневий екстракційний
Elevadores para raíces
المصاعد الجذر

PDL ELEVATORS



RSD2-565
PDL Elevators
Set of 8 Pcs

ROOT ELEVATOR

Элеватор корневой экстракционный
 Elévateurs à racines
 Wurzelheber

Елеватор корневий екстракційний
 Elevadores para raíces
 المصاعد الجذر

RSD2-566

Pdl Elevator Instrument Set of 7 Pcs



- 01 Pdl Elevator
- 02 Pdl Elevator
- 03 Pdl Elevator
- 04 Pdl Elevator
- 05 Pdl Elevator
- 06 Pdl Elevator
- 07 Pdl Elevator

RSD2-567

Pdl Elevator Set of 7 Pcs



- 01 Pdl Elevator
- 02 Pdl Elevator
- 03 Pdl Elevator
- 04 Pdl Elevator
- 05 Pdl Elevator
- 06 Pdl Elevator
- 07 Pdl Elevator



RSD2-568
 PDL Elevator Set of 13 Pcs

- 01 PDL Elevator
- 02 PDL Elevator
- 03 PDL Elevator
- 04 PDL Elevator
- 05 PDL Elevator
- 06 PDL Elevator
- 07 PDL Elevator
- 08 PDL Elevator
- 09 PDL Elevator
- 10 PDL Elevator
- 11 PDL Elevator
- 12 PDL Elevator
- 13 PDL Elevator



RSD2-569
 PDL Elevator Set of 14 Pcs

- 01 PDL Elevator
- 02 PDL Elevator
- 03 PDL Elevator
- 04 PDL Elevator
- 05 PDL Elevator
- 06 PDL Elevator
- 07 PDL Elevator
- 08 PDL Elevator
- 09 PDL Elevator
- 10 PDL Elevator
- 11 PDL Elevator
- 12 PDL Elevator
- 13 PDL Elevator
- 14 PDL Elevator

PERIOSTEALS ELEVATOR 10mm Hollow Handle



RSD2-570
Raspator,
Bone Scraper
Length 16.5cm,
10 mm hollow
handle,
double-ended 6
mm / 4 mm,

RSD2-571
Raspator,
Raspatory
Length 18cm,
10 mm hollow
handle,
double-ended
3 mm | 3 mm,

RSD2-572
Raspator,
Raspatory
Length 18cm,
10 mm hollow
handle,
double-ended
3 mm | 3 mm,

RSD2-573
Raspator,
Papilla
Length 18.5cm,
10 mm hollow
handle,
double-ended
3 mm | 2 mm,

RSD2-574
Raspator,
Bone Scraper
Length 17.5cm,
10 mm hollow
handle,
double-ended,

RSD2-575
Raspator,
Micro-Freer
Length 18cm,
10 mm hollow
double-ended
3 mm | 4 mm

RSD2-576
Raspator,
Freer
Length 18cm,
10 mm hollow
handle,
double-ended
5 mm | 5 mm,

PERIOSTEALS ELEVATOR 10mm Hollow Handle



RSD2-577
Raspator,
Freer
10 mm hollow
handle,
high-quality
medical
double ended

RSD2-578
Raspator,
Molt
Length 18cm,
10 mm Hollow
grip,
double-ended
6 mm

RSD2-579
Raspator,
Micro Freer
Length 18.5cm,
10 mm hollow
handle,
double-ended

RSD2-580
Raspator,
Micro Freer
Length 18.5cm,
10 mm hollow
handle,
double-ended

RSD2-581
Raspator, Micro
Raspator /
Chisel
Length 17cm,
10 mm hollow
handle,
double-ended

RSD2-582
Raspator, Micro
Raspator /
Chisel
Length 17.5cm,
10 mm hollow
handle,
double-ended

RSD2-583
Raspator,
Micro Raspator /
Chisel
Length 17.5cm,
10 mm hollow handle,
double-ended
3.5mm | 3.5 mm

PERIOSTEALS ELEVATOR 10mm Hollow Handle



RSD2-584
Raspatory Micro Molt
Length 17.5cm,
corrosion-resistant

RSD2-585
Raspatory Chipping
Workers:
Chisel 24G and 7.5 mm,
length about 17.5cm,
corrosion-resistant

RSD2-586
Raspatory Molt
Work ends:
5.0 mm and 4.0 mm,
length about 17.5cm,
corrosion-resistant

RSD2-587
Rasp, Molt M9
Length 17.5cm,
corrosion-resistant

RSD2-588
Raspator, Wautier
Length 17.5cm,
corrosion-resistant

RSD2-589
Raspator, Benque
Length 17cm,
workers
5 mm | 7 mm,
corrosion-resistant

RSD2-590
Raspator, Nordberg
Length 18cm,
10 mm hollow handle,
double-ended
6 mm | 4 mm

DIAGNOSTIC, PERIODONTAL & RESTORATIVE INSTRUMENTS



Dental Hand Instruments

The types of hand instruments can be categorized based on their function:

Examining the mouth and teeth Scaling or professional cleaning of teeth Cutting teeth and removing cavities Placing and condensing filling materials Carving and finishing dental fillings Miscellaneous types Instruments used for examining the mouth and teeth

Mouth mirrors are used to examine the condition of the teeth and mouth, as well as to retract the tongue and cheeks upon examination.

Probes— There are many probes available but the main types are the straight probe, Briault probe and periodontal probe.

The straight probe is sharp and it is used to examine cavities in the dentine region as well as to check the margins of restorations or fillings.

The Briault probe is also sharp and it is used to detect the junction between the enamel and dentine. Briault can also be used to detect tartar hidden in the gum pockets.

A periodontal probe can be used to measure the depth of the gum pockets and for other measurements. Instruments for cutting and scraping tissues and deposits

Scalers remove tartar (calculus) on the tooth surface as well as those hidden in the gum pockets. They can also be used to remove other surface deposits and temporary crowns.

Excavators are used to remove softened dentine and temporary fillings. The back of the blade can also be used for placing linings between filling materials. Sometimes excavators can be used for carving metal or amalgam fillings.

Chisels are used to remove unsupported enamel during cavity preparation. Unsupported enamel can be harmful for it can compromise the longevity of the filling if not removed. Hatchets and hoes have similar purposes to chisels. Instruments used for placing and condensing filling materials Plastic instruments convey and shape filling materials and are usually flat blades that do not use heavy pressure. Some instruments have rounded ends for shaping and burnishing the surfaces of fillings. Plastic instruments are usually made of stainless steel for most materials though there are also thin Teflon-coated or titanium nitride instruments which are non-stick for composite fillings.

Condensers or pluggers are used for compressing and forming filling materials. Heavy pressure is applied on use therefore these instruments are usually used for metal or amalgam fillings. Instruments used for carving and finishing fillings Carving and finishing instruments have sharp or semi-sharp blades of various shapes. They are used to carve materials by cutting and scraping.



DIAGNOSTIC, PERIODONTAL & RESTORATIVE INSTRUMENTS

VALUED BUYERS

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours. We regularly manufacturing below articles in showing handles & pattern A to L for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in materials, sizes, shapes and patterns.



A-Hollow Handle Ø 10 mm



B-Hollow Handle Ø 10 mm



C-Hollow Handle Ø 10 mm



D-Hollow Handle Ø 10 mm



E-Hollow Handle Ø 10 mm



F-Hollow Handle Ø 10 mm



G-Hollow Handle Ø 10 mm



H-Hollow Handle Ø 6.4 mm



I-Hollow Handle Ø 8 mm



J-Hollow Handle Ø 10 mm



K-Soild Handle Ø 6 mm



L-Soild Handle Ø 6 mm



M-Hollow Handle Ø 10 mm



N-Hollow Handle Ø 10 mm



O-Hollow Handle Ø 10 mm



P-Hollow Handle Ø 10 mm



Q-Hollow Handle Ø 10 mm



R-Solid Handle Ø 6 mm



S-Solid Handle Ø 6 mm



T-Solid Handle Ø 6 mm



U-Hollow Handle Ø 8 mm



V-Hollow Handle Ø 8 mm



W-Silicon Handle



X-Autoclavable Handle Ø 11 mm



Y-Autoclavable Handle Ø 9.5 mm

MOUTH MIRRORS

зеркал
Espejos de la boca
Mundspiegel

дзеркал
Miroirs à la bouche
مرايا الفم

Plain



RSD3-001 Fig. 3 Ø20mm
RSD3-002 Fig. 4 Ø22mm
RSD3-003 Fig. 5 Ø24mm
RSD3-004 Fig. 6 Ø26mm

Cone-Socket Plain



RSD3-005 Fig. 3 Ø20mm
RSD3-006 Fig. 4 Ø22mm
RSD3-007 Fig. 5 Ø24mm

Magnifying



RSD3-008 Fig. 3 Ø20 mm
RSD3-009 Fig. 4 Ø22 mm
RSD3-010 Fig. 5 Ø24 mm

Cone-Socket Magnifying



RSD3-011 Fig. 3 Ø20mm
RSD3-012 Fig. 4 Ø22mm
RSD3-013 Fig. 5 Ø24mm

Rhodium



RSD3-014 Fig. 3 Ø20mm
RSD3-015 Fig. 4 Ø22mm
RSD3-016 Fig. 5 Ø24 mm

Cone-Socket Rhodium



RSD3-017 Fig. 3 Ø20mm
RSD3-018 Fig. 4 Ø22mm
RSD3-019 Fig. 5 Ø24mm

Front Surface



RSD3-020 Fig. 3 Ø20mm
RSD3-021 Fig. 4 Ø22mm
RSD3-022 Fig. 5 Ø24mm
RSD3-023 Fig. 6 Ø26mm

Cone-Socket Front Surface



RSD3-024 Fig. 3 Ø20mm
RSD3-025 Fig. 4 Ø22mm
RSD3-026 Fig. 5 Ø24mm

HANDLES FOR MOUTH MIRRORS

 Ручки для зеркал
 Manches pour miroirs à bouche
 Mundspiegelgriffe

 Ручки для дзеркал
 Mangos para espejos de boca
 مقابض لمرآيا الفم



RSD3-027
Hollow Mirror
Handle Round
Ø 10 mm



RSD3-028
Hollow Mirror
Hexagonal
Handle
Ø 8 mm



RSD3-029
Hollow
Ø 8 mm



RSD3-030
Hollow Mirror
Hexagonal
Handle
Ø 10 mm



RSD3-031
Hollow Mirror
Handle
Ø 8 mm



RSD3-032
Hollow Mirror
Handle
Ø 10 mm



RSD3-033
Hollow Mirror
Handle
Ø 8 mm



RSD3-034
Mirror Handle
Autoclavable
Material
Available in
Various Colours
Ø 9.5 mm



RSD3-035
Mirror Handle
Autoclavable
Material
Available in
Various Colours
Ø 9.5 mm



RSD3-036
Mirror Handle
Silicon Grip
Autoclavable
Available in
Various Colours



RSD3-037
Mirror Handle
Solid
Ø 6 mm



RSD3-038
Mirror Handle
Solid
Ø 6 mm



RSD3-039
Mirror Handle
Solid Octagonal
Ø 6 mm



RSD3-040
Mirror Handle
Solid
Stem with
Rule 35mm



RSD3-041
Cone-Socket
Solid Octagonal
Ø 6 mm



RSD3-042
Cone-socket
Solid
Ø 6 mm



RSD3-043
Cone-socket
Hexagonal
Handle
Ø 6 mm



RSD3-491
Mouth Mirror

Cone-Socket Mirror Handle

EXAMINATION SETS

■ Комплект для обследования
 ■ Ensemble d'examen
 ■ Prüfungssatz

■ Экзаменаційний набір
 ■ Conjunto de examen
 ■ مجموعة الفحص

RSD3-412

First Aid Examination kit Set of 3 Pieces



RSD3-412-A

Examination Set of 3 Pcs



RSD3-413

First Aid Examination Set of 4 Pieces



RSD3-414

First Aid Examination Set of 3 Pieces



RSD3-415

First Aid Examination Set of 3 Pieces




RSD3-416

First Aid Examination Set of 3 Pieces



RSD3-417

First Aid Examination Set of 3 Pieces



EXAMINATION SETS

■ Комплект для обследования
 ■ Ensemble d'examen
 ■ Prüfungssatz

■ Экзаменаційний набір
 ■ Conjunto de examen
 ■ مجموعة الفحص

RSD3-418

First Aid Examination kit Set of 3 Pieces



RSD3-418-A

Examination Set of 3 Pcs



RSD3-419

First Aid Examination Set of 4 Pieces



RSD3-420

First Aid Examination Set of 3 Pieces



RSD3-421

First Aid Examination Set of 3 Pieces



RSD3-422


First Aid Examination Set of 3 Pieces



RSD3-423

First Aid Examination Set of 3 Pieces






RSD3-1303 Adult Palatal /
RSD3-1304 Child Palatal

RSD3-1305 Buccal

RSD3-1306 Adult Palatal /
RSD3-1307 Child Palatal

RSD3-1308 Adult Palatal /
RSD3-1309 Adult Buccal




RSD3-1310 Lingual upper /
RSD3-1311 Lingual Lower

RSD3-1312 Adult Buccal /
RSD3-1313 Child Buccal

RSD3-1314 Child Palatal /
RSD3-1315 Child Buccal

RSD3-1316 Palatal-Adult



RSD3-1317 Palatal-Adult Large

RSD3-1318 Palatal-Child

RSD3-1319 Buccal-Small

RSD3-1320 Buccal

RSD3-1321 Lingual

RSD3-1322 Lingual-Large

TWEEZERS

 Пинцеты
 Précelles - Bruxelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



RSD3-044
London College
15cm



RSD3-045
Meriam
16cm



Semken Taylor
RSD3-046 12.5cm
RSD3-046A 15cm
Straight



Semken Taylor
RSD3-047 12.5cm
RSD3-047A 15cm
Curved



Semken Taylor
RSD3-048 12.5cm
RSD3-048A 15cm
1x2 Straight



Semken Taylor
RSD3-049 12.5cm
RSD3-049A 15cm
1x2 Curved



RSD3-050
Flagg
15cm



RSD3-051
Meriam
16cm



RSD3-052
Amsterdam
15.5cm



RSD3-053
Perry
13cm
Curved



RSD3-054
Self-Locking
15cm
Angled



RSD3-055
Endo-Locking
Grooved
15cm

TWEEZERS

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пинцети
 Pinzas para algodón
 ملاقط



RSD3-056 9cm
RSD3-057 12cm
RSD3-058 15cm
RSD3-059 18cm
 Taylor



RSD3-060 9cm
RSD3-061 12cm
RSD3-062 15cm
RSD3-063 18cm
 Taylor



RSD3-064 0,6mm
RSD3-065 0,9mm
 Gerald
 17.5cm



RSD3-066 0,6mm
RSD3-067 0,9mm
 Gerald
 17.5cm



RSD3-068
 Gerald
 17.5cm



RSD3-069
 Gerald
 17.5cm



RSD3-070
 15cm
 College



RSD3-071
 Meriam
 15cm



RSD3-072
 De Backy
 Atraumatic
 15cm
 Straight



RSD3-073
 Meriam Tweezers
 smooth 16cm



RSD3-074
 London College
 smooth 16cm



Adson With Teeth
RSD3-075 12cm
RSD3-076 15cm



TWEEZERS

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пинцети
 Pinzas para algodón
 ملاقط



RSD3-077
Perry
10.5 cm



RSD3-078
Perry
10.5 cm



RSD3-079
Fig. 1
Perry
13 cm



Potts-Smith Dressing Forceps
RSD3-080 18 cm
RSD3-081 20 cm
RSD3-082 25 cm
 straight serrated



Potts-Smith Dressing Forceps
RSD3-083 18 cm
RSD3-084 20 cm
RSD3-085 25 cm
 curved serrated



Potts-Smith Tissue Forceps 1x2 teeth
RSD3-086 18 cm
RSD3-087 20 cm
RSD3-088 25 cm
 straight



Potts-Smith Tissue Forceps 1x2 teeth
RSD3-089 18 cm
RSD3-090 20 cm
RSD3-091 25 cm
 curved



De-Bakey With Tungsten Carbide (TC) 17.5cm
RSD3-092



Adson
RSD3-093 12cm
RSD3-094 15cm



Adson 1x2
RSD3-095 12cm
RSD3-096 15cm



Adson Brown
RSD3-097 12cm
RSD3-098 15cm



RSD3-099
Crane Kaplan
15cm



RSD3-100
Hemingway
15cm

TWEEZERS

 Пинцеты
 Précelles - Bruxelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



RSD3-101
Crane Kaplan
15cm



RSD3-102
Crane Kaplan
15cm



RSD3-103 10.5cm
RSD3-104 11.5cm
RSD3-105 13.0cm
RSD3-106 14.0cm
RSD3-107 16.0cm
RSD3-108 18.0cm
RSD3-109 20.0cm
RSD3-110 25.0cm
RSD3-111 30.0cm

Dressing Forceps



RSD3-112 10.5cm
RSD3-113 11.5cm
RSD3-114 13.0cm
RSD3-115 14.0cm
RSD3-116 16.0cm
RSD3-117 18.0cm
RSD3-118 20.0cm
RSD3-119 25.0cm
RSD3-120 30.0cm

Tissue Forceps
1x2 teeth



RSD3-121 11.5 cm
RSD3-122 13.0 cm
RSD3-123 14.0 cm
RSD3-124 16.0 cm
RSD3-125 18.0 cm
RSD3-126 20.0 cm
RSD3-127 25.0 cm
RSD3-128 30.0 cm

Tissue Forceps
2x3 Teeth




Lane Tissue Forceps




RSD3-129
12.5 cm
1x2 teeth



RSD3-130
15 cm
1x2 teeth



RSD3-131
17.5 cm
1x2 teeth




RSD3-132
15 cm
2x3 teeth




RSD3-133
17.5 cm
2x3 teeth



Adson-Brown



RSD3-134
12 cm
Tissue, 7:7 Teeth



RSD3-135
12 cm
Tissue, 9:9 Teeth



Adson-Brown




RSD3-136
12 cm
Teeth 9:9
Curved

TWEEZERS

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط


RSD3-137
 15 cm
 Straight, 9:9 Teeth



Brown


RSD3-138
 12 cm
 Dressing, Serrated



Ewald / Hudson


RSD3-140
 12 cm
 dressing, serrated



Adson /
Micro-Adson


RSD3-141
 12 cm
 dressing,
 Micro serrated


RSD3-143
 12 cm
 Tissue, Micro
 1:2 Teeth


RSD3-139
 12 cm
 Tissue, 1:2 Teeth


RSD3-142
 12 cm
 Tissue,
 1:2 Teeth


RSD3-144
 12 cm
 Dressing,
 cross serration



Adson
T.C


RSD3-145
 12 cm
 cross serration-
 1:2 Teeth


RSD3-146
 15 cm
 cross serration-
 1:2 Teeth


RSD3-147
 Adson
 1 mm
 12 cm


RSD3-148
 Adson T.C
 1,7 mm
 12 cm


RSD3-149
 Adson
 1,7 mm
 15 cm




RSD3-150
 Micro-Adson
 very fine
 0,8 mm
 15 cm



RSD3-151
 Adson
 1x2, 1,5 mm
 12 cm



TWEEZERS

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



RSD3-152
 Micro-Adson
 1x2, 1 mm
 15 cm
 very fine



RSD3-153
 Micro-Adson
 1x2, 1 mm
 15 cm
 curved
 very fine



RSD3-154
 10 cm



RSD3-156
 10 cm



RSD3-155
 10 cm



RSD3-157
 10 cm
 1x2 teeth



RSD3-159
 10 cm
 1x2 teeth



RSD3-158
 10 cm
 1x2 teeth



RSD3-160
 15 cm
 1x2 teeth

Graefe

Graefe Tissue Forceps

Mc Indoe



RSD3-161
 Cushing
 1,2 mm
 18 cm



RSD3-162 13 cm
 1x2, 1,5 mm
RSD3-163 14 cm
 1x2, 1,7 mm



RSD3-164
 Taylor
 1x2, 1,5 mm
 17,5 cm



RSD3-165 13 cm
 1x2, 2,3 mm
RSD3-166 14 cm
 1x2, 2,3 mm



RSD3-167
 Micro-Cushing
 1x2, 0,9 mm
 17,5 cm
 very fine



RSD3-168
 Suture Forceps
 15 cm

TWEEZERS

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пинцети
 Pinzas para algodón
 ملاقط



Delicate Tweezers



RSD3-169
11.5 cm



RSD3-170
13 cm



RSD3-171
14.5 cm



Gerard



RSD3-172
17.5 cm



RSD3-173
23 cm



RSD3-174
17.5 cm



Gerard 1x2



RSD3-175
17.5 cm



RSD3-176
23 cm



RSD3-177
17.5 cm



Selman
RSD3-178
20 cm



Ø 4 mm



Brown



RSD3-179
15 cm



RSD3-180
20 cm



Strassburger Modell



RSD3-181
15 cm



RSD3-182
20cm

TWEEZERS

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



Taylor Dressing Forceps
RSD3-183 18 cm
 straight serrated



Taylor Tissue Forceps
RSD3-184 18 cm
 1x2 teeth, Straight



Taylor Tissue Forceps
RSD3-185 18 cm
 1x2 teeth, Curved



Cushing Dressing Forceps
RSD3-186 18 cm
 Serrated Straight



Cushing Dressing Forceps
RSD3-187 18 cm
 Serrated Curved



Cushing Tissue Forceps
RSD3-188 20 cm
 1x2 teeth Straight



Cushing Tissue Forceps
RSD3-189 18 cm
 1x2 teeth Curved



Brophy Dressing Forceps
RSD3-190 20 cm
 straight serrated



Brophy Dressing Forceps
RSD3-191 20 cm
 curved serrated



Brophy Tissue Forceps
RSD3-192 20 cm
 1x2 teeth straight serrated



Brophy Tissue Forceps
RSD3-193 20 cm
 1x2 teeth curved serrated

TWEEZERS

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пинцети
Pinzas para algodón
ملقط



Gillies Tissue Forceps
1X2 Teeth
RSD3-194 15 cm
straight serrated



Gillies Tissue Forceps
1X2 Teeth
RSD3-195 15 cm
curved serrated



Iris Tissue Forceps
RSD3-196 10 cm
1x2 Teeth, Straight



Treves Tissue Forceps
1x2 Teeth
RSD3-199 12.5 cm
RSD3-200 15.0 cm
RSD3-201 18.0 cm



Waugh Tissue Forceps
1x2 Teeth
RSD3-202 15 cm
RSD3-203 18 cm



Charnley Tissue Forceps
1x2 Teeth
RSD3-204 14.5 cm
RSD3-205 18.0 cm
block on jaws.



Dissecting Dressing
Forceps
RSD3-206 12.5 cm
RSD3-207 15.0 cm
RSD3-208 18.0 cm



Dissecting Dressing
Forceps
RSD3-209 12.5 cm
RSD3-210 15.0 cm
RSD3-211 18.0 cm
RSD3-212 20.0 cm
block end, serrated
jaws.

TWEEZERS

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط




Ferris-Smith Tissue Forceps 1x2 Teeth
 RSD3-213 12 cm
 RSD3-214 14 cm
 RSD3-215 16 cm
 RSD3-216 18 cm
 RSD3-217 20 cm




Ferris-Smith Tissue Forceps 2x3 Teeth
 RSD3-218 12 cm
 RSD3-219 14 cm
 RSD3-220 16 cm
 RSD3-221 18 cm
 RSD3-222 20 cm



Bonney Tissue Forceps
 RSD3-223 17 cm



RSD3-224
 18 cm
 1x2 Teeth



RSD3-225
 18 cm
 2x3 Teeth



Durante Tissue Grasping Forceps
 RSD3-226 14.5cm
 RSD3-227 16 cm
 RSD3-228 20 cm




Russian-Modell Tissue Forceps
 RSD3-229 15 cm
 RSD3-230 20 cm
 RSD3-231 25 cm



Singley-Tuttle Forceps
 RSD3-232 23 cm



Collin-Duval Forceps



TWEEZERS

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



Lucae

RSD3-236
 14 cm
 Dressing, Serrated
 RSD3-237
 16 cm
 Dressing, Serrated
 RSD3-238
 18 cm
 Dressing, Serrated
 RSD3-239
 20 cm
 Dressing, Serrated



Lucae

1x2 teeth
 RSD3-240
 14 cm
 RSD3-241
 16 cm
 RSD3-242
 18 cm
 RSD3-243
 20 cm



Adson Bayonet
 Dressing Forceps
 RSD3-244 21 cm
 serrated



Yasargil
 RSD3-245
 20 cm



Hunt
 RSD3-246
 20 cm



Hunt
 RSD3-247
 20 cm

TWEEZERS

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



Jansen
RSD3-248
 18.5 cm



Jansen
RSD3-249
 16 cm



Adson
RSD3-250
 23 cm



Troeltsch
RSD3-251
 12 cm



Troeltsch
RSD3-252
 11.5 cm
RSD3-253
 14 cm
RSD3-254
 17 cm



Landolt
RSD3-255
 21.5 cm

TWEEZERS

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملقط



2.7 mm

Debakey Atraumatic Dressing Forceps
RSD3-256 15cm
RSD3-257 20cm
RSD3-258 24cm
RSD3-259 30cm



3.5 mm

Debakey Atraumatic Dressing Forceps
RSD3-260 15cm
RSD3-261 20cm
RSD3-262 24cm
RSD3-263 30cm



2 mm

Cooley Atraumatic Dressing Forceps
RSD3-264 15cm
RSD3-265 20cm
RSD3-266 24cm
RSD3-267 30cm

1.2 mm

RSD3-268
18 cm
RSD3-269
21 cm

2 mm

RSD3-270
18 cm
RSD3-271
21 cm



Debakey Dressing Forceps



Debakey Forceps



RSD3-272 16 cm 1.5 mm
RSD3-273 20 cm 1.5 mm
RSD3-274 24 cm 1.5 mm



RSD3-275 16 cm 2 mm
RSD3-276 20 cm 2 mm
RSD3-277 24 cm 2 mm
RSD3-278 30 cm 2 mm




De Bakey Forceps




2 mm

RSD3-287
16 cm
Angular




2 mm

RSD3-288
20 cm
Angular



2 mm

RSD3-289
24 cm
Angular



2 mm

RSD3-290
30 cm
Angular

TWEEZERS

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



Feilichenfeld Forceps



RSD3-291
 7.5 cm
RSD3-292
 9 cm
RSD3-293
 11 cm
 Straight



Feilichenfeld



RSD3-294
 7.5 cm



RSD3-295
 9.5 cm



RSD3-296
 11.5 cm



Hunter Forceps



RSD3-297
 10.5 cm
 Straight



Hunter Forceps



RSD3-298
 10.5 cm
 Curved



Carmalt Forceps



RSD3-299
 10.5 cm
 Straight



Carmalt Forceps



RSD3-300
 10.5 cm
 Curved

TWEEZERS

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



RSD3-301
 15 cm
 Fig. 1
 Thin Pattern



RSD3-302
 15 cm
 Fig. 2
 Wide Pattern



RSD3-303
 16 cm



RSD3-304
 15 cm
 Fig. 1



RSD3-305
 15 cm
 Fig. 2



RSD3-306
 11.5 cm
 Fig. 1

Micro Surgery Instrument



RSD3-307
 18 cm
 Fig. 1



RSD3-308
 18 cm
 Fig. 2



RSD3-309
 18 cm
 Fig. 3



RSD3-310
 18 cm
 Fig. 4

SOLDERING PLIERS - TWEEZERS & SPLINTER FORCEPS

Пинцеты для удаления осколков

Pinces à souder - Pincas et pincas à échardes

Löt-Zange - Pinzette & SplinterZange

Пинцети для видалення осколків

Alicates Soldadura - Pinzas y pinzas de Splinter

لحام كماشة - ملاقط منشقة الملقط



RSD3-311
Soldering Pliers
16cm



RSD3-312
Soldering Pliers
16cm



RSD3-313
Posigrip Articulating
Paper Forceps
15cm



RSD3-314
Standard Articulating
Paper Forceps, Economy Grade
15cm



RSD3-315
Miller
15cm
Curved



Beer Cilia Forceps
RSD3-316
9 cm



Kuehne Cilia Forceps
RSD3-317
10.5 cm



Kuehne Cilia Forceps
RSD3-318
10.5 cm
angled sideways

SOLDERING PLIERS - TWEEZERS & SPLINTER FORCEPS

 Пинцеты для удаления осколков
 Pincettes à souder - Pincettes et pincettes à écharde
 Löt-Zange - Pinzette & SplinterZange

 Пінцети для видалення осколків
 Alicates Soldadura - Pinzas y pinzas de Splinter
 لحام كماشة - ملاقط منشقة الملقط



Splinter Forceps
 RSD3-319 8 cm
 RSD3-320 9 cm
 RSD3-321 10.5 cm
 RSD3-322 11.5 cm
 RSD3-323 12.5 cm
 RSD3-324 14.5 cm
 Straight serrated with
 guide pin



Splinter Forceps
 RSD3-325 8 cm
 RSD3-326 9 cm
 RSD3-327 10.5 cm
 RSD3-328 11.5 cm
 RSD3-329 12.5 cm
 RSD3-330 14.5 cm
 Curved serrated with
 guide pin



Splinter Forceps
 RSD3-331 8 cm
 RSD3-332 9 cm
 RSD3-333 10.5 cm
 RSD3-334 11.5 cm
 RSD3-335 12.5 cm
 RSD3-336 14.5 cm
 Straight plain with
 guide pin



Splinter Forceps
 RSD3-337 8 cm
 RSD3-338 9 cm
 RSD3-339 10.5 cm
 RSD3-340 11.5 cm
 RSD3-341 12.5 cm
 RSD3-342 14.5 cm
 Curved plain with
 guide pin



**Stieglitz Forceps
 with Magnifying
 Mirrors**
 RSD3-343
 12.5 cm

TWEEZERS

 Пинцеты
 Précelles - Bruxelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



Micro Adson



RSD3-350
12 cm



RSD3-351
15 cm



Adson



RSD3-352
12 cm



RSD3-353
15 cm



Adson



RSD3-354
12 cm



Micro Adson
Light Weight



RSD3-355
12 cm



RSD3-356
15 cm



Adson
Light Weight



RSD3-357
12 cm



RSD3-358
15 cm



Adson
Light Weight



RSD3-359
12 cm

TWEEZERS

 Пинцеты
 Précelles - Bruxelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



Micro Adson Tissue Forceps
1x2 teeth



RSD3-360
12 cm



RSD3-361
15 cm



Adson Tissue Forceps
1x2 teeth



RSD3-362
12 cm



RSD3-363
15 cm



Adson Brown Forceps
3x4 teeth



RSD3-364
12 cm



Micro Adson Light Weight
1x2 teeth



RSD3-365
12 cm



RSD3-366
15 cm



Adson Light Weight
1x2 teeth



RSD3-367
12 cm



RSD3-368
15 cm



Adson-Brown Light Weight
3x4 teeth



RSD3-369
12 cm

TWEEZERS

Пинцеты
 Précelles - Bruxelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



RSD3-370
14 cm

0,4 mm

RSD3-371
16 cm

T.C Standard Dressing Forceps



T.C Cushing Forceps
RSD3-372 18 cm



T.C Cushing Dressing Forceps
RSD3-373
18 cm Cross Serrated



Cushing Dressing T.C Forceps Angled
RSD3-375
18 cm



T.C Dissecting Tissue Forceps
RSD3-376 18 cm



T.C Resano Dressing Forceps
RSD3-377 25 cm



T.C Cushing Dressing Forceps
RSD3-378 19 cm

TWEEZERS

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



Wangensteen Dressing Forceps T.C.



RSD3-379
15.5 cm




RSD3-380
23 cm



RSD3-381
25.5 cm



T.C DeBakey Micro Needle
Tissue Forceps
RSD3-382 18 cm



T.C Bonney Dissecting
& Tissue Forceps



RSD3-383
1x2 Teeth
18 cm
Straight
X-Serrations



RSD3-384
1x2 Teeth
18 cm
Straight, Serrated



RSD3-385
2x3 Teeth
18 cm
Straight, Serrated



T.C POTTS-SMITH
Dressing Forceps
RSD3-386 15 cm
RSD3-387 18 cm
RSD3-388 20 cm
RSD3-389 23 cm
RSD3-390 25 cm



T.C POTTS-SMITH
Dressing Forceps
RSD3-391 16 cm
RSD3-392 18 cm
RSD3-393 21 cm
RSD3-394 24 cm
RSD3-395 26 cm



T.C Evans Pelvic Needle
Dressing Forceps
RSD3-396 25.5 cm



T.C Evans Dressing Forceps
RSD3-397 25 cm

TWEEZERS

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пинцети
 Pinzas para algodón
 ملاقط



T.C Mc Indoe Dressing Forceps
RSD3-398 15 cm



T.C De Bakey Dressing Forceps
RSD3-399 15 cm
RSD3-400 19 cm
RSD3-401 24 cm



T.C Gillies Tissue Forceps
RSD3-402 15 cm



RSD3-403
15 cm

Semken
with dissector end



RSD3-404
16 cm

Semken
1x2 teeth



RSD3-405
18 cm



T.C Oehler Tissue Forceps



RSD3-406
14 cm



RSD3-407
16 cm



RSD3-408
18 cm



RSD3-409
20 cm



RSD3-410
25 cm



RSD3-411
30 cm

TISSUE MEMBRANE FORCEPS



RSD6-216
Tissue Graft Forcep
Serrated Teeth
15cm



RSD6-217
Self Locking
Grooved
15cm



RSD6-218
Tissue Graft Forcep
Serrated Teeth
15cm



RSD6-219
Titanium Metal
Locking Plier



RSD6-220
TC Gripper Carbide
Tipped Tissue Forceps
Serrated Carbide Tips
14.5cm



RSD6-221
Adson Carbide
Tipped Tissue Forceps
Serrated Carbide Tips
12cm



RSD6-222
Chalazion Mucocoele
Forceps
9cm Length
12mm x 23mm
Inside Ring



RSD6-223
Micor Adson
Tissue Forcep 12cm
1x2 Teeth -
0.5mm Tips



RSD6-224
Adson 7x7
Tissue Forceps
7x7 Teeth
12cm



RSD6-225
Adson 7x7 Curved
Tissue Forceps
7x7 Teeth 12cm



RSD6-226
Suture Pliers
Suture Forceps
with
1.6mm Ring Tip
16cm



RSD6-227
Membrane Forcep
With Delicate
Ring Tip
10 cm



RSD6-228
Debakey Atraumatic
Tissue Forceps
Straight
Serrated Teeth
16cm



RSD6-229
Debakey Atraumatic
Tissue Forceps Curved
Serrated Teeth
16cm

TISSUE MEMBRANE FORCEPS



RSD6-230
Russian Long
Tissue Forceps
Serrated Tip
20cm



RSD6-231
Russian
Tissue Forceps
Serrated Tip
15cm



RSD6-232
Tissue Punch Forceps
To Remove Tissue
After Using Punch
15cm



RSD6-233
Narrow
13.5cm



RSD6-234
10.5cm

Membrane Tweezers Kuehne



RSD6-235
wide
13.5cm



RSD6-236
Bishop Harman
Micro Tissue Forcep
8.6cm 1x2 Teeth
- 0.3mm Tips



RSD6-237
Pin-Tweezers
length 15cm
with lock



RSD6-238
without lock
15cm



RSD6-239
Haltepinzette with lock
15.5cm, for pins, towards
the tips, etc



RSD6-240
Carry Forceps
15cm



RSD6-241
Carry Forceps
15cm




RSD6-242
Carry Forceps
15cm




SOLDERING PLIERS - TWEEZERS & SPLINTER FORCEPS

 Пинцеты для удаления осколков
 Pincettes à souder - Pincettes et pincettes à écharde
 Löt-Zange - Pinzette & SplinterZange

 Пінцети для видалення осколків
 Alicates Soldadura - Pinzas y pinzas de Splinter
 لحام كمانشة - ملاقط منشقة الملقط



Stieglitz Splinter Forceps
 RSD3-344 14 cm Str.
 RSD3-345 14 cm 45° Angled



Stieglitz Splinter Forceps
 RSD3-346
 14 cm 90° Angled



Peet Splinter Forceps
 RSD3-347
 11.5 cm




Ralk Splinter Forceps
 RSD3-348 15 cm Str.
 RSD3-349 15 cm Cvd.

EXPLORERS

Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين

Single Ended Probe Solid Handle Ø 6 mm, Length 150mm



Single Ended Probe Hollow Handle Ø 8 mm Length 150mm




EXPLORERS

Зонды
 Sondes
 Zahnsonden

Зонди
 Sondas
 المستكشفين

EXPLORERS INSTRUMENTS AVAILABLE IN
 Hollow Handle Ø 8 mm, Ø 10mm Solid Handle Ø 6mm
 Length 170MM



PROBE

Зонд
Sonde
Sonde

Зонд
Sonda
مستبر

Available in Hollow Handle
Ø 8 mm, Ø 10mm Solid Handle Ø 6mm




EXPLORERS

Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين


Explorers Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Pointed Probes Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Pointed Probes Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة

RSD3-490

Periodontal Set of 10 Pieces Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD3-492
EXP,23/17-A

RSD3-493
Exavators

RSD3-494
Filling
Instruments

RSD3-495
EXP,23/CP-18

RSD3-496
EXP,23/CP-QW

RSD3-497
EXP,23/CP-11

RSD3-498
EXP,23/CP-10

RSD3-499
EXP,23/CP-12

RSD3-500
EXP,23/CP-11.5-B

PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك للثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

RSD3-501
Columbia Set

RSD3-505
McCall Set

Columbia 2L/2R



RSD3-502

Columbia 4L/4R



RSD3-503

Columbia 13/14



RSD3-504

McCalls 14S/18S



RSD3-506

McCalls 17/18



RSD3-507

McCalls 13S/14S



RSD3-508

McCalls 13/14




RSD3-509

RSD3-510

William Probe Set Double Ended With Hollow Handle Ø 10 Mm, Length 170mm

Plasma Coated Handle



RSD3-511 EXP, 23/CO,QW	RSD3-512 EXP, 23/CP,10	RSD3-513 EXP, 23/CP,18	RSD3-514 EXP, 23/CP,12	RSD3-515 EXP, 23/CP,2	RSD3-516 EXP, 23/UNC,15	RSD3-517 EXP, 23/CP,11	RSD3-518 EXP, 23/CP, 12 ,S
----------------------------------	----------------------------------	----------------------------------	----------------------------------	---------------------------------	-----------------------------------	----------------------------------	--------------------------------------


EXPLORERS

Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين

Explorers Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

PLUGGERS



PERIODONTAL PROBES


Зонды пародонтологические
Sondes parodontales
sonde parodontali

зонди пародонтологічні
Sondas parodontales
المساطر اللثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SPREADERS

HEAT CARRIERS




PERIODONTAL PROBES

Зонды пародонтологические
Sondes parodontales
sonde parodontali

зонди пародонтологічні
Sondas parodontales
المساطر اللثة

Endodontics Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

ANTERIOR SINGLE ENDED




PERIODONTAL PROBES

Зонды пародонтологические
Sondes parodontales
sonde parodontali

зонди пародонтологічні
Sondas parodontales
المساطر اللثة

Endodontics Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
POSTERIOR SINGLE ENDED

RETRO FILLING
PLUGGERS



PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Инструменти пародонтологічні
instrumentos para periodoncia
المصكوك للثة

Heidemann Spatula

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Инструменти пародонтологічні
instrumentos para periodoncia
المسكوكات للثة

Excavators

Excavators Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD3-671



RSD3-672
Fig. 3



RSD3-673
Fig. 4



RSD3-674
Exc # 14



RSD3-675
Exc # 17



RSD3-676
Exc # 18



RSD3-677
Exc # 19



RSD3-678
Exc # 22


PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة

Composite Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Инструменти пародонтологічні
instrumentos para periodoncia
المصكوك النثة

Crane Kaplan

Available in Hollow Handle Ø 8 mm




RSD3-690



RSD3-691
Jaquette



RSD3-692
Columbia 13/14



RSD3-693
Kirkland 15/16



RSD3-694
Goldman Fox



RSD3-695
Gracey 11/12



RSD3-696
Gracey 13/14



RSD3-697
Columbia 2R/2L



RSD3-698
Orban Knife 1/2

Lucas Cures

Available in Hollow Handle Ø 8 mm



RSD3-699
Ø 2 mm



RSD3-700



RSD3-701
Ø 2.5 mm



RSD3-702
Ø 3 mm

PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك للثة

Sinus Lift Instruments Blue Plasma Coating

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD3-703
Fig. 1




RSD3-704
Fig. 2



RSD3-705
Fig. 3



RSD3-706
Fig. 4



RSD3-707
Fig. 5



RSD3-708
Fig. 6



RSD3-709
Fig. 7



RSD3-710
Fig. 8

Gingival Retractors

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD3-711



RSD3-712




RSD3-713

PERIODONTAL INSTRUMENTS

 Инструменты пародонтологические
 Instruments pour la parodontologie
 Instrumente für die Parodontologie

 Інструменти пародонтологічні
 instrumentos para periodoncia
 الصكوك اللثة

Blue Plasma Coating & Operative instruments Available in Hollow Handle Ø 8 mm



EXPLORERS

Зонды	Зонди
Sondes	Sondas
Zahnsonden	المستكشفين

PROBES



RSD3-822
Probe



RSD3-823
Probe



RSD3-824
Probe



RSD3-825
Probe



RSD3-826
Probe



RSD3-827
Probe



RSD3-828
Probe



RSD3-829
Probe

PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة



RSD3-830
Excavator




RSD3-831
Excavator



RSD3-832
Excavator



RSD3-833
Filling Instruments



RSD3-834
Filling Instruments



RSD3-835
Filling Instruments



RSD3-836
Filling Instruments



RSD3-837
Filling Instruments






RSD3-838
Cement Spatula

PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك للثة

MODELLING INSTRUMENTS



PERIODONTAL INSTRUMENTS

Инструменты пародонтологические

Instruments pour la parodontologie

Instrumente für die Parodontologie

Інструменти пародонтологічні


Instrumentos para periodoncia

المصكوك اللثة

Periodontal Instruments Set



Cotton Meriam Tweezer



Dental Mirror Magnify No.5 with Handle



Probe Fig. 9



Excavator Tip Size Ø2mm-1,5mm




Excavator Tip Size Ø1,2mm-1,5mm



Plugger Tip Size Ø1,5mm-2mm



Plugger Tip Size Ø2mm-2,5mm



Ironing Crescent Tip Size Ø3mm-3,2mm



Ironing Crescent Tip Size Ø2,5mm-2,5mm



Cement Spatula Tip Size 6mm



RSD3-847
Set of 10 Pieces.

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة


Periodontal Files

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

RSD3-848

Sugarman File Set

Excellent to remove and to crush strong calculus deposits



RSD3-849
9/10
Schluger
Curved File

RSD3-850
11/12
Buck
Straight File

RSD3-851
1S/2S
Sugarman File

RSD3-852
3S/4S
Sugarman File

RSD3-853
3-7
Hirschfeld

RSD3-854
5-11
Hirschfeld

RSD3-855
9-10
Hirschfeld

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة


Periodontal Files

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

RSD3-856

Sugarman File Set

Excellent to remove and to crush strong calculus deposits



RSD3-857
9/10
Schluger
Curved File

RSD3-858
11/12
Buck
Straight File

RSD3-859
1S/2S
Sugarman File

RSD3-860
3S/4S
Sugarman File

RSD3-861
3-7
Hirschfeld

RSD3-862
5-11
Hirschfeld


RSD3-863
9-10
Hirschfeld

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Инструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة

Sickle Scalars Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



PERIODONTAL INSTRUMENTS

 Инструменты пародонтологические
 Instruments pour la parodontologie
 Instrumente für die Parodontologie


 Інструменти пародонтологічні
 instrumentos para periodoncia
 المصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

JACQUETTE SCALERS

CRANE-KAPLAN



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SCALER CURETTE




Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Инструменти пародонтологічні
instrumentos para periodoncia
المصكوك للثة

SICKLE SCALERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SICKLE SCALERS



PERIODONTAL INSTRUMENTS

 Инструменты пародонтологические
 Instruments pour la parodontologie
 Instrumente für die Parodontologie


 Інструменти пародонтологічні
 instrumentos para periodoncia
 الصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

JACQUETTE SCALERS

CRANE-KAPLAN




Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SCALER CURETTE



PERIODONTAL INSTRUMENTS

 Инструменты пародонтологические
 Instruments pour la parodontologie
 Instrumente für die Parodontologie

 Інструменти пародонтологічні
 instrumentos para periodoncia
 الصكوك اللثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

PERIODONTAL KNIVES



RSD3-1142
1
Reddy Thin

RSD3-1143
2
Reddy Thin

RSD3-1144
3/4
Buck

RSD3-1145
5/6
Buck

RSD3-1146
1/2
Solt

RSD3-1147
1/2
Sanders

RSD3-1148
USC Towner
19/20

RSD3-1149
CK-3
Crane Kaplan

RSD3-1150
Dr. Allen
End-Cutting
Intrasulcular Knife

PERIODONTAL INSTRUMENTS

■ Инструменты пародонтологические
 ■ Instruments pour la parodontologie
 ■ Instrumente für die Parodontologie

■ Інструменти пародонтологічні
 ■ instrumentos para periodoncia
 ■ المصكوك اللثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

PERIODONTAL KNIVES



PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك للثة

PERIODONTAL CHISELS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

OCHSENBEIN

FEDI



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
instrumentos para periodoncia
المصكوك اللثة

PERIODONTAL CHISELS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD3-1176
1/2
Kramer-Nevins

RSD3-1177
36/37
Rhodes
Back-Action

RSD3-1178
44/45
Remington
Hoe

RSD3-1179
TG

RSD3-1180
TG-0
Modified TG with
Ochsendein cutouts

RSD3-1181
3/4
Solt

RSD3-1182
13K/TG

RSD3-1183
13K/13KL

FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




FILLING INSTRUMENTS

Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك

Amalgam Instruments and Burnishers Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




FILLING INSTRUMENTS

Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك

Amalgam Instruments and Burnishers Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Composite Gold Coated Instruments



RSD3-640
Markley 2



RSD3-641
Westcott N.21



RSD3-642
NYU/B2



RSD3-643
Filling Instruments



RSD3-644
Ball Filling Instruments



RSD3-645
Filling Instruments



RSD3-646
Filling Instruments



RSD3-647
Filling Instruments



RSD3-648
Plastic Filling Instruments



RSD3-649
Heidemann



RSD3-650
Plastic Filling Instruments

FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Composite Gold Tip Instruments



RSD3-651

Markley 2



RSD3-652

Westcott N.21



RSD3-653

NYU/B2



RSD3-654

Filling instruments



RSD3-655

Ball Filling Instruments



RSD3-656

Filling Instruments



RSD3-657

Filling Instruments



RSD3-658

Filling Instruments



RSD3-659

Plastic Filling Instruments



RSD3-660

Heidemann



RSD3-661

Plastic Filling Instruments

FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)


Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

CARVERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




FILLING INSTRUMENTS

Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك

CARVERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації


Instrumentos para obturaciones

ملء الصكوك

CARVERS INSTRUMENTS


Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

CLEOIDS



Nystrom

CHISELS



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації


Instrumentos para obturaciones

ملء الصكوك

ENAMEL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

MARGINTRIMMERS




FILLING INSTRUMENTS

Инструменты для реставрации
Instruments pour plombage
Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації
Instrumentos para obturaciones
ملء الصكوك

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
CHISELS



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації


Instrumentos para obturaciones

ملء الصكوك

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

Heidemann Spatulas

SPATULAS



FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

P.K. THOMAS WAXING INSTRUMENTS

P.K. THOMAS WAX CARVER SET OF 5 PCS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD3-1199 PKT-1 RSD3-1200 PKT-2 RSD3-1201 PKT-3 RSD3-1202 PKT-4 RSD3-1203 PKT-5

P.K. THOMAS WAX CARVER SET OF 5 PCS


Aluminum Anodizing Colours Handle Ø 8 mm



RSD3-1205 PKT-1 RSD3-1206 PKT-2 RSD3-1207 PKT-3 RSD3-1208 PKT-4 RSD3-1209 PKT-5



RSD3-1204
P.k. Thomas Wax Carvers Set of 5 Pcs



RSD3-1210
P.k. Thomas Wax Carvers Set of 5 Pcs
with Sterilized Cassette



RSD3-1210-A
P.k. Thomas Wax Carvers Set of 5 Pcs

FILLING INSTRUMENTS

Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC



RSD3-1211

RSD3-1212

RSD3-1213

RSD3-1214

RSD3-1215

RSD3-1216



RSD3-1217
Composite Instruments Set of 6 Pcs
Blue Titanium with Cassettes



RSD3-1217-A
Composite Instruments Set of 6 Pcs


EXCAVATORS

Экскаваторы
 Excavateurs
 Excavatoren
 Excavatori
 Excavadores
 حفارات

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm


Oval Spoon

Spoon



DENTINS REMOVALS


Spoon



EXCAVATORS


Экскаваторы
 Excavateurs
 Excavatoren
 حفارات
 Excavatori
 Excavadores

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Spoon Excavators Round Shape

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



IMPLANTOLOGY

Инструменты для имплантологии

Implantologie

Implantologie

Инструменти для имплантології

Implantologia

زراعة الأسنان

Curettes Hollow Handle 10mm Length 170mm



RSD6-207 Miller



RSD6-208 Miller




RSD6-209 LUCAS #84



RSD6-210 LUCAS #85



RSD6-211 LUCAS #86



RSD6-212 LUCAS #87



RSD6-213 LUCAS #88



RSD6-214 LUCAS86S
SERRATED




RSD6-215 LUCAS86AS
SERRATED

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

SICKLE SCALERS

Designed for the removal of supragingival calculus



RSD3-864
H6/H7
Hygienist

RSD3-865
H5/33
Hygienist-
Jacquette

RSD3-866
H5/U15
Hygienist-
Towner

RSD3-867
6/7
Goldman

RSD3-868
U15
Towner

RSD3-869
U15/30
Towner-
Jacquette


RSD3-870
U15/33
Towner-
Jacquette

CURETTES

Кюрети Кюрети
Curetas Curette
Küretten کرتس


SCALERS

JACQUETTE SCALERS



CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس



CURETTES

Кюрети	Кюрети
Curetas	Curette
Küretten	کرتس

SICKLE SCALERS

POSTERIOR

ANTERIOR



RSD3-890
135



RSD3-891
REM 3/4
Remington



RSD3-892
129



RSD3-893
T2/T3
Taylor



RSD3-894
137



RSD3-895
0/00
Morse

CURETTES

Кюрети	Кюрети
Curetas	Curette
Küretten	کرتس



RSD3-896
137
Mini Five



RSD3-897
OD1/33
Odont/
Jacquette



RSD3-898
128
Nebraska



RSD3-899
USC-128



RSD3-900
WH-2
Whiteside




RSD3-901
YG-15
Younger-
Good

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

BALL SCALERS

CRANE-KAPLAN



RSD3-902
CI-1

RSD3-903
CI-2/3

RSD3-904
152
Cumine

RSD3-905
B-1

RSD3-906
B-4

RSD3-907
B-2/3

RSD3-908
B-5/6

RSD3-909
CK-6

RSD3-910
CK-6D


RSD3-911
333/334
Mini Kaplan

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten كرتس

SCALER / CURETTE SETS

Mc CALL SET



RSD3-912
11/12
Scaler

RSD3-913
11A/12A
Scaler

RSD3-914
11/12S
Modified
McCall

RSD3-915
13S/14S
Pointed
Curette

RSD3-916
MTC 13/14
Pointed
Curette

RSD3-917
17/18
Curette

RSD3-918
17S/18PT
Curette

RSD3-919
19/20
Curette


RSD3-920
Goldman-Fox

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

SCALER / CURETTE SETS

GOLDMAN-FOX SET



CURETTES


Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

SCALER / CURETTE SETS

GOLDMAN-FOX SET

CATTONI SET

COLUMBIA




CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

UNIVERSAL CURETTES

BARNHART



CURETTES

■ Кюрети ■ Кюрети
■ Curetas ■ Curette
■ Kûretten ■ کرتس



RSD3-948 1/2 (shown) RSD3-949 1/2 Rigid



RSD3-950 3/4 (shown) RSD3-951 3/4 Rigid



RSD3-952 5/6 (shown) RSD3-953 5/6 Rigid RSD3-954 5/6 Extra Rigid



RSD3-955 7/8 (shown) RSD3-956 7/8 Rigid RSD3-957 7/8 Extra Rigid



GRACEY 1/2



GRACEY 7/8



GRACEY 11/12



GRACEY 13/14



GRACEY 15/16



GRACEY 17/18



RSD3-958 9/10 (shown) RSD3-959 9/10 Rigid



RSD3-960 11/12 (shown) RSD3-961 11/12 Rigid RSD3-962 11/12 Extra Rigid



RSD3-963 13/14 (shown) RSD3-964 13/14 Rigid RSD3-965 13/14 Extra Rigid



RSD3-966 15/16 (shown) RSD3-967 15/16 Rigid



RSD3-968 17/18 (shown) RSD3-969 17/18 Rigid


CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

GRACEY CURETTES

MICRO MINI

GRACEY PROPHY CURETTES



CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

RSD3-986-A

Set of 9 Pcs Gracey Curettes Solid Handle Ø 6 mm, Length 170mm




RSD3-986-B
Gracey Scaler Set of 7 Pcs

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کړتس

RSD3-995-A

Set of 9 Pcs Gracey Curettes Hollow Handle Ø 10 mm, Length 170mm




RSD3-995-B
Gracey Scaler Set of 7 Pcs



CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس




RSD3-1003
 Gracey Curette Set of 7 Pcs
 Gracey Curette 1-2
 Gracey Curette 3-4
 Gracey Curette 5-6
 Gracey Curette 7-8
 Gracey Curette 9-10
 Gracey Curette 11-12
 Gracey Curette 13-14



RSD3-1011
 Gracey Curette
 (Black & Gold Coating)
 Set of 7 Pcs
 Gracey Curette 1-2
 Gracey Curette 3-4
 Gracey Curette 5-6
 Gracey Curette 7-8
 Gracey Curette 9-10
 Gracey Curette 11-12
 Gracey Curette 13-14

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس



RSD3-1019
Gracey Curette Set of 7 Pcs
Gracey Curette 1-2
Gracey Curette 3-4
Gracey Curette 5-6
Gracey Curette 7-8
Gracey Curette 9-10
Gracey Curette 11-12
Gracey Curette 13-14




RSD3-1027
Gracey Curette Set of 7 Pcs
with Silicon Handle
Gracey Curette 1-2
Gracey Curette 3-4
Gracey Curette 5-6
Gracey Curette 7-8
Gracey Curette 9-10
Gracey Curette 11-12
Gracey Curette 13-14

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

RSD3-1035

SET OF 7 PCS HOLLOW HANDLE PLASMA BLUE COATED Ø 10 MM, LENGTH 170MM



RSD3-1035-A
Gracey Curettes Set of 7 Pcs
 Gracey Curettes 1/2
 Gracey Curettes 3/4
 Gracey Curettes 5/6
 Gracey Curettes 7/8
 Gracey Curettes 9/10
 Gracey Curettes 11/12
 Gracey Curettes 13/14



CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

HOLLOW HANDLE Ø 10 MM, LENGTH 170MM WITH SILICON RING




RSD3-1043

Gracey Curettes Set of 7 Pcs

Gracey Curettes 1/2
 Gracey Curettes 3/4
 Gracey Curettes 5/6
 Gracey Curettes 7/8
 Gracey Curettes 9/10
 Gracey Curettes 11/12
 Gracey Curettes 13/14

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس



RSD3-1051
Gracey Curettes Set of 7 Pcs
with Silicon Grip
 Gracey Curettes 1/2
 Gracey Curettes 3/4
 Gracey Curettes 5/6
 Gracey Curettes 7/8
 Gracey Curettes 9/10
 Gracey Curettes 11/12
 Gracey Curettes 13/14

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

Silicone Handles



RSD3-1052
Explorer 1



RSD3-1053
Explorer 2



RSD3-1054
Explorer 3



RSD3-1055
Root Canal Spreader 1



RSD3-1056
Root Canal Spreader 2



RSD3-1057
Root Canal Spreader 3



RSD3-1058
Root Canal Spreader 4



RSD3-1059
Amalgam
Instruments 1



RSD3-1060
Amalgam
Instruments 2



RSD3-1061
Amalgam
Instruments 3



RSD3-1062
Amalgam
Instruments 4



RSD3-1063
Amalgam
Instruments 5



RSD3-1064
Amalgam
Instruments 6

CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس



RSD3-1065
Filling Instruments 1



RSD3-1066
Filling Instruments 2



RSD3-1067
Filling Instruments 3



RSD3-1068
Filling Instruments 4



RSD3-1069
Filling Instruments 5



RSD3-1070
Lucas Curettes 1



RSD3-1071
Lucas Curettes 2



RSD3-1072
Lucas Curettes 3



RSD3-1073
Lucas Curettes 4



RSD3-1074
Lucas Curettes 5

CURETTES


Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

RSD3-1075-A
Set of 6 Pcs

GRACEY MINI FIVE
HOLLOW HANDLE Ø 10 MM, LENGTH 170MM
CERAMIC COATED HANDLE

RSD3-1081-A
Set of 6 Pcs

GRACEY AFTER FIVE
HOLLOW HANDLE Ø 10 MM, LENGTH 170MM
CERAMIC COATED HANDLE




RSD3-1075 Gracey 1/2	RSD3-1076 Gracey 3/4	RSD3-1077 Gracey 5/6	RSD3-1078 Gracey 7/8	RSD3-1079 Gracey 11/12	RSD3-1080 Gracey 13/14	RSD3-1081 Gracey 1/2 AF	RSD3-1082 Gracey 3/4 AF	RSD3-1083 Gracey 5/6 AF	RSD3-1084 Gracey 7/8 AF	RSD3-1085 Gracey 11/12 AF	RSD3-1086 Gracey 13/14 AF
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-------------------------------------	-------------------------------------	--------------------------------------	--------------------------------------	--------------------------------------	--------------------------------------	--	--

SCALERS

 Стоматологический набор скейлеров
 Instruments a detartrer
 Zahnreiniger

 Стоматологічний набір скейлерів
 Limpiadores de sarro
 المتسلقون

SCALERS Solid Handle Ø 6 mm



RSD3-1127

RSD3-1128

RSD3-1129


RSD3-1130

RSD3-1131

RSD3-1132

RSD3-1133

Tartar Remover Instruments Solid Handle Ø 6 mm




CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

BONE CURETTES

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



CURETTES

Кюрети Кюрети
 Curetas Curette
 Küretten کرتس

LUCAS SPOON SHAPE

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

LUCAS



RSD3-1190
84

RSD3-1191
85

RSD3-1192
86

RSD3-1193
87

RSD3-1194
88

RSD3-1195
85/86
Lucas
Surgical Curette

RSD3-1196
86/87
Lucas
Surgical Curette

RSD3-1197
Lucas

RSD3-1198
Lucas

BONE FILE

Инструменты для костной ткани
Limes a os
Knochenfeilen

Инструменти для кісткової тканини
Limas para huesos
ملفات العظام



RSD3-1166
Fig. 1
PoloKoff



RSD3-1167
Fig. 33



RSD3-1168
Fig. 45
Miller



RSD3-1169
Fig. 64
Miller



RSD3-1170
Fig. 10



RSD3-1171
Fig. 1
Wahl



RSD3-1172
Fig. 2
Wahl



RSD3-1173
Fig. 12
Howard



RSD3-1174
Fig. 12A



RSD3-1175
Fig. 12CA

CEMENT SPATULAS

Шпатели для цемента
Spatules à ciment
Zementspatel

Шпателі для цементу
Espátulas para cemento
ملاعق الأسمنت

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




RSD3-1218
Fig. 1




RSD3-1219
Fig. 2



RSD3-1220
Fig. 3



RSD3-1221
Fig. 00



RSD3-1222
Fig. 0



RSD3-1223
Fig. 1



RSD3-1224
Fig. 2




RSD3-1225
Fig. 3



RSD3-1226
Fig. 4



RSD3-1227
Weston



RSD3-1228



RSD3-1229

CEMENT SPATULAS

Шпатели для цемента
Spatules à ciment
Zementspatel

Шпателі для цементу
Espátulas para cemento
ملاعق الأسمنت


Wax and Modelling Carvers



CARVERS INSTRUMENTS

Инструменты для реставрации
Instruments a modeler
Modellierinstrumente

Инструменти для реставрації
Instrumentos para modelar
أدوات تقطيع



RSD3-1246-A
Carvers Set of 10 Pcs
01 Wax Knife
02 Cement Spatulas
03 Cement Spatulas
04 Cement Spatulas
05 Cement Spatulas
06 Cement Spatulas
07 Cement Spatulas
08 Cement Spatulas
09 Cement Spatulas
10 Cement Spatulas



RSD3-1238
Graded Handle
Hylin Carver



RSD3-1239
Hylin



RSD3-1240
Hylin-Le Cron



RSD3-1241
Hylin-Ceramica



RSD3-1242
Hylin



RSD3-1243



RSD3-1244
Le Cron



RSD3-1245
Le Cron



RSD3-1246
Le Cron

WAX KNIVES

Ножи для воска
Couteaux à cire
Wachsmesser

Ножі для воску
Cuchillos para cera
السكاكين الشمع

WAX KNIVES



RSD3-1247
Fig. 1
Fahnenstock
125mm

RSD3-1248
Fig. 2
Fahnenstock
170mm

RSD3-1249
Fig. 1
Lessmann
125mm

RSD3-1250
Fig. 2
Lessmann
170mm

RSD3-1251
Fig. 1
with small spoon
125mm

RSD3-1252
Fig. 2
with large spoon
170mm

WAX KNIVES

Ножи для воска
Couteaux à cire
Wachsmesser

Ножі для воску
Cuchillos para cera
السكاكين الشمع

PLASTER KNIVES



RSD3-1253
160mm



RSD3-1254
160mm



RSD3-1255
Buffalo



RSD3-1256
Buffalo



RSD3-1257
Plaster Knife



RSD3-1258
170mm

WAX KNIVES

Ножи для воска
Couteaux à cire
Wachsmesser

Ножі для воску
Cuchillos para cera
السكاكين الشمع

WAX KNIVES



RSD3-1259
Sachs



RSD3-1260
130mm
Ermert



RSD3-1261
Hook Knives



RSD3-1262
125mm
Antkoviak



RSD3-1263



Special Inox
RSD3-1264 170mm
RSD3-1265 200mm

ALGINATE & PLASTER SPATULAS

Шпатель для замешивания
Spatules a alginate et a platre
Alginat- und Gipsspatel

Шпатель для змішування
Espatulas para alginato y para yeso
ملاعق الجص الجينات

WAX SPATULAS



RSD3-1266
Fig. 1
18.5 cm

RSD3-1267
Fig. 2
18 cm

RSD3-1268
Wooden Handle
19.5cm curved

RSD3-1269
Wooden Handle
20cm curved

RSD3-1270
Wooden Handle
22.5cm curved

RSD3-1271
Wooden Handle
19.5cm

RSD3-1272
Wooden Handle
21.5cm

RSD3-1273
Wooden Handle
20.5cm





Maxiofacial Surgery Instruments

A Bard-Parker handle with a No. 15 blade is the usual scalpel used to make incisions on the gums where a flap can be raised to gain access to tooth.

Retractors

Tongue, cheek and flap retractors are needed to provide adequate access to the tooth and are legion in number.

Needle Holders and Sutures

Sutures are used to replace and secure the gum flap in its former position, to limit the size of the wound and therefore restrict subsequent contamination with food debris, and to assist in controlling bleeding. The usual suture size for intra-oral work is 4/0 and the material may be absorbable or non-absorbable Scissors Dissecting scissors and suture-cutting scissors are usually used.

Artery Forceps

These are designed to hold soft tissue without damaging it. College tweezers are used to lift up sutures prior to removing them.

Periosteal Elevators

A number are available and are similarly used like a normal elevator but periosteal elevators have two heads – one at each end. One end is sharp-tipped to expand the tooth socket while the other has a blunt end to remove granulation tissues or raise a gum flap.

Chisels and Mallet

Chisels are more appropriate in the operating theater and are particularly useful for disto-angular third molars and upper wisdom teeth extraction.

Bone Curettes

The Mitchells trimmer is probably the most valuable instrument in the category and is used to remove granulation tissue in the tooth socket.

Bone Ronguer and Bone File


Used to smoothen any sharp bone edges after tooth extraction is done.

SCALPEL HANDLES

Ручки для скальпелей
 Manches de bistouris
 Skalpellgriffe

Ручки для скальпелів
 Mangos de bisturi
 مشرط مقابض

Scalpel Handles



SCALPEL BLADE & SCALPEL KNIVES

Ножи для скальпеля и скальпеля

Blapel Blade & Scalpel Knives

Scalpel Blade & Scalpel Messer

Ножи скальпельного леза та скальпеля

Cuchillos de Escalpo y Escalpo

مشرط شفرة و مشرط السكاكين

STERIL SCALPEL BLADES



RSD4-011
10



RSD4-012
11



RSD4-013
20



RSD4-014
21



RSD4-015
22



RSD4-016
12



RSD4-017
15



RSD4-018
23



RSD4-019
24



RSD4-020
16.5cm



Kit of scalpel blade
remover and orientable
blade scalpel handle

SCALPEL KNIVES



DIEFFENBACH
12.5cm



RSD4-021
Fig. 1
12.5cm



RSD4-022
Fig. 2
13cm



RSD4-023
Fig. 4
13.5cm



RSD4-024
Fig. 5
14cm



RSD4-025
Fig. 7
12.5cm



RSD4-026
Fig. 8
13cm



RSD4-027
Fig. 9
13.5cm



RSD4-028
Fig. 10
14cm

SCALPEL KNIVES

Скальпельные ножи Ножі скальпеля
 Couteaux à scalpel Cuchillos Scalpel
 Skalpellmesser سكاكين مشرط



Dieffenbach
16cm



RSD4-029
Fig. 1
16cm



RSD4-030
Fig. 2
16cm



RSD4-031
Fig. 3
16cm



RSD4-032
Fig. 4
16cm



RSD4-033
Fig. 5
16cm



RSD4-034
Fig. 6
16cm



RSD4-035
Fig. 7
16cm



RSD4-036
Fig. 8
16cm



RSD4-037
Fig. 9
16cm



RSD4-038
Fig. 10
16cm



RSD4-039
Fig. 15
16cm



RSD4-040
Fig. 16
16cm

RETRACTORS - MOUTH GAGS

Роторасширители

Ecarteurs - Ouvre-bouches

Vestibulumabhalter / Mundspanner

Роторозширювачі

Separadores - Abre-bocas

الكاشطات - الكعصات الفم



RSD4-041 10.5cm
RSD4-042 13.5cm
Weitlaner
Fig. 1



RSD4-043 10.5cm
RSD4-044 13.5cm
Weitlaner
Fig. 2



RSD4-045
Stainless Steel
Double
Ended Cheek
Retractors



RSD4-046
Stainless Steel
Single
Ended Cheek
Retractors



RSD4-047 11cm
RSD4-048 13cm
RSD4-049 15cm
Weitlander



RSD4-050 11cm
RSD4-051 13cm
RSD4-052 15cm
Jennings



RSD4-053
Denhard
12cm




RSD4-054
Roser-Konig
16cm

RSD4-055
Roser-Konig
19cm

RETRACTORS - MOUTH GAGS

- | | |
|--|--|
|  Роторасширители |  Роторозширювачі |
|  Ecarteurs - Ouvre-bouches |  Separadores - Abre-bocas |
|  Vestibulumabhalter / Mundspanner |  الكاشطات - الكعصات الفم |




RSD4-057
Mouth Gag Molt
11cm for (Children)

RSD4-058
Mouth Gag Molt
14cm for (Adults)

RETRACTORS

Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
الكامشات



Heiss Retractor
RSD4-106 10 cm
Sharp 4:4
RSD4-107 10 cm
Blunt 4:4

Finsen Retractor
RSD4-108 7 cm
3 x 4 teeth, sharp
RSD4-109 7 cm
3 x 4 teeth, Blunt

Finsen Retractor
RSD4-110 5 cm
2 x 3 teeth sharp/blunt
RSD4-111 7 cm
2 x 3 teeth sharp/blunt

Buckley Retractor
RSD4-112
6 cm



ALM Retractor
RSD4-113 7 cm
4 x 4 teeth sharp/blunt
RSD4-114 10 cm
4 x 4 teeth sharp/blunt

Jansen Retractor
RSD4-115 10 cm
Sharp 3x3
RSD4-116 10 cm
Blunt 3x3

**Markham-Meyerdig
Hemilaminectomy Retractor**
RSD4-117 18cm
Side Blade 2.5cm wide

Gelpi-Loktite Retractor

RSD4-118
9 cm
6 mm Sharp
RSD4-119
11 cm
10 mm Sharp
RSD4-120
14 cm
10 mm Sharp
RSD4-121
16 cm
10 mm Sharp
RSD4-122
19cm
10 mm Sharp


RETRACTORS

Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
الكامشات



Mollison Retractor
RSD4-123 15 cm
sharp



Milligan Retractor
RSD4-124 13 cm
18.5 mm x 20 mm
Blunt, 3x3



Johnson Retractor
RSD4-125 12.5cm



Miskimon Retractor
RSD4-126 22 cm
33 mm x 38 mm
Sharp, 7x7



Travers Retractor
RSD4-127
21 cm Blunt



Righby Retractor
RSD4-128 35 x 20 mm
RSD4-129 50 x 20 mm
RSD4-130 65 x 15 mm
18 cm



Cone Retractor

RSD4-131 31 cm
3x4 Blunt prongs,
13 mm x 14.5 mm
RSD4-132 16.5 cm
17 mm x 19 mm
3x4 Sharp prongs,
RSD4-133 16.5 cm
17 mm x 19 mm
3x4 Blunt prongs,
RSD4-134 20 cm
21 mm x 210 mm
3x4 Blunt prongs,

RETRACTORS - MOUTH GAGS

Роторасширители
 Ecarteurs - Ouvre-bouches
 Vestibulumabhalter / Mundspanner

Роторозширювачі
 Separadores - Abre-bocas
 الكاشطات - الكعصات الفم



RSD4-059
14cm

RSD4-060
15cm

Cheek Retractor



RSD4-061
Wangen Holder
15cm, with or without
lip protection



RSD4-062
Wangen Holder
15cm, with
or without lip protection



RSD4-063
Wangen Holder
Wangen holder with dual
function, 17cm



RSD4-064
Retractor with
Ergo Grip
24cm
ergonomic handle



RSD4-065
Cheek Retractor
Retractor Ivory
Width adjustable
& flexible, 20cm

RETRACTORS

 Ранорасширители
 Ecarteurs
 Wundhaken

 Ранорозширювачі
 Separadores
 الكاشحات


RETRACTORS



RSD4-066
Columbia
14.5 cm



RSD4-067
Sternberg
14 cm
Double Ended



RSD4-068
Sternberg
16 cm



RSD4-069
Sternberg
14 cm
Single Ended



RSD4-070
Oringer Medium
11x8 cm



RSD4-071
Brunings
Cheek Retractor
19cm



RSD4-072
Minnesota
Small
14cm



RSD4-073
Minnesota
Large
16cm



RSD4-074
Cheek Retractor
20cm
Easy to Bend
by Hand



RSD4-075
Cheek Retractor
and Tongue
Pushers Bishop with
Comfortable Margin,
16cm



RSD4-076
Cheek Retractor
Combination
25cm

RETRACTORS

Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
الكامشات



RSD4-077
Cheek Retractor
Long



RSD4-078
Cheek Retractor
Long



RSD4-079
Cheek Retractor
Long



RSD4-080
Unique Retractor
Length 17.5cm
angled 45°



RSD4-081
Unique Retractor
Length 17.5cm
angled 90°



RSD4-082
Retractor
16cm



RSD4-083
Retractor
19cm



RSD4-084
Mini-Langenbeck
Retractor
16cm, Hooks 22 mm



RSD4-085
Cheek Retractor
Combination
21cm, 2 Pcs



RSD4-086
Senn Miller Retractor
16cm Hooks 22 mm



RSD4-087
75x75
Austin Tissue
Retractors

8cm
Available with
round tip on
request

RSD4-088
75x75
Austin Tissue
Retractors




Joll (Friedman)
RSD4-089 15.5 cm
Thyroid- and Vaginal Retractor

RETRACTORS

Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
الكامشات



GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren


Ножиці для ясен

Tijeras para encias

مقص صمغ



RSD4-135
Micro Iris
9 cm
Straight



RSD4-136
Micro Iris
9 cm
Curved



RSD4-137
Iris
11.5 cm
Curved



Iris
RSD4-138 9cm
RSD4-139 10.5cm
RSD4-140 11.5cm
RSD4-141 12cm
Straight



Iris
RSD4-142 9cm
RSD4-143 10.5cm
RSD4-144 11.5cm
RSD4-145 12cm
Curved



RSD4-146
Angular, Saw edge
11.5 cm



RSD4-147
Angular
11.5 cm

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias

مقص صمغ



RSD4-148
Goldman-Fox
13 cm
Saw edge,
straight



RSD4-149
Goldman-Fox
13 cm
Saw edge,
curved



RSD4-150
Goldman-Fox
13 cm
Saw edge,
double curved



RSD4-151
La Grange
11.5 cm
Saw edge,
double curved



RSD4-152
Neumann
12 cm
Saw edge,
double curved



RSD4-153
Kilner (Ragnell)
12 cm
Straight



RSD4-154
Kilner (Ragnell)
12 cm
Curved



RSD4-155
Joseph
14 cm
Straight

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias

مقص صمغ



Joseph
RSD4-156 14 cm
Curved



Kelly
RSD4-157 16cm
RSD4-158 18cm
Straight



Kelly
RSD4-159 16cm
RSD4-160 18cm
Curved



Locklin
RSD4-161 16 cm
Saw edge,
angular



Locklin
RSD4-162 12 cm
RSD4-163 16 cm
Saw edge,
angular



RSD4-164
Dean
17 cm
Saw edge,
angular



RSD4-165
Straight
12 cm



RSD4-166 11.5cm
RSD4-167 13cm
RSD4-168 14cm
Straight



RSD4-169 11.5cm
RSD4-170 13cm
RSD4-171 14cm
Curved



Mayo
RSD4-172 14 cm
RSD4-173 16 cm
RSD4-174 18 cm
Straight
RSD4-175 14 cm
RSD4-176 16 cm
RSD4-177 18 cm
Curved

GUM SCISSORS

Ножницы десневые


Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias

مقص صمغ



RSD4-178
Universal
12 cm
Saw edge



Beebee
RSD4-179 10cm
RSD4-180 12cm
Saw edge, straight



Beebee
RSD4-181 10cm
RSD4-182 12cm
Saw edge, curved



RSD4-183
15 cm



RSD4-184
Dental scissors
large 17.5 cm



RSD4-185 14cm
RSD4-186 18cm
universal
autoclavable up to
143 °C

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножици для ясен

Tijeras para encias

مقص صمغ



Standard Operating Scissors

Sharp/Sharp Straight

RSD4-187	11.5 cm
RSD4-188	13.0 cm
RSD4-189	14.0 cm
RSD4-190	15.0 cm
RSD4-191	16.0 cm
RSD4-192	17.5 cm
RSD4-193	18.5 cm
RSD4-194	20.0 cm



Standard Operating Scissors

Sharp/Sharp Curved

RSD4-195	11.5 cm
RSD4-196	13.0 cm
RSD4-197	14.0 cm
RSD4-198	15.0 cm
RSD4-199	16.0 cm
RSD4-200	17.5 cm
RSD4-201	18.5 cm
RSD4-202	20.0 cm



Standard Operating Scissors

Blunt/Blunt Straight

RSD4-203	11.5 cm
RSD4-204	13.0 cm
RSD4-205	14.0 cm
RSD4-206	15.0 cm
RSD4-207	16.0 cm
RSD4-208	17.5 cm
RSD4-209	18.5 cm
RSD4-210	20.0 cm



Standard Operating Scissors

Blunt/Blunt Curved

RSD4-211	11.5 cm
RSD4-212	13.0 cm
RSD4-213	14.0 cm
RSD4-214	15.0 cm
RSD4-215	16.0 cm
RSD4-216	17.5 cm
RSD4-217	18.5 cm
RSD4-218	20.0 cm



Standard Operating Scissors

Sharp/Blunt Straight

RSD4-219	11.5 cm
RSD4-220	13.0 cm
RSD4-221	14.0 cm
RSD4-222	15.0 cm
RSD4-223	16.0 cm
RSD4-224	17.5 cm
RSD4-225	18.5 cm
RSD4-226	20.0 cm



Standard Operating Scissors

Sharp/Blunt Curved

RSD4-227	11.5 cm
RSD4-228	13.0 cm
RSD4-229	14.0 cm
RSD4-230	15.0 cm
RSD4-231	16.0 cm
RSD4-232	17.5 cm
RSD4-233	18.5 cm
RSD4-234	20.0 cm

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias

مقصص صمغ



Mayo

RSD4-235 14 cm RSD4-238

RSD4-236 16 cm RSD4-239

RSD4-237 18 cm RSD4-240

Straight

Curved



Metzenbaum

RSD4-241 14 cm RSD4-245

RSD4-242 16 cm RSD4-246

RSD4-243 18 cm RSD4-247

RSD4-244 20 cm RSD4-248

Straight

Curved



RSD4-249 9cm

RSD4-250 11cm

RSD4-251 13cm

RSD4-252 14cm



Spencer



Lister Bandage
Plaster Scissors

RSD4-253 14cm

RSD4-254 16cm

RSD4-255 18cm

Angled, Probe End



Fergusson Operating
Dissecting Scissors

RSD4-256 16 cm

RSD4-257 18 cm

B/B, Angled

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножици для ясен

Tijeras para encias

مقص صمغ



Wertheim Operating Dissecting Scissors

RSD4-258 14.5 cm
RSD4-259 22.5 cm
B/B, Cvd.



Metzenbaum Dissecting Scissors

RSD4-260 20 cm
RSD4-261 23 cm
RSD4-262 28 cm
B/B, Cvd.



Metzenbaum-Jamison Scissors

RSD4-263 15 cm
S/S, Str.



Metzenbaum-Baby Operating Dissecting Scissors

RSD4-264
15 cm
B/B Str.

RSD4-265
15 cm
B/B Cvd.



Mayo Dissecting & Gynecological Scissors
T/C Insert



RSD4-266
14.5 cm
RSD4-267
15 cm
RSD4-268
17 cm
RSD4-269
23 cm
B/B, Str.



RSD4-270
14.5 cm
RSD4-271
15 cm
RSD4-272
17 cm
RSD4-273
23 cm
B/B, Cvd.



Mayo-stille Dissecting & Gynecological Scissors
T/C Insert



RSD4-274
15 cm
RSD4-275
17 cm
RSD4-276
19.5 cm
RSD4-277
21 cm
B/B, Str.



RSD4-278
15 cm
RSD4-279
17 cm
RSD4-280
19.5 cm
RSD4-281
21 cm
B/B, Cvd.

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножици для ясен

Tijeras para encias

مقص صمغ



RSD4-282
22 cm
B/B, Str.

RSD4-283
22 cm
B/B, Cvd.



RSD4-284
17 cm
B/B, Str.

RSD4-285
17 cm
B/B, Cvd.



Metzenbaum-Fino
Operating Scissors
T/C Insert



RSD4-286 14.5 cm
RSD4-287 18 cm
RSD4-288 20 cm
RSD4-289 23 cm
RSD4-290 25 cm
RSD4-291 28 cm
B/B Str.



RSD4-292 14.5 cm
RSD4-293 18 cm
RSD4-294 20 cm
RSD4-295 23 cm
RSD4-296 25 cm
RSD4-297 28 cm
B/B Cvd.



Metzenbaum-Nelson
Operating Scissors
T/C Insert



RSD4-298 18 cm
RSD4-299 20 cm
RSD4-300 23 cm
RSD4-301 25 cm
RSD4-302 28 cm
RSD4-303 30 cm
B/B Str.



RSD4-304 18 cm
RSD4-305 20 cm
RSD4-306 23 cm
RSD4-307 25 cm
RSD4-308 28 cm
RSD4-309 30 cm
S/B Str.



RSD4-310 18 cm
RSD4-311 20 cm
RSD4-312 23 cm
RSD4-313 25 cm
RSD4-314 28 cm
RSD4-315 30 cm
S/S Str.



Metzenbaum-Nelson
Operating Scissors
T/C Insert



RSD4-316 18 cm
RSD4-317 20 cm
RSD4-318 23 cm
RSD4-319 25 cm
RSD4-320 28 cm
RSD4-321 30 cm
B/B Cvd.



RSD4-322 18 cm
RSD4-323 20 cm
RSD4-324 23 cm
RSD4-325 25 cm
RSD4-326 28 cm
RSD4-327 30 cm
S/B Cvd.



RSD4-328 18 cm
RSD4-329 20 cm
RSD4-330 23 cm
RSD4-331 25 cm
RSD4-332 28 cm
RSD4-333 30 cm
S/S Cvd.

GUM SCISSORS

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



**Metzenbaum-Fino
Operating Scissors T/C**

RSD4-334
14.5 cm
B/B Str.

RSD4-335
14.5 cm
S/S Str.

RSD4-336
14.5 cm
B/B Cvd.

RSD4-337
14.5 cm
S/S Cvd.



**Metzenbaum-Fino
Scissors T/C**

RSD4-338
18 cm
RSD4-339
20 cm
RSD4-340
23 cm
B/B Str.

RSD4-341
18 cm
RSD4-342
20 cm
RSD4-343
23 cm
B/B Cvd.

RSD4-344
18 cm
RSD4-345
20 cm
RSD4-346
23 cm
B/B Slightly
Cvd.



**Kaye Face Lift Scissors
T/C Inserted**
RSD4-347 11.5cm
RSD4-348 14 cm
RSD4-349 15 cm
RSD4-350 18 cm
RSD4-351 20 cm
B/B, Str.

Kaye Face Lift Scissors
RSD4-352 11.5 cm
RSD4-353 14 cm
RSD4-354 15 cm
RSD4-355 18 cm
RSD4-356 20 cm
B/B, Cvd.



Gorney Scissors T/C Inserted
RSD4-357 18 cm
RSD4-358 18 cm serrated
RSD4-359 23 cm serrated
B/B, Str.

RSD4-360 18 cm
RSD4-361 23 cm
B/B, Cvd, Serrated

RSD4-362 18 cm
B/B, Cvd.

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischschere

Ножици для ясен

Tijeras para encias

مقص صمغ

T.C GOLD



RSD4-363
Iris
11 cm
Straight

RSD4-364
Iris
11 cm
Curved

RSD4-365
Goldman-Fox
13 cm
Saw edge, straight

RSD4-366
Goldman-Fox
13 cm
Saw edge, Curved

RSD4-367
Iris Double Bend
11 cm



RSD4-368
La Grange
11.5 cm
Saw edge,
double Curved

Mayo
RSD4-369 14 cm
RSD4-370 17 cm
RSD4-371 23 cm
Straight

Mayo
RSD4-372 14 cm
RSD4-373 17 cm
RSD4-374 23 cm
Curved

Dean Scissors TC
RSD4-375
17cm curve

Metzenbaum
RSD4-376 14.5 cm
RSD4-377 18.0 cm
RSD4-378 20.0 cm
Straight

Metzenbaum
RSD4-379 14.5 cm
RSD4-380 18.0 cm
RSD4-381 20.0 cm
Curved

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias

مقص صمغ

T.C GOLD




Kelly
RSD4-382
16 cm
Straight



Kelly
RSD4-383
16 cm
Curved



RSD4-385
Beebee
11 cm
Saw edge,
straight



RSD4-386
Beebee
11 cm
Saw edge,
curved



RSD4-387
Beebee
12 cm
Saw edge,
straight



RSD4-388
Beebee
12 cm
Saw edge,
curved



RSD4-389
Universal
12 cm



RSD4-390
Super Cut Goldman Fox
13.2cm Curved
Serrated
For Slip-Free
Cutting



RSD4-391
Super Cut Iris
11cm Straight
Serrated For
Slip-Free Cutting



RSD4-392
Super Cut Iris
11cm Curved
Serrated For
Slip-Free Cutting



RSD4-393
Super Cut LaGrange
11.5cm
Double Curved
Serrated For
Slip-Free Cutting

Double
Curved

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias


مقص صمغ




Potts Tenotomy Scissors
RSD4-394 15 cm
RSD4-395 17 cm
S/S



Potts-Smith Vascular Scissors
S/S Angled




RSD4-396
19 cm



RSD4-397
19 cm



RSD4-398
19 cm



RSD4-399
19 cm



Potts-Smith Scissors
RSD4-400 19.5 cm
B/B Angled.



Cardiovascular Scissors
RSD4-401 12 cm
Probe End



Enterotomy Scissors
RSD4-402 21 cm



Enterotomy Scissors
RSD4-403 21 cm



Enterotomy Scissors
RSD4-404 21 cm



Spencer Stitch Scissors
RSD4-405 9 cm
RSD4-406 11 cm
RSD4-407 14 cm
Str.

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias

مقص صمغ



Artery Scissors
Probe end Angled



RSD4-408
18 cm



RSD4-409
18 cm




RSD4-410
18 cm




RSD4-411
18 cm



RSD4-412
18 cm




RSD4-413
18 cm



RSD4-414
18 cm



RSD4-415
18 cm



RSD4-416
18 cm



RSD4-417
18 cm

Quality Symbol

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias


مقص صمغ




Martin-Cartilage Scissors
RSD4-418 20 cm
S/S Serrated



Mc Indoe Scissors
RSD4-419 19 cm
S/S Str.



Potts-Smith Scissors
RSD4-420 19 cm
B/B Cvd.



Craford Scissors
RSD4-421 30 cm
B/B Cvd.



RSD4-422
10 cm
RSD4-423
11 cm
S/S Str.



RSD4-424
10 cm
RSD4-425
11 cm
S/S Cvd.



RSD4-426
10 cm
RSD4-427
11 cm
B/B Str.



RSD4-428
10 cm
RSD4-429
11 cm
B/B Cvd.



Graefe Scissors
RSD4-430 10.5 cm
B/B Angled

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножици для ясен

Tijeras para encias

مقص صمغ



Cottle Scissors
RSD4-431 16 cm
B/B, Angled




Fomon Scissors
RSD4-432 14.5 cm
B/B, Angled



Sailer Scissors
RSD4-433 17 cm
B/B, Angled



Harvey Scissors
RSD4-434 12.5 cm
Saw Edge



Schoemaker-Loth Scissors
RSD4-435 12.5 cm



Eiselsberg Scissors
RSD4-436 11 cm

GUM SCISSORS

Ножницы десневые

Ciseaux a gencives

Zahnfleischscheren

Ножиці для ясен

Tijeras para encias

مقص صمغ




Scissors
RSD4-437 9 cm
RSD4-438 10.5 cm
S/S, Str.



Scissors
RSD4-439 13 cm
RSD4-440 15 cm
S/B, Str.



Scissors
RSD4-441 9 cm
RSD4-442 10.5 cm
S/S, Str.



Scissors
RSD4-443 9 cm
RSD4-444 10.5 cm
S/S, Cvd.



Cuticle Scissors



RSD4-445
10 cm
RSD4-446
11.5 cm
S/S, Str.



RSD4-447
10 cm
RSD4-448
11.5 cm
S/S, Cvd.



Cuticle Scissors



RSD4-449
9 cm
RSD4-450
10 cm
RSD4-451
11.5 cm
S/S, Str.



RSD4-452
9 cm
RSD4-453
10 cm
RSD4-454
11.5 cm
S/S, Cvd.

HAEMOSTATIC FORCEPS

Зажимы артериальные Затискачі артеріальні
 Pincers hemostatiques Pinzas hemostaticas
 Arterienklemmen ملفط هيموستاتيک



HAEMOSTATIC FORCEPS

■ Зажимы артериальные ■ Затискачі артеріальні
 ■ Pincers hemostaticques ■ Pinzas hemostaticas
 ■ Arterienklemmen ■ ملفط هيموستاتيک



Heiss Forceps
20 cm



RSD4-817

RSD4-818

RSD4-819




RSD4-820



RSD4-821



Rochester-Carmalt Forceps



RSD4-822

16 cm

RSD4-823

18 cm

RSD4-824

20 cm

RSD4-825

23 cm

RSD4-826

16 cm

RSD4-827


18 cm

RSD4-828


20 cm

RSD4-829

23 cm



Crile Haemostatic Forceps



RSD4-830

14 cm

RSD4-831

16 cm

RSD4-832

18 cm

RSD4-833

20 cm


RSD4-834

22 cm

RSD4-835

24 cm

straight



RSD4-836

14 cm

RSD4-837

16 cm

RSD4-838

18 cm

RSD4-839

20 cm


RSD4-840

22 cm

RSD4-841

24 cm

curved



Crile Haemostatic
Tissue Forceps



RSD4-842

14 cm

RSD4-843

16 cm

RSD4-844

18 cm

RSD4-845

20 cm

RSD4-846

22 cm

RSD4-847

24 cm

1x2 teeth

(Kocher) straight



RSD4-848

14 cm

RSD4-849

16 cm

RSD4-850

18 cm

RSD4-851

20 cm

RSD4-852

22 cm

RSD4-853

24 cm

1x2 teeth



(Kocher) curved

HAEMOSTATIC FORCEPS


Зажимы артериальные Затискачі артеріальні
 Pincers hemostaticques Pinzas hemostaticas
 Arterienklemmen ملفط هيموستاتيک




**Pean Haemostatic Forceps
Slender Type**




**Spencer-Wells Haemostatic Forceps
Slender Type**



**Kocher Haemostatic
Forceps Strong Pattern
1x2 teeth**



**Rochester-Carmalt
Haemostatic Forceps**



HAEMOSTATIC FORCEPS

Зажимы артериальные Затискачі артеріальні
 Pincers hemostaticques Pinzas hemostaticas
 Arterienklemmen ملفط هيموستاتيک



**Jones
Towel Forceps**
 RSD4-874 9 cm
 Slight Curve



**Engl. Modell
Towel Forceps**
 RSD4-875 9 cm



**Schaedel Towel
Forceps**
 RSD4-876 9 cm



Mayo
 RSD4-877 14cm



Backhaus Clip Towel Forceps
 RSD4-878 11 cm
 with clip for tubing or cable



Roeder Towel Forceps
 RSD4-879 13 cm
 RSD4-880 15 cm
 Curved




Bernhard Towel Forceps
 RSD4-881 16.5 cm

HAEMOSTATIC FORCEPS

Зажимы артериальные Затискачі артеріальні
 Pincers hemostaticques Pinzas hemostaticas
 Arterienklemmen ملفظ هيموستاتيک



Pratt Haemostatic Forceps
RSD4-882 15 cm



Pennington / Siezing Haemostatic Forceps
RSD4-883 15.5 cm
RSD4-884 20 cm
Straight



Pean Forceps

RSD4-885
12.5 cm

RSD4-886
14 cm

RSD4-887
16 cm



Pean Forceps

RSD4-888
16 cm
RSD4-889
18 cm



Lister Sinus Forceps

RSD4-890
13cm

RSD4-891
15.5cm



RSD4-892
17.5cm

HAEMOSTATIC FORCEPS

Зажимы артериальные Затискачі артеріальні
 Pincers hemostaticues Pinzas hemostaticas
 Arterienklemmen ملفط هيموستاتيک



**Foerster-Ballenger Sponge Forceps
& (bullet holding forceps)**



RSD4-893 18 cm
RSD4-894 20 cm
RSD4-895 25 cm
 Straight Smooth

RSD4-896 18 cm
RSD4-897 20 cm
RSD4-898 25 cm
 Straight Serrated

RSD4-899 18 cm
RSD4-900 20 cm
RSD4-901 25 cm
 Curved Smooth

RSD4-902 18 cm
RSD4-903 20 cm
RSD4-904 25 cm
 Curved Serrated

STERILIZING FORCEPS



RSD4-905 20 cm
RSD4-906 25 cm
RSD4-907 30 cm



Rogge
RSD4-908 20 cm
RSD4-909 30 cm



Cheatel
 Box Joint
RSD4-910
 27 cm




Cheatel
 Screw Joint
RSD4-911
 27 cm


HAEMOSTATIC FORCEPS

Зажимы артериальные
Pincers hemostaticques
Arterienklemmen

Затискачі артеріальні
Pinzas hemostaticas
ملقط هايموستاتيک




Backhaus
RSD4-912
Not-perforating and
not-damaging 14cm



Backhaus
RSD4-913 11cm
RSD4-914 13cm
RSD4-915 14cm
RSD4-916 15cm



Backhaus
RSD4-917
9cm



Gross Maier Dressing Forceps
RSD4-918 16cm straight
RSD4-919 20cm straight
RSD4-920 25cm straight



RSD4-921
16cm
RSD4-922
20cm
RSD4-923
25cm
curved



RSD4-924
Young
Tongue Forceps
17cm

Spare rubber




RSD4-925
Collins
Tongue Forceps
17cm

STERILIZING & TISSUE CLAMPS

 Стерилизационные и тканевые зажимы
 Stérilisation et pincettes pour tissus
 Esterilizadores y sujetadores de tejidos


 Стерилізаційні та тканинні затискачі
 Esterilização e grampos de tecido
 تعقيم & المشابك الأنسجة






NEEDLE HOLDERS

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة





Castroviejo



RSD4-940
Halsey
13 cm



Baumgartner



Derf



RSD4-950
Delicate
15cm



RSD4-951
Adson
18cm



Mayo-Hegar
RSD4-952 13cm
RSD4-953 14cm
RSD4-954 15cm
RSD4-955 16cm
RSD4-956 18cm
RSD4-957 20cm
RSD4-958 23cm
RSD4-959 24cm
RSD4-960 26cm
RSD4-961 30cm
Serrated Grooved




Crile-Wood
RSD4-962 15cm
RSD4-963 18cm
RSD4-964 20cm
RSD4-965 22cm
RSD4-966 25cm

NEEDLE HOLDERS


Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة




Webster
Needle Holder

RSD4-967 11 cm
Smooth
RSD4-968 12 cm
Smooth
RSD4-969 13 cm
Smooth
RSD4-970 14 cm
Smooth
RSD4-971 16 cm
Smooth
RSD4-972 18 cm
Smooth
RSD4-973 20 cm
Smooth
RSD4-974 23 cm
Smooth



Heaney Needle Holders

RSD4-975
20 cm Cvd.
RSD4-976
27 cm Cvd.




Crile-Murray
Needle Holders

RSD4-977
15 cm
RSD4-978
18.5 cm
Slender
pattern



Doyen Needle Holders

RSD4-979
12.5 cm



Olsen-Hegar Needle
Holder With Scissors
14cm

RSD4-980
12 cm
RSD4-981
14 cm
RSD4-982
17 cm
RSD4-983
19 cm

RSD4-984
12 cm
smooth
RSD4-985
14 cm
smooth
RSD4-986
17 cm
smooth
RSD4-987
19 cm
smooth



Olsen-Hegar Needle
Holder With Scissors
14cm

RSD4-988
12 cm
RSD4-989
14 cm
RSD4-990
17 cm
RSD4-991
19 cm

RSD4-992
12 cm
smooth
RSD4-993
14 cm
smooth
RSD4-994
17 cm
smooth
RSD4-995
19 cm
smooth



RSD4-996
Gillies Needle
Holder With Scissors
16cm

NEEDLE HOLDERS

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة



Kilner Needle Holders

RSD4-997
13 cm
serrated



Crile
RSD4-998
15 cm



Toennis
RSD4-999
18 cm



Mathieu

RSD4-1000 RSD4-1001
14 cm 17 cm

RSD4-1002
20 cm



Mathieu

RSD4-1003 RSD4-1004
14 cm 17 cm

RSD4-1005
20 cm



Mathieu

RSD4-1006 RSD4-1007
14 cm 17 cm

RSD4-1008
20 cm



Mathieu

RSD4-1009 RSD4-1010
14 cm 17 cm

RSD4-1011
20 cm

NEEDLE HOLDERS

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة

RSD4-1012
14 cm
Curved, Smooth
without lock
curved

RSD4-1013
14 cm
Straight, Smooth
without lock
straight

RSD4-1014
14 cm
Curved
Smooth
with lock

RSD4-1015
14 cm
Curved
Smooth
without lock

RSD4-1016
14 cm
Straight
Smooth
with lock

RSD4-1017
14 cm
Straight
Smooth
without lock

RSD4-1018
18 cm
Curved
Smooth
with lock

RSD4-1019
18 cm
Curved
Smooth
without lock

RSD4-1020
18 cm
Straight
Smooth
with lock

RSD4-1021
18 cm
Straight
Smooth
without lock

Barraquer Needle Holders

Castroviejo Needle Holders

RSD4-1024
13 cm
with lock
straight

RSD4-1026
13 cm
with lock
curved

RSD4-1025
13 cm
without lock
straight

RSD4-1027
13 cm
without lock
curved

RSD4-1028
14 cm
Curved, Smooth
with lock curved

RSD4-1029
14 cm
Straight, Smooth
with lock straight

Castroviejo Olsen Heager /
Needle Holders / Scissors
RSD4-1022 14 cm
RSD4-1023 15 cm

Castroviejo Needle Holders

Barraquer Needle Holders

NEEDLE HOLDERS

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة



RSD4-1030 16 cm
RSD4-1031 20 cm
RSD4-1032 23 cm
RSD4-1033 16 cm
RSD4-1034 20 cm
RSD4-1035 23 cm
Smooth Str. Smooth Cvd.

Micro Needle Holders



RSD4-1036
16 cm
Smooth Cvd.

Micro Needle Holders



RSD4-1037
12 cm
Smooth Str
with lock

Castroviejo Needle
Holders



RSD4-1038
12 cm
Smooth Cvd
with lock

Castroviejo Needle
Holders



RSD4-1039
10.5 cm
Smooth
Str.

Barraquer-Troutman
Needle Holders

RSD4-1040
10.5 cm
Smooth
Cvd.

RSD4-1041
11.5 cm
Smooth
Str.

RSD4-1042
11.5 cm
Smooth
Cvd.



RSD4-1043
13.5 cm
Serrated

Kalt Needle Holders



RSD4-1044
13.5 cm
Serrated
Str.

Arruga Needle Holders



RSD4-1045
13.5 cm
Serrated
Cvd.



RSD4-1046
14 cm

Mathieu-Kocher



RSD4-1047
17 cm



RSD4-1048
Boynton
12 cm

NEEDLE HOLDERS

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة

T.C GOLD




NEEDLE HOLDERS


Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة


T.C GOLD





Mathieu-Kocher




Mathieu



Olsen-Hegar Needle Holder With Scissors T.C



Olsen-Hegar Needle Holder With Scissors



RSD4-1098
Gillies Needle Holder With Scissors
16cm



Castroviejo
RSD4-1099 14cm straight
RSD4-1100 18cm straight
RSD4-1101 14cm curved
RSD4-1102 18cm curved



Castroviejo
RSD4-1103 14cm straight
RSD4-1104 18cm straight
RSD4-1105 14cm curved
RSD4-1106 18cm curved

NEEDLE HOLDERS

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة



RSD4-1107
Lichtenberg-
Ryder T.C
14 cm



RSD4-1108
Lichtenberg-
Ryder T.C
14 cm



RSD4-1109
Lichtenberg-
Ryder T.C
17 cm



RSD4-1110
Lichtenberg-
Ryder T.C
17 cm



RSD4-1111
Barraquer
18 cm



RSD4-1112
Barraquer
18 cm



RSD4-1113
Barraquer T.C
18 cm



RSD4-1114
Barraquer T.C
18 cm



RSD4-1115
Tönnis
18 cm



RSD4-1116
Castroviejo
13 cm
Straight



RSD4-1117
Kalt Needle Holders T.C
14 cm Straight



RSD4-1118
Peet Splinter Forceps
11.5 cm



RSD4-1119
La Grange
Curved



RSD4-1120
All-Purpose Pliers


NEEDLE HOLDERS

Иглодержатели
Porte-aiguilles
Nadelhalter


Голкотримачи
Porta-agujas
اصحاب ابرة




Kilner
Needle Holders T.C



Wertheim Needle Holders T.C
RSD4-1122 20 cm
RSD4-1123 24 cm




Halsey Needle Holders T.C



Baumgartner Needle Holders T.C



Castroviejo Needle Holders T.C



Castroviejo Needle Holders T.C



Mathieu-Ryder Needle
Holders T.C



Crile

BONE ROUNGER

Костные кусачки
Hohlmeißelzangen
Pincas-goubas

Кусачки для кісток
Pinzas-gubias
رونچيوس العظام



RSD4-1450
Blumenthal
15cm -90°



RSD4-1451
Blumenthal
15cm -45°



RSD4-1452
Blumenthal
15cm -30°



RSD4-1453
Friemdan
14cm



RSD4-1454
Friemdan 14.5cm
Thin
Fine



RSD4-1455
Micro-Friedman 15cm
A form "S"
"S" Shaped



RSD4-1456
Micro-Friedman
16cm



RSD4-1457
Friedman
14 cm
30°



RSD4-1458
Friedman
14 cm
45°



RSD4-1459
Micro-Friedman
14cm -90°

BONE ROUNGER

Костные кусачки
Hohlmeißelzangen
Pincès-gouges

Кусачки для кісток
Pinzas-gubias
رونجيورس العظام



RSD4-1460
Kleinert-Kutz
16cm



RSD4-1461
Cohen 10cm
Tissue Nippers



RSD4-1462
Cohen 14cm
Tissue Nippers



RSD4-1463
Beyer
18.5cm



RSD4-1464
Clevelanda
14cm

BONE ROUNGER

Костные кусачки
Hohlmeißelzangen
Pincers-gouges

Кусачки для кісток
Pinzas-gubias
رونچيوس العظام



RSD4-1465
Rongeurs
16.5cm



RSD4-1466
Rongeurs
Side Cutting
16.5cm



RSD4-1467
Round Beaks
Double Action
18 cm



RSD4-1468
Ruskin-Mini
15.5 cm
Bone rongeurs,
multiple action



RSD4-1469
Bohler
16 cm
Straight

BONE ROUNGER

Костные кусачки
 Hohlmeißelzangen
 Pincès-gouges

Кусачки для кісток
 Pinzas-gubias
 روتجیورس العظام



RSD4-1470
Beyer
18 cm



RSD4-1471
Böhler
15 cm
Bone cutting forceps,
multiple action



RSD4-1472
Fig. 50
Liston
14cm



RSD4-1473
Fig. 51
Liston
19cm



Kerrison
15 cm

RSD4-1474
2 mm

RSD4-1475
3 mm

RSD4-1476
5 mm




RSD4-1477
Tomati
12cm

BONE ROUNGER

- | | |
|------------------|--------------------|
| Костные кусачки | Кусачки для кісток |
| Hohlmeißelzangen | Pinzas-gubias |
| Pinces-gouges | رونچيوس العظام |



RSD4-1478
Universal



RSD4-1479
Mead 1 A
16.5 cm



RSD4-1480
Mead 2 A
16.5 cm



RSD4-1481
Luer
14 cm



RSD4-1482
Luer
14 cm



RSD4-1483
Luer 15.5cm




RSD4-1484
Luer, 17.5cm

BONE ROUNGER

Костные кусачки	Кусачки для кісток
Hohlmeißelzangen	Pinzas-gubias
Pinces-gouges	رونچيورس العظام



RSD4-1485
Adson



RSD4-1486
Cleveland 3
14 cm



RSD4-1487
Cleveland 4
14 cm



RSD4-1488
Cleveland 4 A
17 cm



RSD4-1489
Cleveland 5
Side cutting
17 cm



RSD4-1490
Cleveland 5 A
Side Cutting
17 cm



RSD4-1491
Cleveland 5 S
15 cm
Side Cutting

BONE ROUNGER

Костные кусачки
Hohlmeißelzangen
Pincers-gouges

Кусачки для кісток
Pinzas-gubias
رونچيوس العظام



RSD4-1492
Cleveland 10
17 cm



RSD4-1493
Cleveland-Gardner
end Cutting
17.5 cm



RSD4-1494
Jansen
17.5 cm



RSD4-1495
Jansen
18 cm



RSD4-1496
Bane
Slight Curved
17.5cm






RSD4-1497
Lempert Delicate
16cm

RUBBER DAM CLAMPS FORCEPS

Инструменты для системы коффердам
Pinces porte-crampons a digue
Kofferdam Klammerzangen

Инструменты для системы коффердам
Pinzas de engrapar para dique de caucho
السداد المطاط الملقط المشابيك



RUBBER DAM CLAMPS FORCEPS

Инструменты для системы коффердам
Pincas porte-crampons a digue
Kofferdam Klammerzangen

Инструменты для системы коффердам
Pinzas de engrapar para dique de caucho
السداد العطار المعلق المشابك



Rubber Dam Kits





RSD4-1148

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Rubber Dam Frame L
- 05 Rubber Dam Frame M
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7




RSD4-1149

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Ivory
- 02 Ivory light weight
- 03 Clamps Stand
- 04 Rubber Dam Frame M
- 05 Rubber Dam Frame L
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 209
- 09 Rubber Dam Clamp Fig. 212
- 10 Rubber Dam Clamp Fig. 211
- 11 Rubber Dam Clamp Fig. 8
- 12 Rubber Dam Clamp Fig. 8A
- 13 Rubber Dam Clamp Fig. 14
- 14 Rubber Dam Clamp Fig. 14A
- 15 Rubber Dam Clamp Fig. 0
- 16 Rubber Dam Clamp Fig. 00
- 17 Rubber Dam Clamp Fig. 2
- 18 Rubber Dam Clamp Fig. 2A
- 19 Rubber Dam Clamp Fig. 7



RSD4-1150
Rubber Dum Template



RSD4-1150-A
Rubber Dam Set of 12 Pcs

- 01 Anisworth
- 02 Rubber Dam Frame
- 03 Ivory
- 04 Clamp
- 05 Clamp
- 06 Clamp
- 07 Clamp
- 08 Clamp
- 09 Clamp
- 10 Clamp
- 11 Clamp
- 12 Clamp

RUBBER DAM CLAMPS

Резиновая дамба зажимы комплекты
 Trousses de caoutchouc
 Gummidammklammern

Гумова дамба затискачі комплекти
 Conjuntos de abrazaderas de retención de caucho
 المطاط السد المشابك مجموعات



RSD4-1151

Stainless Steel Gold Plated Rubber Dam Clamp Tray with 11 Different Clamps
 Fig. 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7




RSD4-1152

Stainless Steel Rubber Dam Clamp Tray with 12 Different Clamps
 Fig. 209, 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



RSD4-1153

Stainless Steel Rubber Dam Clamp Tray with 8 Different Clamps
 Fig. 1, 26, 29, 51, 60, 202, 205, 210



RSD4-1154

Stainless Steel Rubber Dam Clamp Tray
 (without clamps)



RSD4-1155

Anisworth 17cm
 Ivory 17cm
 Brewer 17.5cm
 Rubber Dam Frame (Children)
 Rubber Dam Frame (Adult)
 Rubber Dam Clamps 11 Pcs
 Fig. 00, Fig.0, Fig.1, Fig.2A, Fig.3, Fig.7,
 Fig.8, Fig.8A, Fig.9, Fig.212, Fig.14

RUBBER DAM CLAMPS

 Кламмера для коффердама
 Crampons pour digue
 Kofferdam Klammern

 Кламмера для коффердама
 Grapas para dique de caucho
 المشابك المد المطاط



RSD4-1156
Fig. 15



RSD4-1157
Fig. 18



RSD4-1158
Fig. 24



RSD4-1159
Fig. 25



RSD4-1160
Fig. 26



RSD4-1161
Fig. 28



RSD4-1162
Fig. 30



RSD4-1163
Fig. 31



RSD4-1164
Fig. 51



RSD4-1165
Fig. 138



RSD4-1166
Fig. 139



RSD4-1167
Fig. 200



RSD4-1168
Fig. 201



RSD4-1169
Fig. 202



RSD4-1170
Fig. 203



RSD4-1171
Fig. 204



RSD4-1172
Fig. 205



RSD4-1173
Fig. 3



RSD4-1174
Fig. W3



RSD4-1175
Fig. 4



RSD4-1176
Fig. 56A



RSD4-1177
Fig. 7



RSD4-1178
Fig. 7A



RSD4-1179
Fig. 8



RSD4-1180
Fig. 8A



RSD4-1181
Fig. W8A



RSD4-1182
Fig. 10



RSD4-1183
Fig. 11



RSD4-1184
Fig. 12A



RSD4-1185
Fig. 13A



RSD4-1186
Fig. 14



RSD4-1187
Fig. 14A



RSD4-1188
Fig. 14AN



RSD4-1189
Fig. 22



RSD4-1190
Fig. 27



RSD4-1191
Fig. 207



RSD4-1192
Fig. 206



RSD4-1193
Fig. 208



RSD4-1194
Fig. 209



RSD4-1195
Fig. 1



RSD4-1196
Fig. 1A



RSD4-1197
Fig. 2



RSD4-1198
Fig. 2A



RSD4-1199
Fig. W2A



RSD4-1200
Fig. W2A



RSD4-1201
Fig. 0



RSD4-1202
Fig. 00



RSD4-1203
Fig. 210



RSD4-1204
Fig. 211



RSD4-1205
Fig. 212



RSD4-1206
Fig. 6



RSD4-1207
Fig. 9



RSD4-1208
Fig. N06



RSD4-1209
Fig. U67



RSD4-1210
Fig. P1



RSD4-1211
Fig. P2



RSD4-1212
Fig. G-1



RSD4-1213
Fig. G-2

CROWN REMOVER


Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменти для видалення коронок
Arranca-corona
مزيلات التاج

Complete With 3 Attachment



RSD4-1214
Wynman Crown Gripper
No Lock T.C 16 cm



RSD4-1215
Set of 4 Pcs



RSD4-1216
Crown Spreader, Crown Remover
Special Instrument for crown spreading
after separating, simple and safe
operation.



RSD4-1217
Crown Remover
Set of 3 Pcs



RSD4-1218
Fig. 7



RSD4-1219
Fig. 8



RSD4-1220
Fig. 9



RSD4-1221
Fig. 10

For Molars and Premolars



RSD4-1222
Fig. 11



RSD4-1223
Fig. 12



RSD4-1224
Fig. 13
Tofflemire

CROWN REMOVER

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменти для видалення коронок
Arranca-corona
مزيلات الناج



RSD4-1225
Crown Spreaders



RSD4-1226
Crown Spreaders
with Spring



RSD4-1227
Crown Spreader



RSD4-1228
Fig. 1



RSD4-1229
Fig. 2



RSD4-1230
Crown Spreaders
Fig. 3



RSD4-1230-A
Crown Remover
Set of 5 Pcs



RSD4-1231
Crown Remover
Set of 3 Pcs



RSD4-1232
Crown Remover
Set of 5 Pcs



RSD4-1233
Crown Remover
Set of 6 Pcs



CROWN REMOVER

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменты для видаления коронок
Arranca-corona
مزيلات الشاج



RSD4-1498
16cm

Crown Plier with plastic inserts Straight
For ceramical and plastic crowns and bridges, with
plastic inserts (sterilization up to 180° C).



RSD4-1499
14.5cm

Crown Plier
with plastic inserts
Angled



RSD4-1500
14cm

Universal plier for the removal
of ceramic crowns and bridges



RSD4-1501
14cm

Universal plier for the removal
of ceramic crowns and bridges



RSD4-1502
Furrer
Copper Ring
Removing Plier.
15cm
Curved



RSD4-1503 12.5 cm
RSD4-1504 15.0 cm
Böhm
Pliers for telescope
crowns.
Diamond points
coated.



RSD4-1505
Nerve-Canal Pliers
straight



RSD4-1506
Nerve-Canal Pliers
Curved

MATRIX RETAINERS

■ Матрицедержатели ■ Матрицеутримувачи
 ■ Porte Matrices ■ Porta Matrices
 ■ Matrizenspanner ■ مصفوفة الخدم



RSD4-1249
Fig. 5
Ivory 8
5 mm



RSD4-1250
Fig. 6
Ivory 8
6 mm



RSD4-1251
Fig. 7
Ivory 8
7 mm



RSD4-1252
Ivory 8 N



RSD4-1253
Didner 8 P



RSD4-1254
Fig. 1
Tofflemire-Universal



RSD4-1255
Fig. 2
Tofflemire-Junier



RSD4-1256
Fig. 3
Tofflemire-Senior

Amalgam Well



RSD4-1257
Stainless Steel
Amalgam Well
Medium



RSD4-1258
Stainless Steel
Amalgam Well
large

MATRIX RETAINERS

■ Матрицедержатели ■ Матрицеутримувачи
 ■ Porte Matrices ■ Porta Matrices
 ■ Matrizenspanner ■ مصفوفة الخدم



RSD4-1259
Tofflemire-Junior
45 mm



RSD4-1260
Tofflemire-Universal
60 mm



RSD4-1261
Matrix Ivory 8
6 mm



RSD4-1262
Didner 8 P
70 mm



RSD4-1263
Nyström
Left



RSD4-1264
Nyström
Right



RSD4-1265
Narrow



RSD4-1266
Wide



RSD4-1267
Ivory
with shields



RSD4-1268
Ivory
without shields



RSD4-1269
Ivory



RSD4-1270
Elliot

AMALGAM INSTRUMENTS

Инструменты для амальгамы
Instruments pour amalgame
Amalgamisinstrumente

Инструменты для амальгами
Instrumentos para amalgama
الصكوك الملغم

Ø 2mm



RSD4-1271
Fig. 10
16.5cm
English pattern

Ø 2mm



RSD4-1272
Fig. 1
18.5cm
German pattern

Ø 2mm



RSD4-1273
Fig. 1
15cm
Swiss pattern

Ø 2mm



RSD4-1274
Fig. 1
15cm

Ø 2mm



RSD4-1275
Fig. 2
15cm, 45°

Ø 2mm



RSD4-1276
Fig. 3
14cm, 90°



RSD4-1277
Chip Blower



RSD4-1278



RSD4-1279



RSD4-1280



RSD4-1281
Water Syringe

AMALGAM INSTRUMENTS

■ Инструменты для амальгамы
 ■ Instruments pour amalgame
 ■ Amalgamisinstrumente

■ Инструменты для амальгамы
 ■ Instrumentos para amalgama
 ■ الصكوك الملعقم

AMALGAM CARRIERS

TIP SIZES FOR DISTAL OR SINGLE-END

1.5 mm

2.0 mm

2.8 mm

1.5 mm

2.0 mm



RSD4-1282
Mini/
Single-End



RSD4-1283
Regular
Single-End



RSD4-1284
Large
Single-End



RSD4-1288
Mini
Distal



RSD4-1289
Regular
Distal

TIP SIZES FOR DOUBLE-END

1.5 mm

2.8 mm

3.2 mm



RSD4-1285
Mini/
Regular



RSD4-1286
Large/
Regular



RSD4-1287
Jumbo/
Regular



RSD4-1290
Amalgam
Well

Single-End

TIP SIZES FOR DOUBLE-END

RSD4-1291 - 1.5mm-2.0mm **RSD4-1294** - 1.0mm-2mm
RSD4-1292 - 2.0mm-2.8mm **RSD4-1295** - 2.0mm-3mm
RSD4-1293 - 2.0mm-3.2mm **RSD4-1296** - 2.0mm-4mm

RSD4-1297 - 1.5mm-2.0mm **RSD4-1300** - 1.0mm-2mm
RSD4-1298 - 2.0mm-2.8mm **RSD4-1301** - 2.0mm-3mm
RSD4-1299 - 2.0mm-3.2mm **RSD4-1302** - 2.0mm-4mm

RSD4-1303 - 1.5mm-2.0mm **RSD4-1306** - 1.0mm-2mm
RSD4-1304 - 2.0mm-2.8mm **RSD4-1307** - 2.0mm-3mm
RSD4-1305 - 2.0mm-3.2mm **RSD4-1308** - 2.0mm-4mm

RSD4-1309 - 1.5mm-2.0mm **RSD4-1312** - 1.0mm-2mm
RSD4-1310 - 2.0mm-2.8mm **RSD4-1313** - 2.0mm-3mm
RSD4-1311 - 2.0mm-3.2mm **RSD4-1314** - 2.0mm-4mm

RSD4-1315 - 1.5mm
RSD4-1316 - 2.0mm
RSD4-1317 - 2.0mm
RSD4-1318 - 3.0mm

ROOT PLANNING CHISELES & BONE FILES



BONE FILES



CHISELS



Periodontal Chisels

Bone Splitting Chisels



PERIODONTAL INSTRUMENTS

 Инструменты пародонтологические
 Instruments pour la parodontologie
 Instrumente für die Parodontologie

 Інструменти пародонтологічні
 Instrumentos para periodoncia
 الصكوك اللثة



RSD4-1348
1
Ochsenbin



RSD4-1349
2
Ochsenbin



RSD4-1350
Wakefield 1W
17cm



RSD4-1351
Wakefield 2W
17cm



Partsch 14 cm
RSD4-1357 3mm
RSD4-1358 4mm
RSD4-1359 5mm
RSD4-1360 6mm



Partsch 17 cm
RSD4-1361 3mm
RSD4-1362 4mm
RSD4-1363 5mm
RSD4-1364 6mm



Partsch 14 cm
RSD4-1365 3mm
RSD4-1366 4mm
RSD4-1367 5mm
RSD4-1368 6mm



Partsch 17 cm
RSD4-1369 3mm
RSD4-1370 4mm
RSD4-1371 5mm
RSD4-1372 6mm

PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

PERIOSTEEL ELEVATORS



RSD4-1373
Buser



RSD4-1374
Ohi 10



RSD4-1375
Howard 9H



RSD4-1376
Allen 9A



RSD4-1377
Prichard 3



RSD4-1378
Seldin 23



RSD4-1379



RSD4-1380
Raspator, Raspatory
Length 17.5cm



RSD4-1381



RSD4-1382




RSD4-1383
Raspator, Williger
Length 17cm, different widths
(4mm, 5mm or 7mm)



RSD4-1384



RSD4-1385
Trelat
14cm



PERIODONTAL INSTRUMENTS

 Инструменты пародонтологические
 Instruments pour la parodontologie
 Instrumente für die Parodontologie

 Інструменти пародонтологічні
 Instrumentos para periodoncia
 الصكوك للثة



RSD4-1386
Freer



RSD4-1387
Freer 15



RSD4-1388
Freer 16



RSD4-1389
Goldman-Fox
14



RSD4-1390
P 24G



RSD4-1391
Molt
17.5cm



RSD4-1392
Molt
17.5cm



RSD4-1393
Meat
17.5cm




RSD4-1394
Meat
17.5cm



RSD4-1395
Meat
17.5cm

BLACK TITANIUM COLOUR COATED INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD6-041
Fig. 1

RSD6-042
Fig. 2

RSD6-043
Fig. 3
Extra-Flex

RSD6-044
Micro-Mini

RSD6-045
Mini 1

RSD6-046
Mini 4
Extra-Flex


RSD6-047
Fig. 4
Extra-Flex

RSD6-048
Fig. 5
Flexi-Thin

RSD6-049
Fig. 6
Flexi-Thin

RSD6-050
Mini 3
Extra-Flex

BLACK TITANIUM COLOUR COATED INSTRUMENTS Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD6-051
AB1
Boghosian

RSD6-052
AB2
Boghosian

RSD6-053
W3

RSD6-054
8A

RSD6-055
4/5
Gregg

RSD6-056
Interproximal
Carver

RSD6-057
Interproximal
Carver, Long

RSD6-058
A6

RSD6-059
Contouring
instrument

BLACK TITANIUM COLOUR COATED INSTRUMENTS Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

BURNISHERS

FREEDMAN BURNISHERS



RSD6-060
21B



RSD6-061
27/29



RSD6-062
2
Ladmore



RSD6-063
3
Ladmore



RSD6-064
Freedman
"Duckhead"



RSD6-065
Freedman
Small Contact
Forming




RSD6-066
Freedman
Large Contact
Forming



RSD6-067
Small/
Medium Ball
Burnisher

BLACK TITANIUM COLOUR COATED INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD6-068
A
Cosmetic
Contouring

RSD6-069
B
Cosmetic
Contouring

RSD6-070
C
Cosmetic
Contouring

RSD6-071
D
Cosmetic
Contouring

RSD6-072
E
Cosmetic
Contouring


RSD6-073
F
Cosmetic
Contouring

RSD6-074
Marginal Ridge
& Embrasure
Shaping Instrument

RSD6-075
Occlusal
Anatomy
Instrument


RSD6-076
Composite
Packing
Instrument

BLACK TITANIUM COLOUR COATED INSTRUMENTS Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




TITANIUM PLASMA COLOUR COATED SINUS LIFT INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



TITANIUM PLASMA COLOUR COATED SINUS LIFT INSTRUMENTS Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD6-092
Sinus Lift Elevator


RSD6-093
Sinus Lift Elevator

RSD6-094
Sinus Lift Elevator

RSD6-095
Sinus Lift Elevator

TITANIUM PLASMA COLOUR COATED IMPLANT SINUS LIFT INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD6-096
Acute 6577

RSD6-097
Sinus1

RSD6-098
Sinus2

RSD6-099
Sinus3

RSD6-100
Sinus4


RSD6-101
Sinus5

RSD6-102
Sinus6

RSD6-103
Sinus7

TITANIUM PLASMA COLOUR COATED SINUS LIFT INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



RSD6-104
Sinuslift 1

RSD6-105
Sinuslift 2

RSD6-106
Sinuslift 3

RSD6-107
Sinuslift 4

RSD6-108
Sinuslift 5

RSD6-109
Sinuslift 6

RSD6-110
Membrane
3Positioning 1

RSD6-111
Condenser 1


TITANIUM PLASMA COLOUR COATED INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

Sinus Lift Instrument,
Flexible Sinus Lift instruments
reduce the risk of perforation of the
"Schneider's membrane"

Sinus Lift Elevator,
9.5mm Straight grip
Length 175 mm, Sinus lift elevators, dull,
for gentle peeling and lifting the
Schneider membrane

Sinus lift instrument
9.5mm Straight grip Length 180 mm,
Buser (B1, B2, B3) raspatories,
quality medical stainless steel,
corrosion-resistant



RSD6-112
Fig. 1

RSD6-113
Fig. 2

RSD6-114
Fig. 3

RSD6-115
Fig. 1

RSD6-116
Fig. 2

RSD6-117
Fig. 3

RSD6-118
B01


RSD6-119
B02

RSD6-120
B03

TITANIUM PLASMA COLOUR COATED INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

Sinus lift instrument 9.5mm Straight grip Sinus lift elevators, blunt or sharply curved (H01, H02, K03, K04, K78, K77)



RSD6-121
Sinus lift
instrument
H01

RSD6-122
Sinus lift
instrument
H02

RSD6-123
Sinus lift
instrument
K03

RSD6-124
Sinus lift
instrument
K04

RSD6-125
Sinus lift
instrument
K78


RSD6-126
Sinus lift
instrument
K77

TITANIUM PLASMA COLOUR COATED INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SINUS LIFT INSTRUMENT

Sinus Lift Elevators, 9.5mm Hollow Handle length 170mm Blunt Curved Set of 6 Pcs



RSD6-127

RSD6-128

RSD6-129

RSD6-130

RSD6-131

RSD6-132

IMPLANTOLOGY

Инструменты для имплантологии

Implantologie

Implantologie

Инструменти для імплантології

Implantologia

زراعة الأسنان

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

PERIODONTAL SURGERY SURGICAL CHISELS

KRAMER NEVINS SINUS LIFT



RSD6-133 TG



RSD6-134 SOLT3/4



RSD6-135 TGO



RSD6-136 29K



RSD6-137 FEDI4



RSD6-138 OCH3



RSD6-139 OCH4



RSD6-140 FEDI1



RSD6-141 FEDI2



RSD6-142 FEDI3



RSD6-143 Gf36

IMPLANTOLOGY

Инструменты для имплантологии

Implantologie

Implantologie

Инструменти для имплантології

Implantologia

زراعة الأسنان

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
IMPLANT SURGICAL CHISELS



RSD6-144 13K/TG



RSD6-145 13K/TG Titanium Plasma Coated Point



RSD6-146 13K/TG



RSD6-147 13KL/R Titanium Plasma Coated Point



RSD6-148 13KL 4/6mm



RSD6-149 Pt2 Titanium Plasma Coated Point
A/V Regular



RSD6-150 Pt1 Titanium Plasma Coated Point
A/V Regular



RSD6-151 Pt3 Titanium Plasma Coated Point Angel S/E



RSD6-152 Pt4 Titanium Plasma Coated Point STR S/E



RSD6-153 Pt5 Titanium Plasma Coated Point TG CVD S/E

PERIODONTAL INSTRUMENTS

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
المسكوك للثة



RSD4-1352
Mead
19cm

ø 28 mm
Teflon
interchangeable
inserts
Autoclavable.



RSD4-1353
Mead
17cm

ø 26 mm
With
interchangeable
inserts
320 gr – 11 ¼ oz



RSD4-1354
Partson
18cm

ø 22 mm
Lead filled
200 gr – 7 oz



RSD4-1355
Vickers
19 cm



RSD4-1356
Herbert
19 cm

DENTAL SYRINGES

Шприцы стоматологические карпульные
Seringues à usage dentaire
Dentalspritzen


Шприці стоматологічні карпульні
Jeringas dentales
المحاقن الأسنان




Type A
RSD4-1396 1.8ml
RSD4-1397 2.2ml



Type C
RSD4-1398 1.8ml
RSD4-1399 2.2ml



Type C
RSD4-1400 1.8ml
RSD4-1401 2.2ml



Type Euro
RSD4-1402 1.8ml
RSD4-1403 2.2ml




3-Ring
RSD4-1404 1.8ml
RSD4-1405 2.2ml



2-Ring
RSD4-1406 1.8ml
RSD4-1407 2.2ml



Euro
Breech loading
RSD4-1408 1.8ml
RSD4-1409 2.2ml



Euro
Breech loading
RSD4-1410 1.8ml
RSD4-1411 2.2ml



Manual Aspirating Syringe
RSD4-1412 1.8ml
RSD4-1413 2.2ml



Quatro
RSD4-1414 1.8ml
RSD4-1415 2.2ml



Trio
RSD4-1416 1.8ml
RSD4-1417 2.2ml



RSD4-1417-A

DENTAL SYRINGES

 Шприцы стоматологические карпульные
 Seringues à usage dentaire
 Dentalspritzen

 Шприці стоматологічні карпульні
 Jeringas dentales
 المحاقن الأسنان



RSD4-1418 1.8ml.
RSD4-1419 2.2ml.
 Pistol-Ject
 Injects 0.20ml with every click
 Lever action for high injection pressure

Autoclavable



RSD4-1420 1.8ml
RSD4-1421 2.2ml
 Fig. 1001
 Intraligamental
 Anesthesia Syringe in
 Black Colour Auto
 Cleavable

Autoclavable



Ject Ultima
RSD4-1422 1.8ml
RSD4-1423 2.2ml
 Injects 0.05 ml with every click
 Supplied with angled and straight tubes



Swedish Type
RSD4-1424 1.8ml
RSD4-1425 2.2ml



Ject angled
RSD4-1426 1.8ml
RSD4-1427 2.2ml



Ject straight
RSD4-1428 1.8ml
RSD4-1429 2.2ml



Disc System
RSD4-1430 1.8ml
RSD4-1431 2.2ml



Disc System
RSD4-1432 1.8ml
RSD4-1433 2.2ml



Disc System
RSD4-1434 1.8ml
RSD4-1435 2.2ml



Projection System
RSD4-1436 1.8ml
RSD4-1437 2.2ml



Projection System
RSD4-1438 1.8ml
RSD4-1439 2.2ml

DENTAL SYRINGES

 Шприцы стоматологические карпульные
 Seringues à usage dentaire
 Dental spritzen

 Шприці стоматологічні карпульні
 Jeringas dentales
 المحاقن الأسنان



Aluminium - Dark
RSD4-1440 1.8 ml
RSD4-1441 2.2 ml



Aluminium - Silver
RSD4-1442 1.8 ml
RSD4-1443 2.2 ml



Aluminium - Orange
RSD4-1444 1.8 ml
RSD4-1445 2.2 ml






Aluminium - Green
RSD4-1446 1.8 ml
RSD4-1447 2.2 ml



Aluminium - Blue
RSD4-1448 1.8 ml
RSD4-1449 2.2 ml

ORTHODONTIC PLIERS & INSTRUMENTS





Orthodontics Pliers in Clinical Dentistry Orthodontics And Dentofacial Orthodontics

Driven by the demand for tooth alignment, general practitioners are increasingly being asked to provide orthodontic treatment.

Dentists keen to participate in this area of expertise need to demonstrate the skills to identify suitable patients, the ability to provide conventional and contemporary orthodontic treatment and be aware of its limitations.

This is exactly what this programme delivers, along with the vital core competencies to excel in this unique postgraduate area.

The objective of the Orthodontics and Dentofacial Orthopaedics is to provide dental surgeons with the core knowledge and clinical skills essential for the delivery of competent, high-quality dental care.

It encompasses a blended learning approach and consists of a combination of lectures, seminars, workshops, clinical case discussions and laboratory sessions.

The flexible modular structure enables the programme to meet a wide range of student needs.



ORTHODONTIC INSTRUMENTS

Инструменты ортодонтические
Instrumentos Ortodonticos
Orthodontische Instrumente

Инструменти ортодонтични
Instruments Orthodontiques
اوتات تلويم الأسنان



RSD5-001
Height Gauge 0.18
Aluminum



RSD5-002
Height Gauge 0.22
Aluminum



RSD5-003
Height Gauge 0.18
Stainless Steel



RSD5-004
Height Gauge 0.22
Stainless Steel



RSD5-005
Moveable
Height Gauge



RSD5-007
Endo Finger
Ruler Span
Measure Scale



RSD5-006
Bone Bracket Gauge



RSD5-008
Endo Finger
Ruler Span
Measure Scale



RSD5-009
Endo Finger
Ruler Span
Measure Scale



RSD5-010
Bone Thickness Gauge
40mm, 18cm



RSD5-011
Osteometer
Measuring Forceps For Jaw
Ridge Width measurement, scale
on both sides, rustfree, autoclavable



RSD5-012
Positive
Postioner



RSD5-013
Strip Saw



RSD5-014
Compasso
De Willis



RSD5-015
Endo Finger Ruler
Span Measure Scale



RSD5-016
Endo Finger Ruler
Span Measure Scale

Bracket Holder Tweezers



RSD5-017
Deluxe Bracket
Holder Tweezers



RSD5-018
Deluxe Bracket
Holder Tweezers



RSD5-019
Direct Bond Bracket
Holder Tweezers



RSD5-020
Bracket Tweezers
Slot Aligner



RSD5-021
Direct Bond
Bracket Holder
Tweezers, Small Tip



RSD5-022
Exact Bracket
Holder Tweezers



RSD5-023
Bracket Placing
Tweezers



RSD5-024
Easy Access
Tweezers



RSD5-025
Bracket Placer
Tweezers
(Extra Fine Tip)



RSD5-026
Deluxe Bracket
Holder Tweezers
13.5cm



RSD5-027
Easy Eccess
Tweezers



RSD5-028
Slipper
Tweezers



RSD5-029
Direct Bonding
Bracket Holders
13.5cm



RSD5-030
Direct Bonding
Bracket Holders
13.6cm



RSD5-031
Direct Bonding
Bracket Holders
13.5cm



RSD5-032
Deluxe Bracket
Holder



RSD5-033
Direct Bonding
Bracket Holders
14cm

ORTHODONTIC INSTRUMENTS

Инструменты ортодонтические
Instrumentos Ortodonticos
Orthodontische Instrumente

Инструменти ортодонтични
Instruments Orthodontiques
ادوات تلويم الأسنان



RSD5-034
Elastic instruments



RSD5-035
Twiel-on elastic module instrument



RSD5-036
Schure Band Seater



RSD5-037
Ortho light scaler



RSD5-038
Ortho light scaler



RSD5-039
Ortho instrument



RSD5-040
Ortho instrument



RSD5-041
Ligature instruments



RSD5-042
Ligature instruments

ORTHODONTIC INSTRUMENTS

Инструменты ортодонтические
Instrumentos Ortodonticos
Orthodontische Instrumente

Инструменти ортодонтични
Instruments Orthodontiques
اڤوات تلويم الاڤستان



RSD5-043
Ortho Light double ended
ligature director



RSD5-044
Ortho Light double ended ligature scaler director
remove modules with one instruments



RSD5-045
Ortho band pusher scaler



RSD5-046
Ortho ligature director single ended



RSD5-047
Ortho Light double ended ligature scaler/director
remove modules with one instruments



RSD5-048
Ortholite Ligature Director Scaler



RSD5-049
Ortho light ligature director double ended



RSD5-050
Double ended Tampers

ORTHODONTIC INSTRUMENTS

Инструменты ортодонтические
Instrumentos Ortodonticos
Orthodontische Instrumente

Инструменти ортодонтични
Instruments Orthodontiques
اوتات تلويم الأسنان



RSD5-051
Band seater
auto cleavable



RSD5-052
Band seater
auto cleavable
Triangle Shape



RSD5-053
Band seater
auto cleavable
Square Shape



RSD5-054
Band seater
auto cleavable
Round Shape



RSD5-055
Mershon Band
Pusher
14cm



RSD5-056
Band seating &
removing
instruments



RSD5-057
Ortho Light double
ended
ligature director
scaler



RSD5-058
Ortho Light Band
Pusher Scaler



RSD5-059
Band Seating &
Removing
Silicon Grip



RSD5-060
Ligature Instrument
Silicon Grip



RSD5-061
Elastic Instrument
Silicon Grip



RSD5-062
Distal Bender
Silicon Grip



RSD5-063
Distal Bender
Silicon Grip



RSD5-064
Distal Bender
Silicon Grip

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس



RSD5-065
 Aderer
 Large wire bending
 and loop forming pliers
 13cm

RSD5-066
 Aderer Plier
 For the manufacturing of appliances and
 activation of wires for wires up to .030" fine
 tip and up to .036" mean tip.
 Length 13cm

RSD5-067
 Angle Plier
 Fig. 1
 13.5cm



RSD5-068
 Angle Plier
 Fig. 2
 14cm

RSD5-069
 Angle Plier
 Fig. 3
 13cm

RSD5-070
 Angle Plier
 Fig. 4
 14cm

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروتيسيس



RSD5-071
 Adams Plier
 Adams pliers for bending wire
 arches and bends right angles
 for wires up to .022", .025 tip
 to beak. Length 13cm

RSD5-072
 Flat Nose Pliers
 14cm

RSD5-073
 Johnson Contouring
 Pliers 14cm



RSD5-074
 Falcon Beak
 Contouring Pliers.
 13cm

RSD5-075
 How Plier
 Universally usable tool for the
 multiband technology for wires
 up to .028".
 Length 14.5cm

RSD5-076
 How Plier curved
 Universally usable tool for
 the multiband technology for
 wires up to .028".
 Length 14.5cm

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كمشة لتقويم الأسنان وبروستيسيس




RSD5-077
De La Rosa
Fig. 1
13.5cm



RSD5-078
Arch Counting Plier Grooved
for double back and triple back
bands in 0.30" / 0.76mm & 0.36"
/ 0.91mm lingual arch wires
13.5cm



RSD5-079
Contouring Plier
Fig. 2
13.5cm



RSD5-080
Nance Pliers
13.5cm



RSD5-081
Tweed Loop
forming Pliers
14cm



RSD5-082
Nance Loop
forming Pliers
14cm

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонтії та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيزيس




RSD5-083
Tweed Loop Plier
13cm



RSD5-084
Band Removing Plier
The curved blade of the band removing plier can strip easily and gently remove directly to the tooth. Length 14cm



RSD5-085
Crown & Shell Plier
13cm



RSD5-086
Young Wire Bending &
Loop Forming Pliers
13.5cm



RSD5-087
Ligature Cutters T.C
Hard Wire
13cm



RSD5-088
Mini Ligature Cutters T.C
Light Tip
13cm


PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس



RSD5-089
Mini Ligature cutters T.C
fine tip
13cm



RSD5-090
Large
Fluch-cut distal and cutting Pliers T.C
with safety hold fitted with
un-breakable, corrosion
resistant and durable
stainless steel tips
13cm



RSD5-091
Small
Fluch-cut distal and cutting Pliers T.C
with safety hold fitted with
un-breakable, corrosion
resistant and durable
stainless steel tips
13cm



RSD5-092
Bracket Removing Plier
Bracket removing plier for comfortable
and quick removal of the brackets.
Length 13.5cm



RSD5-093
Bracket Removing Plier curved
Bracket Removing Plier curved for
easy removal of brackets in the
posterior region. Length 14.5cm



RSD5-094
Grooved Fine Plier
14.5cm

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروتستيس



RSD5-095
Aderer
Facebow loop forming Pliers
14.5cm



RSD5-096
V Bending Pliers
13.5cm



RSD5-097
Tweed / angle rectangular
arch forming & Torquing pliers
13.5cm



RSD5-098
Weingart Plier medium
Weingart plier that are ideal due to
their curved tip for setting sheets and
is capable of bending the ends.
Length 14.5cm



RSD5-099
Universal Plier
16.5cm



RSD5-100
Jarabak Plier
14.5cm

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
Pinces pour l'orthodontie et la prothèse
Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
Alicates para ortodoncia y prótesis
كشاة لتقويم الأسنان بروستيزيس



RSD5-101
Wald sachs Wire Bending
and Cutting Pliers
15cm



RSD5-102
Universal Dolphin Style
Wire Bending & cutting pliers



RSD5-103
Arrow Clasp Forming Pliers
12.5cm

RSD5-104
Schwarz
13cm



RSD5-105
Flat Nose Pliers
both sides smooth
14cm



RSD5-106
Flat Nose Pliers
both sides serrated
14cm



RSD5-107
Peeso
Wire Bending Pliers
13.5cm



RSD5-108
Goslee
Flat Nose Pliers
14cm



RSD5-109
Flat Nose Pliers
both sides serrated
with TC 14cm inserts

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس



RSD5-110
Universal Pliers
for soft and hard wire
16cm

Max. Ø 0,7mm
 Max. Ø 0,9mm



RSD5-111
Thick grooved plier
14.5cm



RSD5-112
Bird's beak claw with cutter
13cm



RSD5-113
Elastomeric Separator
Placing Pliers 16cm



Separators placing pliers
for easy placement of radiopaque
and non-radiopaque separators



RSD5-114
for incisors with
convex beaks
16cm



RSD5-115
For bicusps
and molars
with concave beaks
16cm



PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس


CUTTERS



RSD5-116
Hard
Wire Cutter 15



RSD5-120
Mini Pin &
Ligature Cutter



RSD5-117
Distal End Cutter
Flush Cut



RSD5-121
Hard Wire Cutter
Straight



RSD5-118
Universal Cut & Hold
Distal End Cutter



RSD5-122
Micro-Mini Pin &
Ligature Cutter



RSD5-119
Pin & Ligature
Cutter



RSD5-123
Cutter

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس

CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS




RSD5-124
Weingart Pliers




RSD5-129
Adhesive
Removing Pliers



RSD5-125
Slim Weingart
Pliers



RSD5-130
Posterior Band Removing
Pliers - Long



RSD5-126
How Pliers
Straight




RSD5-131
Posterior Band Removing
Pliers - Short




RSD5-127
How Pliers Offset



RSD5-132
Bird Beak Pliers



RSD5-128
Band seating Pliers



RSD5-133
Omega Loop
Forming Pliers

Replacement tip and
hex key included.

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان وبروستيسيس




RSD5-134
Three Jaw Pliers



RSD5-138
Jarabak Pliers



RSD5-135
Light Wore Bird Beak
Pliers with Groove



RSD5-139
Hollow Chop
Contouring Pliers




RSD5-136
Tweed loop
Forming Pliers



RSD5-140
Lingual Arch
Forming Pliers



RSD5-137
Arch Bending Pliers



RSD5-141
Rectangular Arch
Bending Pliers

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كمشة لتقويم الأسنان وبروستيس

L-KEY JOINT PLIERS T.C.



RSD5-142
T.C Distal end cutter with safety hold, mini head, long handles.
fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



RSD5-146
T.C Soft wire and ligature cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

⌋ .25 mm to .38 mm (.010" to .015")



RSD5-143
T.C Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" to .025")



RSD5-147
T.C Soft wire and ligature cutter, 45° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

⌋ .25 mm to .38 mm (.010" to .015")



RSD5-144
T.C Soft wire and ligature cutter, Mini fitted with un-breakable, corrosion resistant and durable stainless steel tips

⌋ .25 mm to .38 mm (.010" to .015")



RSD5-148
T.C Hard wire cutter, straight Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



RSD5-145
T.C Soft wire and ligature cutter, standard fitted with un-breakable, corrosion resistant and durable stainless steel tips

⌋ .25 mm to .38 mm (.010" to .015")



RSD5-149
T.C Hard wire cutter, 15° angled Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



RSD5-150
T.C Hard wire cutter, 15° angled Non-autoclavable

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
Pinces pour l'orthodontie et la prothèse
Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
Alicates para ortodoncia y prótesis
كماشة لتقويم الأسنان وبروستيسيس

L-KEY JOINT PLIERS T.C



RSD5-151
T.C Band removing pliers



RSD5-156
T.C Weingart-mini
Utility Pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max 0.56 mm x 0.71 mm (.022" x .028")



RSD5-152
T.C Bracket removing pliers angled



RSD5-157
T.C Weingart
Utility Pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max 0.56 mm x 0.71 mm (.022" x .028")



RSD5-153
T.C Bracket removing pliers straight



RSD5-158
T.C How Straight
Utility Pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max 0.56 mm x 0.71 mm (.022" x .028")



RSD5-154
T.C Bracket removing pliers with plastic tip



RSD5-159
T.C How Angled
Utility Pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max 0.56 mm x 0.71 mm (.022" x .028")



RSD5-155
T.C Adhesive removing pliers



RSD5-160
T.C Wire bending and loop forming pliers
Angled (Bird beak)

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max 0.56 mm x 0.71 mm (.022" x .028")

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
Pincettes pour l'orthodontie et la prothèse
Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
Alicates para ortodoncia y prótesis
كماشة لتقويم الأسنان بروتستيس

L-KEY JOINT PLIERS T.C.



RSD5-161
T.C Cinch Back / tie back
Distal end bending pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



RSD5-162
T.C Detailing step pliers / Bayonet
bending pliers - 1mm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



RSD5-163
T.C Stop pliers V-Bebd - 1mm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.43 mm x 0.56 mm (.017" x .028")



RSD5-164
T.C Tweed / angel
Rectangular arch forming and
torquing pliers

● 0.30 mm to 0.71 mm (.012" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



RSD5-165
T.C Surgical ball hook crimping
pliers, straight

L-KEY JOINT PLIERS Without T.C.



RSD5-166
Multi-purpose pliers, bird beak,
De la rosa and cutter in one

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")



RSD5-167
Hammer head
Distal end bending pliers

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



RSD5-168
Nance
Loop forming pliers

● 0.36 mm to 0.71 mm (.014" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



RSD5-169
Lingual / palatal arch forming plies

● 0.76 mm (.030") / 0.91 mm (.036")



RSD5-170
Torquing pliers



RSD5-171
Surgical ball hook crimping pliers,
angled

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
Pincettes pour l'orthodontie et la prothèse
Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
Alicates para ortodoncia y prótesis
كماشة لتقويم الأسنان بروتستیس

L-KEY JOINT PLIERS



RSD5-172
Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/.76mm.



RSD5-178
Falcon beak contouring plier.



RSD5-173
Bird beak plier for bending round, square and rectangle wire up to .30"/.76mm.



RSD5-179
For placement and removal of arch wires, pins and other auxiliaries.




RSD5-174
Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.



RSD5-180
For placement and removal of arch wires, pins and other auxiliaries in angle design



RSD5-175
Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.



RSD5-181
contouring Plier.



RSD5-176
Adams pliers for hooks and bending wire up to .030"/0.76mm.



RSD5-182
Bending plier. In order to skirt molares bands. With grooved. Recommended wire. 016", 018" and 022" 0.41mm, 0.46mm and 0.56mm



RSD5-183
Loop forming in four length 3,4,5 & 6 mm wire capacity 0.22", 56mm.



RSD5-177
Johnson pliers for crown band contouring



RSD5-184
Arch bending plier design for convenient handling of square or rectangular wire, .020" / 0.51mm.



PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонтії та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس



RSD5-185
Thermo Forming Pliers



RSD5-186
Thermo Forming Pliers



RSD5-187
Thermo Forming Pliers



RSD5-188
Ortho Pliers



RSD5-189
Ortho Pliers



RSD5-190
Ortho Pliers



RSD5-191
Ortho Pliers



RSD5-192
Ortho Pliers

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان وبروستيسيس

LAP & L-KEY JOINT PLIERS



RSD5-193
Pin & Ligature Cutters
TC. Gold.



RSD5-194
Mini Pin & Ligature
Cutter TC. Gold.



RSD5-195
Standard Ligature
Cutter TC. Gold.



RSD5-196
Distal end Cutter
TC. Gold.



RSD5-197
Distal end Cutter
Long Handle
TC. Gold.



RSD5-198
Hook Crimping Plier Angled
TC. Gold.



RSD5-199
Distal end Cutter, Safety Hold Plier
(Ideal for Cutting Wire) TC. Gold.

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincettes pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس

L-KEY JOINT PLIERS



RSD5-200
Direct Bracket Remover
TC. Gold.



RSD5-201
Straight Direct Bond Bracket
Remover T.C



RSD5-202
Pin Holding Plier



RSD5-203
3 Prong Clasp Adjustable
Plier TC. Gold.



RSD5-204
Utility Plier Weingart TC. Gold.



RSD5-205
Nance Plier TC. Gold.

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كماشة لتقويم الأسنان بروستيثيس



RSD5-206
Separating Pliers For Separators
15.5cm



RSD5-207
Elastic
Separating Plier



RSD5-208 Straight
RSD5-209 Curved
Smaha Elastomeric ligature
placing forceps
14cm



RSD5-210
Elastomeric ligature placing
forceps with groove
14cm



RSD5-211
Elastomeric ligature placing
forceps with groove
14cm

PLIERS FOR ORTHODONTICS & PROSTHETICS

Щипцы для ортодонтии и протезирования
 Pincers pour l'orthodontie et la prothèse
 Zangen für Kieferorthopädie und Prothetik

Щипці для ортодонції та протезування
 Alicates para ortodoncia y prótesis
 كمشة لتقويم الأسنان بروستيثيس



RSD5-212
Methieu Plier
With Ring



RSD5-213
Pliers to Place Modules
14cm



RSD5-214
Forceps For Placing Dental astics
12cm



RSD5-215
Mathieu
14cm




RSD5-216
Mathieu Ligature Pliers
14cm



RSD5-217
Ligature Pliers T.C Gold
With Tungsten Carbide Inserts
14cm



Ligature Pliers
14cm



RSD5-218 RSD5-219 RSD5-220

PLIERS FOR ORTHODONTICS & PROSTHETICS

- | | |
|---|--------------------------------------|
| Щипцы для ортодонтии и протезирования | Щипці для ортодонтії та протезування |
| Pincettes pour l'orthodontie et la prothèse | Alicates para ortodoncia y prótesis |
| Zangen für Kieferorthopädie und Prothetik | كمامة لتقويم الأسنان بروستيسيس |




- RSD5-221**
Orthodontic Set of 21 Pcs
- 01 Plier
 - 02 Plier
 - 03 Plier
 - 04 Plier
 - 05 Plier
 - 06 Plier
 - 07 Plier
 - 08 Plier
 - 09 Plier
 - 10 Plier
 - 11 Plier
 - 12 Clamp
 - 13 Plier
 - 14 Mathieu Needle Holder
 - 15 Elevator
 - 16 Gauge
 - 17 Tweezer
 - 18 Band Seater
 - 19 Tweezer
 - 20 Scaler
 - 21 Scaler



- RSD5-222**
Orthodontic Set of 18 Pcs
- 01 Mirror With Handle
 - 02 Plier
 - 03 Scissor
 - 04 Plier
 - 05 Tweezer
 - 06 Plier
 - 07 Matrix Retainer
 - 08 Band Seater
 - 09 Plier
 - 10 Tweezer
 - 11 Plier
 - 12 Probe
 - 13 Plier
 - 14 Plier
 - 15 Plier
 - 16 Needle Holder
 - 17 Plier
 - 18 Plier



- RSD5-223**
Orthodontic Set of 10 Pcs
- 01 Plier
 - 02 Plier
 - 03 Scaler
 - 04 Scaler
 - 05 Mathieu Needle Holder
 - 06 Plier
 - 07 Plier
 - 08 Tweezer
 - 09 Band Seater
 - 10 Plier





**Instruments for dental implants are made of high quality stainless steel.
Our line of instruments include, but is not limited to, the following products:**

Bone Grafting Instruments

- Bone Curettes (SPRATT) – designed for harvesting bone
- Bone Grinder – used in conjunction with a mallet to grind bone into appropriate particle sizes
- Bone Basin – weighted bottom design prevents tipping or spilling while mixing bone or other bone grafting materials
- Bone Compactors – feature serrated tips for compacting bone
- Scalpels – offers convenience of a disposable with razor sharp cutting edge for superior cutting
- Biopsy Punches – designed at a steeper angle to ensure a finely honed, sharper edge
- Bone Spreading Chisels – double or single cut tips used for widening ridges during implant preparations
- Graduated Chisels – feature reference markings for measurement during a bone ridge augmentation and are available in straight and curved patterns

Osteotomes

- Concave Osteotomes – break the floor of the maxillary sinus through the implant cavity and are available in curved (off set) and straight (with or without stops)
- Convex Osteotomes – raise the floor of the maxillary sinus and are available in curved (off set) and straight patterns (with or without stops)
- Progressive Osteotomes – designed to prepare the implant site in the maxillary area where bone density is 3 or 4
- Bone Spreading Osteotomes – sharp or blunt configurations for progressive circumferential expansion of the implant site

Sinus Lift Instruments

- Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

Gauges & Calipers

- Gauges & Calipers – various measuring devices available to determine appropriate length, width, and distance for implant procedures

Collagen Products

- Collagen wound dressing product offering includes products.

Sterilization Trays

- Sterilization Trays – designed to protect and store instruments during sterilization, transportation and storage

IMPLANT GAGES



RSD6-001
Multifunctional Meter
Length 9cm
different sizes,
8-10-11mm,
corrosion-resistant

RSD6-002
Multifunctional Meter
Length 9cm
different sizes,
5-13-15-17 mm,
corrosion-resistant

RSD6-003
Measuring Probe
Length 17cm
(6-8-10-14-16 mm),
corrosion-resistant


RSD6-004
Implant
Depth Gauge

RSD6-005
Implant
Depth Gauge
& Key

RSD6-006
Implant
Depth Gauge

RSD6-007
Implant
Depth Gauge

IMPLANT GAGES



RSD6-008
Double Ended Osteotomy Depth Gauge

RSD6-009
Double Ended Osteotomy Depth Gauge

RSD6-010
Double Ended Osteotomy Depth Gauge

RSD6-011
Depth & Gingival Gauges

Caliper / Calipers Rule



RSD6-012
Castroviejo Calipers
20mm, 45°
9cm



RSD6-013
Castroviejo Calipers
20mm
9cm



RSD6-014
Castroviejo Calipers
20mm
9cm



RSD6-015
Ridge Mapping
Double Caliper



RSD6-016
Fig. 5
Castroviejo
18cm



RSD6-017
Fig. 4
Castroviejo
18cm



RSD6-018
Fig. 3
Castroviejo
18cm



RSD6-019
Clamp, Osteometer
Length about 15.5cm,
ridge width measurement,
two-sided scale,
measuring
range 0 - 25 mm



RSD6-020
Measuring Loop
straight length 9cm,
measuring range
0 - 40 mm,



RSD6-021
Measuring Loop
angled, length 9cm,
measuring range
0 - 40 mm,



RSD6-022
Measuring loop, Castroviejo
angled, length 17.5cm,
measuring range 0 - 40 mm,
with window



RSD6-023
Measuring Loop
angled, length
17.5cm,
measuring range
0 - 40 mm,



RSD6-024
Measuring Loop
straight length
17.5cm,
measuring range
0 - 40 mm,



RSD6-025
Fig. 1
Iwanson
spring caliper
For metal and
porcelain.



RSD6-026
Fig. 2
Iwanson
spring caliper
For wax.



RSD6-027
Bone Thickness Gauge
40mm,
18cm

IMPLANTOLOGY

Инструменты для имплантологии

Implantologie

Implantologie

Инструменти для імплантології

Implantologia


زراعة الأسنان



RSD6-028
Boley Gauge



RSD6-029
Vernier Gauge



RSD6-030
Aluminium Caliper
up to 20cm



RSD6-031
Boley Gauge




RSD6-032
Delux Boley
Gauge



RSD6-033
Beerendonk
13cm
Mini Gauge

Инструменты для имплантологии
 Implantologie
 Implantologie

Инструменти для имплантології
 Implantologia
 زراعة الأسنان



RSD6-034
 Micro Boley
 Gauge Twisted



RSD6-035
 Micro Boley
 Gauge Straight



RSD6-036
 Dentaurum Ruler




RSD6-037
 Ruler Wedge



RSD6-038
 Ruler



RSD6-039
 Ortho Ruler / Cone Socket



RSD6-040
 Dynamometer

TITANIUM PLASMA COLOUR COATED INSTRUMENTS Hollow Handle 10mm MIXING SPATULA WITH SHOVEL



RSD6-154
Bone Carrier
Titanium
18.5cm



RSD6-155
Raspator, Prichard
Length 18cm,
10 mm hollow handle,
double-ended
10.5 mm | 4.5 mm



RSD6-156
Raspator, Prichard
Length 18cm,
10 mm hollow handle,
double ended



RSD6-157
Raspator, Prichard
Length 18cm,
10 mm hollow handle,
double-ended
10.5 mm | 4.5 mm



RSD6-158
Mixing Spatula
with shovel
17cm




RSD6-159
Mixing Spatula
with shovel
18.5cm

BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS

Product offering includes instruments for harvesting, mixing and compacting bone Provides clinicians with the proper instrumentation for bone grafting procedures Weighted bone basin Prevents tipping and spilling of bone and other materials during mixing

Hollow Handle 9.5mm Length 170mm



RSD6-160
Graft Packer,
Concave Tips
DE 4/6mm

RSD6-161
Bone Spoon/
Packer DE 3mm

RSD6-162
Graft Packer,
Serrated Tips
DE 6/8mm

RSD6-163
Graft Packer,
Serrated Tips
DE 3/4mm

RSD6-164
Graft Packer,
Smooth Tips
DE 5/8mm

RSD6-165
Titanium graft packer/
placement site depth
gauge 3.3mm & 4.0mm
Diameter tips Marked
at 8, 11, 13, 15 & 18mm

RSD6-166
Graft Material Spoon
Concave Tip For
Larger Grafts

BONE GRAFTING INSTRUMENTS TITANIUM PLASMA COLOUR COATED Hollow Handle Ø 10 mm



RSD6-167
Blade and pusher
Length 18cm,
for straight or angled,
blade: 15 x 5 mm,
pusher: 4 mm diameter
Straight



RSD6-168
Blade and pusher
Length 18cm,
for straight or angled,
blade: 15 x 5 mm,
pusher: 4 mm diameter
angular



RSD6-169
Spoon and pusher
Length 18cm,
for angled spoon:
10 x 5 mm, pusher:
diameter 3 mm



RSD6-170
Tunnelierungs
Instrument
length 17.5cm,



RSD6-171
Tunnelierungs
Instrument
length 17.5cm,
different designs



RSD6-172
Tampers
Length 17cm,
for bone material,
double-ended,
diameter 3.2 mm
(2-4-6-8-10) and diameter
2.5mm (2-4-6-8-10)

Pluggers

Hollow Handle 9.5mm Length 160 mm, with titanium tips for mixing and pouring of autologous, heterologous bone substitute material, diameter (3/4 or 6/8 mm)

BONE GRAFTING INSTRUMENTS



RSD6-173 Ø 4mm
 RSD6-174 Ø 5mm
 Compactor
 15cm



RSD6-175
 Sinus Window
 Bone Tamper
 4 mm



RSD6-176
 Sinus Window
 Bone Tamper
 7 mm




RSD6-177
 6/8mm



RSD6-178
 3/4mm

SPOON AND PUSHER (TRINOVO)

Hollow Handle Ø 10 mm Length 180mm



RSD6-179
10x5mm
Spoon and
Pusher
with safety
stop

RSD6-180
10x5mm
Spoon and
Pusher
with safety
stop

RSD6-181
10x5mm
Spoon and
Pusher
with safety
stop



RSD6-182
Bone Graft/Carrier
Titanium
2-4-6-8-10mm

RSD6-183
15x10mm
Spoon and
Pusher
with safety
stop

RSD6-184
10x5mm
Spoon and
Pusher
with safety
stop



RSD6-185
10x5mm
Spoon and
Pusher
with safety
stop

OSTEOTOME CONVEXES / CONCAVE / CONIQUES



	CONVEXES STRAIGHT	RSD6-248	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES STRAIGHT	RSD6-254	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES STRAIGHT	RSD6-260	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm

OSTEOTOME CONVEXES / CONCAVE / CONIQUES



	CONVEXES CURVED	RSD6-266	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES CURVED	RSD6-272	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES CURVED	RSD6-278	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm

OSTEOTOME CONVEXES / CONCAVE / CONIQUES

RSD6-278-A

Osteotome Convexes Set of 10 Pcs



- 01 Osteotome Convexes (Str)
- 02 Osteotome Convexes (Str)
- 03 Osteotome Convexes (Str)
- 04 Osteotome Convexes (Str)
- 05 Osteotome Convexes (Str)

- 06 Osteotome Convexes (Cvd)
- 07 Osteotome Convexes (Cvd)
- 08 Osteotome Convexes (Cvd)
- 09 Osteotome Convexes (Cvd)
- 10 Osteotome Convexes (Cvd)

KEYS FOR IMPLANT SCREW



RSD6-313



RSD6-314



RSD6-315



RSD6-316

TORQUE WRENCH & RATCHET



RSD6-317 Star Torque Wrench 6.35mm Universal



RSD6-323 Star Torque Wrench 6.35mm Universal



RSD6-329 Star Ratchet Wrench 6.35mm Universal



RSD6-318 Round Torque Wrench 4.0mm Universal



RSD6-324 Star Torque Wrench 6.35mm Universal



RSD6-330 Star Ratchet Wrench 4.0mm Universal



RSD6-319 Hex Driver With 4.0mm Driver



RSD6-325 Star Torque Wrench 6.35mm Universal



RSD6-331 Round Ratchet Wrench 4.0mm Universal



RSD6-320 Star Torque Wrench 6.35mm Universal



RSD6-326 Star Torque Wrench 6.35mm Universal



RSD6-332 Round Ratchet Wrench 3.35mm Universal



RSD6-321 Hex Torque Wrench 4.00mm Universal



RSD6-327 Star Ratchet Wrench 6.35mm Universal



RSD6-333 Star Ratchet Wrench 4.0mm Universal



RSD6-322 Star Torque Wrench 6.35mm Universal



RSD6-328 Star Ratchet Wrench 6.35mm Universal



RSD6-334 Star Ratchet Wrench 6.35mm Universal

IMPLANTOLOGY

Инструменты для имплантологии

Implantologie

Implantologie


Инструменти для імплантології

Implantologia

زراعة الأسنان



RSD6-337
Bone Crusher



RSD6-335
Bone Crusher
Titanium Metal Cutter
Set of 2 Pcs
Bone Crusher + Mallet

mallet



RSD6-338
Bone Crusher
Bone Grinding Grain
1.5mm, 2mm, 2.5mm



RSD6-336
Bone Crusher
Titanium Metal Cutter
Set of 2 Pcs
Bone Crusher + Mallet

mallet



RSD6-339
Bone Crusher
Bone Grinding Grain
1.5mm, 2mm, 2.5mm

BONE EXPANDER SCREW & SURGICAL KITS



RSD6-340

Bone Universal Kit

Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm

Trephines Bit:

3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.

Driver:

2.5mm, 2.5mm Hand Pcs

1.25mm, 1.25mm Hand Pcs

Short Driver 2.5mm

Long Driver 2.5mm

1.25mm, Short Hex Driver Hand Pcs

1.25mm, Long Hex Driver Hand Pcs

Torque Wrench 4.0mm



RSD6-341

Bone Expander Kit

Expander: 2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm

Ratchet Driver

Parallel Pin

Pilot Drill

Drill 1.4mm

Drill Driver 4mm

Hand Driver 4mm

Ratchet 4.0mm



RSD6-342

Bone Universal Kit

Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm

Trephines Bit:

3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.

Bone Compression:

2mm, 3.3mm, 2.8mm, 4.2mm, 3mm, 3.5mm

Driver:

2.5mm, 2.5mm Hand Pcs

1.25mm, 1.25mm Hand Pcs

Short Driver 2.5mm

Long Driver 2.5mm

1.25mm, Short Hex Driver Hand Pcs

1.25mm, Long Hex Driver Hand Pcs

Torque Wrench 6.35mm 1 Pcs



RSD6-343

Conical Drills: 20 Pcs

Tissue Punches:

3.0mm, 3.50mm, 4.0mm, 4.5mm, 5.0mm, 6.0mm, 6.50mm

Bone Expander:

1.6/3.8mm, 1.9/3.0mm, 2.1/3.6mm,

2.6/3.8mm, 2.9/4.3mm

Drill Adapter: 1 Pcs

1.25mm, 2.50mm Hex Driver Hand Pcs

Pilot Drill: 1.6mm

Drill Bit: 1.4mm, 1.6mm, 1.8mm, 2.0mm, 2.2mm

Ratchet Driver: 4mm/6.35, 1.25/2.50

Rimer: 3.2/4.7mm, 3.7/4.7mm


Parallel Pin: 4.5mm, 2.5mm

Double step Drill: 2.4/2.8mm

Depth Gauge: 1 Pcs

Torque Wrench 6.35mm Universal: 1 Pcs

BONE INJECTORS, COLLECTORS & SUCTION TUBE



RSD6-344
Ø 2,5mm
Straight / Bend
18cm
I.D 2.5mm
O.D 3.5

RSD6-345
Ø 3,5mm
Straight / Bend
18cm
I.D 3.5mm
O.D 4.5

RSD6-346
Ø 4,5mm
Straight / Bend
18cm
I.D 4.5mm
O.D 5.5

RSD6-347
Ø 6,0mm
Straight / Bend
18cm
I.D 6.0mm
O.D 7.0

RSD6-348
Ø 2,5mm
17cm
I.D 2.5mm
O.D 3.5

RSD6-349
Ø 3,5mm
17cm
I.D 3.5mm
O.D 4.5

RSD6-350
Ø 4,5mm
17cm
I.D 4.5mm
O.D 5.5

RSD6-351
Ø 6,0mm
17cm
I.D 6.0mm
O.D 7.0



RSD6-352
Amalgam Gun



RSD6-353
17.5cm
Ø 20mm

With Filter



RSD6-354
18cm
Ø 27mm

With Filter



RSD6-355
coarse
Ø 5mm



RSD6-356
fine
Ø 9mm
Suction tubes
to collect
the extracted
bone material
180mm

Инструменты для имплантологии
Implantologie
Implantologie

Инструменти для імплантології
Implantologia
زراعة الأسنان

1.5mm



RSD6-357
Suction
Tubes
18cm

3mm



RSD6-358
Suction
Tubes
18cm

5mm



RSD6-359
Suction
Tubes
18cm



RSD6-361
Ø 3 mm
length 18cm
RSD6-362
Ø 4 mm
length 18cm
RSD6-363
Ø 5 mm
length 18cm

RSD6-360
Fig. 6FR
Frazier
15cm - Ø 2 mm



RSD6-364
Suction Cannula Kit
with 3 Needles



Suction Pipe



Suction Needle
1mm



Suction Needle
1.5mm



Suction Needle
2mm



IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع


PERFORATED RIM LOCK UPPER



LOWER




All the above impression trays are available in these 5 sizes



UPPER		U 5	U 4	U 3	U 2	U 1
		XS	S	M	L	XL
A	A	52	53	54	60	62
	B	68	69	73	75	85

LOWER		L 5	L 4	L 3	L 2	L 1
		XS	S	M	L	XL
A	A	52	53	54	60	62
	B	68	69	73	75	85



All Dimension are in mm

SOLID RIM LOCK UPPER




LOWER




IMPRESSION TRAYS

- Оттисные ложки Відбиткові ложки
 Porte-empreinte Cubetas para impresiones
 Abdrucklöffel صواني الانطباع


PERFORATED FULL DENTURE UPPER




LOWER



All the above impression trays are available in these 4 sizes



UPPER		S	M	L	XL
	A	60	64	67	70
UPPER		S	M	L	XL
	B	70	73	76	79




All Dimension are in mm

SOLID FULL DENTURE UPPER



LOWER



PERFORATED FULL DENTURE UPPER

SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
RSD7-037

SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-038

SOLID FULL DENTURE UPPER

SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L
RSD7-039

SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-040

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع

SOLID FULL DENTURE

UPPER




LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-049


SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-050

All the above impression trays are available in these 4 sizes



UPPER		S	M	L	XL
	A	60	64	67	70
UPPER	B	70	73	76	79

LOWER		S	S	M	L
	A	57	60	64	64
LOWER	B	70	70	73	86



All Dimension are in mm

PERFORATED FULL DENTURE

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-059

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-060

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع

UPPER




LOWER




SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-067

All the above impression trays are available in these 3 sizes



UPPER		U 3	U 2	U 1
		S	M	L
	A	57	64	67
	B	70	71	81

LOWER		L 3	L 2	L 1
		S	M	L
	A	57	57	60
B	76	79	81	



All Dimension are in mm

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-074

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع

English Pattern
Mesh Type

UPPER




LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-083


SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-084

All the above impression trays are available in these 4 sizes



UPPER		U 1	U 2	U 3	U 4
		S	M	L	XL
	A	12	13	15	15
	B	50	51	55	58
C	67	72	78	80	


LOWER		L 1	L 2	L3	L 4
	A	S	M	L	XL
		14	15	15	17
		B	50	51	55
C	70	77	79	80	



All Dimension are in mm

English Pattern
Punched Type

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-093


SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-094

IMPRESSION TRAYS

- | | |
|---|--|
|  Оттисные ложки |  Відбиткові ложки |
|  Porte-empreinte |  Cubetas para impresiones |
|  Abdrucklöffel |  صواني الانطباع |

ENGLISH PATTERN SOLID TYPE


UPPER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-103


SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-104

All the above impression trays are available in these 4 sizes



UPPER		U 1	U 2	U 3	U 4
		S	M	L	XL
	A	12	13	15	15
	B	50	51	55	58
C	67	72	78	80	

LOWER		L 1	L 2	L3	L 4
		S	M	L	XL
	A	14	15	15	17
B	50	51	55	60	
C	70	77	79	80	



All Dimension are in mm

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع

SOLID FULL DENTURE

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-113

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-114

PERFORATED FULL DENTURE

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-123

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-124

IMPRESSION TRAYS

- | | |
|------------------|--------------------------|
| Оттисковые ложки | Відбиткові ложки |
| Porte-empreinte | Cubetas para impresiones |
| Abdrucklöffel | صواني الانطباع |

ORTHODONTIAL SOLID

UPPER



S
RSD7-125



M
RSD7-126



L
RSD7-127



XL
RSD7-128

LOWER



S
RSD7-129



M
RSD7-130



L
RSD7-131



XL
RSD7-132

SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-133

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-134

ORTHODONTIAL PERFORATED

UPPER



S
RSD7-135



M
RSD7-136



L
RSD7-137



XL
RSD7-138

LOWER



S
RSD7-139



M
RSD7-140



L
RSD7-141



XL
RSD7-142

SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-143

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-144


IMPRESSION TRAYS

Оттисные ложки
 Porte-empreinte
 Abdrucklöffel


Відбиткові ложки
 Cubetas para impresiones
 صواني الانطباع

PERFORATED WITHOUT RIM LOCK

UPPER



LOWER



SET OF 06 PCS
 Upper: S, M, L
 Lower: S, M, L
RSD7-157


SET OF 08 PCS
 Upper: XS, S, M, L
 Lower: XS, S, M, L
RSD7-158

SET OF 08 PCS
 Upper: S, M, L, XL
 Lower: S, M, L, XL
RSD7-159

SET OF 10 PCS
 Upper: XS, S, M, L, XL
 Lower: XS, S, M, L, XL
RSD7-160


SET OF 12 PCS
 Upper: XXS, XS, S, M, L, XL
 Lower: XXS, XS, S, M, L, XL
RSD7-161

All the above impression trays are available in these 6 sizes



UPPER		U1	U2	U3	U4	U5	U6
		XXS	XS	S	M	L	XL
	A	61	64	66	69	71	73
	B	47	51	53	57	60	62

LOWER		L1	L2	L3	L4	L5	L6
		XXS	XS	S	M	L	XL
	A	55	60	72	74	76	78
	B	47	49	54	57	59	61



All Dimension are in mm

IMPRESSION TRAYS

Оттисные ложки Відбиткові ложки
 Porte-empreinte Cubetas para impresiones
 Abdrucklöffel صواني الانطباع

SOLID WITHOUT RIM LOCK

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-174


SET OF 08 PCS
Upper: XS, S, M, L
Lower: XS, S, M, L
RSD7-175

SET OF 08 PCS
Upper: S, M, L, XL
Lower: S, M, L, XL
RSD7-176

SET OF 10 PCS
Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL
RSD7-177


SET OF 12 PCS
Upper: XXS, XS, S, M, L, XL
Lower: XXS, XS, S, M, L, XL
RSD7-178

All the above impression trays are available in these 6 sizes



UPPER		U1	U2	U3	U4	U5	U6
		XXS	XS	S	M	L	XL
	A	61	64	66	69	71	73
	B	47	51	53	57	60	62

LOWER		L1	L2	L3	L4	L5	L6
		XXS	XS	S	M	L	XL
	A	55	60	72	74	76	78
	B	47	49	54	57	59	61



All Dimension are in mm

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع

PARTIAL SOLID



RSD7-179 44mm
RSD7-180 55mm
RSD7-181 60mm
RSD7-182 65mm



RSD7-183



RSD7-184 44mm
RSD7-185 55mm
RSD7-186 60mm
RSD7-187 65mm



RSD7-188



RSD7-189

Set of 5 Pcs RSD7-190

PARTIAL PERFORATED



RSD7-191 44mm
RSD7-192 55mm
RSD7-193 60mm
RSD7-194 65mm



RSD7-195



RSD7-196 44mm
RSD7-197 55mm
RSD7-198 60mm
RSD7-199 65mm



RSD7-200



RSD7-201



RSD7-202



RSD7-203

Set of 7 Pcs. RSD7-204

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع



Perforated



non-Perforated



1/2L

RSD7-205 1/2L/S
RSD7-206 1/2L/M
RSD7-207 1/2L/L

1/2R

RSD7-208 1/2R/S
RSD7-209 1/2R/M
RSD7-210 1/2R/L

1/2L

RSD7-211 1/2L/S
RSD7-212 1/2L/M
RSD7-213 1/2L/L

1/2R

RSD7-214 1/2R/S
RSD7-215 1/2R/M
RSD7-216 1/2R/L



Universal

RSD7-217 UN/1

Front

with swivel handle and fixing screw



Back



Universal

RSD7-218 UN/1

Front



RSD7-219 1

Front

with swivel handle



Back



RSD7-220 1

Front



Stolley

RSD7-221

ST/1

Partial occlusal stainless steel impression trays

They allow you to take an impression of two semi-arches antagonistic



RSD7-222

OCP



RSD7-223

OCG

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel


Відбиткові ложки
Cubetas para impresiones
صواني الانطباع



UPPER



	A mm.	B mm.
U 0	59	52
U 1	66	52
U 2	71	62
U 3	71	64
U 4	80	70



LOWER



	A mm.	B mm.
L 0	65	46
L 1	67	50
L 2	76	54
L 3	79	56
L 4	79	61



SET OF 10 PCS
Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL
RSD7-234

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع

Ehricke Stainless Steel Impression Trays

UPPER



LOWER



SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-241

	A mm.	B mm.	UPPER
U 1	64	46	
U 2	68	53	
U 3	79	61	

	A mm.	B mm.	LOWER
L 1	70	48	
L 2	70	53	
L 3	72	57	

IMPRESSION TRAYS

Оттисные ложки
Porte-empreinte
Abdrucklöffel

Відбиткові ложки
Cubetas para impresiones
صواني الانطباع

UPPER



S
RSD7-242
U1



M
RSD7-243
U2



L
RSD7-244
U3

LOWER



S
RSD7-245
L1



M
RSD7-246
L2




L
RSD7-247
L3


Quality Symbol

SET OF 06 PCS
Upper: S, M, L
Lower: S, M, L
RSD7-248

	A mm.	B mm.
U 1	62	50
U 2	70	56
U 3	78	62



	A mm.	B mm.
L 1	70	55
L 2	78	59
L 3	86	63



IMPRESSION TRAYS

- | | |
|---|--|
|  Оттисные ложки |  Відбиткові ложки |
|  Porte-empreinte |  Cubetas para impresiones |
|  Abdrucklöffel |  صواني الانطباع |



RSD7-249

Impression Trays Set of 6 Pcs

- 01 Impression Tray
- 02 Impression Tray
- 03 Impression Tray
- 04 Impression Tray
- 05 Impression Tray
- 06 Impression Tray



RSD7-250

Impression Tray Set of 6 Pcs

- 01 Impression Tray
- 02 Impression Tray
- 03 Impression Tray
- 04 Impression Tray
- 05 Impression Tray
- 06 Impression Tray



RSD7-251

Impression Tray Set of 6 Pcs

- 01 Impression Tray
- 02 Impression Tray
- 03 Impression Tray
- 04 Impression Tray
- 05 Impression Tray
- 06 Impression Tray





Scalers Tray
RSD8-001 200 x 100 mm



Scalers Trays
RSD8-002 205 x 115 x 17 mm
RSD8-003 205 x 152 x 17 mm
RSD8-004 290 x 115 x 17 mm



Instruments Trays
RSD8-005 200 x 110 x 30 mm
RSD8-006 215 x 115 x 35 mm
RSD8-007 220 x 130 x 30 mm
RSD8-008 230 x 125 x 35 mm
RSD8-009 245 x 140 x 40 mm
RSD8-010 285 x 160 x 45 mm
RSD8-011 315 x 175 x 50 mm
RSD8-012 330 x 185 x 50 mm



Instruments Trays
RSD8-013 150 x 100 x 50 mm
RSD8-014 200 x 100 x 50 mm
RSD8-015 220 x 120 x 50 mm
RSD8-016 240 x 140 x 50 mm
RSD8-017 200 x 150 x 50 mm
RSD8-018 250 x 200 x 50 mm
RSD8-019 300 x 250 x 50 mm
RSD8-020 300 x 200 x 50 mm
RSD8-021 350 x 250 x 50 mm
RSD8-022 400 x 300 x 50 mm
RSD8-023 450 x 300 x 50 mm
RSD8-024 450 x 350 x 50 mm



Instruments Trays With Lid
RSD8-025 150 x 100 x 50 mm
RSD8-026 200 x 100 x 50 mm
RSD8-027 220 x 120 x 50 mm
RSD8-028 240 x 140 x 50 mm
RSD8-029 200 x 150 x 50 mm
RSD8-030 250 x 200 x 50 mm
RSD8-031 300 x 250 x 50 mm
RSD8-032 300 x 200 x 50 mm
RSD8-033 350 x 250 x 50 mm
RSD8-034 400 x 300 x 50 mm
RSD8-035 450 x 300 x 50 mm
RSD8-036 450 x 350 x 50 mm



Mayo Tray Perforated
RSD8-037 150 x 100 x 25 mm
RSD8-038 200 x 100 x 25 mm
RSD8-039 220 x 120 x 25 mm
RSD8-040 240 x 140 x 25 mm
RSD8-041 200 x 150 x 25 mm
RSD8-042 250 x 200 x 25 mm
RSD8-043 300 x 250 x 25 mm
RSD8-044 300 x 200 x 25 mm
RSD8-045 350 x 250 x 25 mm
RSD8-046 400 x 300 x 25 mm
RSD8-047 450 x 300 x 25 mm
RSD8-048 450 x 350 x 25 mm



Scaler Tray
RSD8-049 200 x 115 mm
RSD8-050 205 x 150 mm
RSD8-051 290 x 150 mm
RSD8-052 300 x 200 mm



Instruments Tray Perforated
RSD8-053 150 x 100 x 50 mm
RSD8-054 200 x 100 x 50 mm
RSD8-055 220 x 120 x 50 mm
RSD8-056 240 x 140 x 50 mm
RSD8-057 200 x 150 x 50 mm
RSD8-058 250 x 200 x 50 mm
RSD8-059 300 x 250 x 50 mm
RSD8-060 300 x 200 x 50 mm
RSD8-061 350 x 250 x 50 mm
RSD8-062 400 x 300 x 50 mm
RSD8-063 450 x 300 x 50 mm
RSD8-064 450 x 350 x 50 mm

Sterilizing boxes & instruments Trays



Instruments Box with Hooks

RSD8-065	170 x 70 x 30 mm
RSD8-066	180 x 80 x 40 mm
RSD8-067	200 x 100 x 40 mm
RSD8-068	200 x 100 x 50 mm
RSD8-069	220 x 120 x 50 mm
RSD8-070	250 x 150 x 50 mm
RSD8-071	250 x 200 x 50 mm
RSD8-072	300 x 150 x 60 mm
RSD8-073	300 x 200 x 50 mm
RSD8-074	350 x 250 x 50 mm
RSD8-075	400 x 200 x 100 mm
RSD8-076	450 x 250 x 50 mm



Instruments Box

RSD8-077	170 x 70 x 30 mm
RSD8-078	180 x 80 x 40 mm
RSD8-079	200 x 100 x 40 mm
RSD8-080	200 x 100 x 50 mm
RSD8-081	220 x 120 x 50 mm
RSD8-082	250 x 150 x 50 mm
RSD8-083	250 x 200 x 50 mm
RSD8-084	300 x 150 x 60 mm
RSD8-085	300 x 200 x 50 mm
RSD8-086	350 x 250 x 50 mm
RSD8-087	400 x 200 x 100 mm
RSD8-088	450 x 250 x 50 mm



Instruments Box Perforated

RSD8-089	170 x 70 x 30 mm
RSD8-090	180 x 80 x 40 mm
RSD8-091	200 x 100 x 40 mm
RSD8-092	200 x 100 x 50 mm
RSD8-093	220 x 120 x 50 mm
RSD8-094	250 x 150 x 50 mm
RSD8-095	250 x 200 x 50 mm
RSD8-096	300 x 150 x 60 mm
RSD8-097	300 x 200 x 50 mm
RSD8-098	350 x 250 x 50 mm
RSD8-099	400 x 200 x 100 mm
RSD8-100	450 x 250 x 50 mm




Gloves Box

RSD8-101 300 x 150 mm



Catheters Tray with lid

RSD8-102	200 x 60 x 40 mm
RSD8-103	200 x 90 x 40 mm
RSD8-104	320 x 95 x 50 mm
RSD8-105	328 x 90 x 50 mm
RSD8-106	365 x 90 x 50 mm
RSD8-107	440 x 115 x 50 mm



Box With Filter

RSD8-108	100 x 50 mm
RSD8-109	250 x 150 mm
RSD8-110	150 x 65 mm
RSD8-111	200 x 100 mm
RSD8-112	300 x 200 mm



Instruments Trays With Cover

RSD8-113	210 x 110 x 50 mm
RSD8-114	240 x 160 x 50 mm
RSD8-115	300 x 180 x 50 mm
RSD8-116	360 x 210 x 50 mm
RSD8-117	400 x 230 x 55 mm

Instruments Trays Without Cover

RSD8-118	210 x 110 x 50 mm
RSD8-119	240 x 160 x 50 mm
RSD8-120	300 x 180 x 50 mm
RSD8-121	360 x 210 x 50 mm
RSD8-122	400 x 230 x 55 mm



Instruments Box with lid

RSD8-123	170 x 70 x 30 mm
RSD8-124	180 x 80 x 40 mm
RSD8-125	200 x 100 x 40 mm
RSD8-126	200 x 100 x 50 mm
RSD8-127	220 x 120 x 50 mm
RSD8-128	250 x 150 x 50 mm
RSD8-129	250 x 200 x 50 mm
RSD8-130	300 x 150 x 60 mm
RSD8-131	300 x 200 x 50 mm
RSD8-132	350 x 250 x 50 mm
RSD8-133	400 x 200 x 100 mm
RSD8-134	450 x 250 x 50 mm



Kidney Tray
RSD8-135 160 mm
RSD8-136 200 mm
RSD8-137 250 mm
RSD8-138 280 mm
RSD8-139 300 mm



Kidney Tray Shallow
RSD8-140 170 mm
RSD8-141 220 mm
RSD8-142 240 mm



Kidney Tray With Cover
RSD8-143 160 mm
RSD8-144 200 mm
RSD8-145 250 mm
RSD8-146 280 mm
RSD8-147 300 mm



Kidney Tray
RSD8-148 160 mm
RSD8-149 200 mm
RSD8-150 250 mm
RSD8-151 280 mm
RSD8-152 300 mm



Forceps Jars
RSD8-153 30 x 90 mm
RSD8-154 50 x 100 mm
RSD8-155 50 x 130 mm
RSD8-156 50 x 175 mm
RSD8-157 75 x 100 mm
RSD8-158 75 x 130 mm
RSD8-159 75 x 175 mm
RSD8-160 100 x 175 mm
RSD8-161 100 x 200 mm



Methot Jar Without Lid
RSD8-162



Dressing Jar / Cotton Container
RSD8-163 80 x 80 mm
RSD8-164 100 x 60 mm
RSD8-165 100 x 80 mm
RSD8-166 100 x 100 mm
RSD8-167 120 x 80 mm
RSD8-168 120 x 100 mm
RSD8-169 150 x 100 mm
RSD8-170 150 x 150 mm
RSD8-171 180 x 150 mm
RSD8-172 180 x 180 mm
RSD8-173 200 x 200 mm
RSD8-174 240 x 240 mm



Sputum Cup
RSD8-175



Dressing Jar With Lid & Knob
RSD8-176 80 x 80 mm
RSD8-177 90 x 90 mm
RSD8-178 100 x 100 mm
RSD8-179 120 x 120 mm
RSD8-180 100 x 130 mm
RSD8-181 150 x 150 mm
RSD8-182 125 x 170 mm



Round Bowls
RSD8-183 80 x 35 mm
RSD8-184 125 x 60 mm



Spirit Lamp
RSD8-185



Spirit Lamp
RSD8-186



Spirit Lamp
RSD8-187



Cotton Holding Jar
RSD8-188



Cotton Holder
RSD8-189



Cotton Holder
RSD8-190



**Cotton Wool Dispenser Acc.
to Methot**
RSD8-191

Quality Symbol



Needle Sterillization Tub
RSD8-192



Round Endodontic Boxes
complete with lid
RSD8-193 Ø mm. 94 - h. mm. 70



Endodontic Box
RSD8-194



RSD8-195
100 x 180 x 11mm
Mini size



RSD8-196
100 x 180 x 11mm
Mini size



RSD8-197
100 x 180 x 11mm
Mini size



RSD8-198
100 x 180 x 11mm
Mini size



RSD8-199
100 x 180 x 11mm
Mini size



RSD8-200
103 x 182 x 10mm
Mini size

Lid
Coperchio



RSD8-201
Silver



RSD8-202 Blue
RSD8-203 Red



RSD8-204 Blue
RSD8-205 Red



RSD8-206
Silver



RSD8-207
Green



RSD8-209
Blue



RSD8-208
Golden



RSD8-210
Red



RSD8-211
Silver



RSD8-212
Midi
125 x 105 x 54mm



RSD8-213
Maxi
204 x 105 x 54mm



RSD8-214
Maxi
204 x 105 x 54mm



Round Sterilizer Drum

RSD8-215	mm. Ø 80 x 80 H
RSD8-216	mm. Ø 125 x 80 H
RSD8-217	mm. Ø 125 x 100 H
RSD8-218	mm. Ø 125 x 125 H
RSD8-219	mm. Ø 150 x 80 H
RSD8-220	mm. Ø 150 x 100 H
RSD8-221	mm. Ø 150 x 130 H
RSD8-222	mm. Ø 150 x 150 H
RSD8-223	mm. Ø 150 x 190 H
RSD8-224	mm. Ø 190 x 125 H
RSD8-225	mm. Ø 190 x 145 H
RSD8-226	mm. Ø 190 x 165 H
RSD8-227	mm. Ø 190 x 190 H
RSD8-228	mm. Ø 190 x 240 H
RSD8-229	mm. Ø 190 x 290 H
RSD8-230	mm. Ø 240 x 120 H
RSD8-231	mm. Ø 240 x 145 H
RSD8-232	mm. Ø 240 x 165 H
RSD8-233	mm. Ø 240 x 190 H
RSD8-234	mm. Ø 240 x 240 H
RSD8-235	mm. Ø 240 x 290 H
RSD8-236	mm. Ø 240 x 340 H
RSD8-237	mm. Ø 240 x 390 H
RSD8-238	mm. Ø 290 x 145 H
RSD8-239	mm. Ø 290 x 165 H
RSD8-240	mm. Ø 290 x 190 H
RSD8-241	mm. Ø 290 x 240 H
RSD8-242	mm. Ø 290 x 290 H
RSD8-243	mm. Ø 290 x 340 H
RSD8-244	mm. Ø 290 x 390 H
RSD8-245	mm. Ø 340 x 145 H
RSD8-246	mm. Ø 340 x 165 H
RSD8-247	mm. Ø 340 x 190 H
RSD8-248	mm. Ø 340 x 240 H
RSD8-249	mm. Ø 340 x 290 H
RSD8-250	mm. Ø 340 x 340 H
RSD8-251	mm. Ø 340 x 390 H
RSD8-252	mm. Ø 390 x 140 H
RSD8-253	mm. Ø 390 x 165 H
RSD8-254	mm. Ø 390 x 190 H
RSD8-255	mm. Ø 390 x 240 H
RSD8-256	mm. Ø 390 x 290 H
RSD8-257	mm. Ø 390 x 340 H
RSD8-258	mm. Ø 390 x 390 H
RSD8-259	mm. Ø 390 x 490 H



Round Sterilizer Drums "Filter"

RSD8-260	mm. Ø 80 x 80 H
RSD8-261	mm. Ø 125 x 80 H
RSD8-262	mm. Ø 125 x 100 H
RSD8-263	mm. Ø 125 x 125 H
RSD8-264	mm. Ø 150 x 100 H
RSD8-265	mm. Ø 150 x 130 H
RSD8-266	mm. Ø 150 x 150 H
RSD8-267	mm. Ø 150 x 190 H
RSD8-268	mm. Ø 190 x 145 H
RSD8-269	mm. Ø 190 x 165 H
RSD8-270	mm. Ø 190 x 190 H
RSD8-271	mm. Ø 190 x 240 H
RSD8-272	mm. Ø 240 x 145 H
RSD8-273	mm. Ø 240 x 165 H
RSD8-274	mm. Ø 240 x 190 H
RSD8-275	mm. Ø 240 x 240 H
RSD8-276	mm. Ø 240 x 290 H
RSD8-277	mm. Ø 240 x 340 H
RSD8-278	mm. Ø 240 x 390 H
RSD8-279	mm. Ø 290 x 145 H
RSD8-280	mm. Ø 290 x 165 H
RSD8-281	mm. Ø 290 x 190 H
RSD8-282	mm. Ø 290 x 240 H
RSD8-283	mm. Ø 290 x 290 H
RSD8-284	mm. Ø 290 x 340 H
RSD8-285	mm. Ø 290 x 390 H
RSD8-286	mm. Ø 340 x 145 H
RSD8-287	mm. Ø 340 x 165 H
RSD8-288	mm. Ø 340 x 190 H
RSD8-289	mm. Ø 340 x 240 H
RSD8-290	mm. Ø 340 x 290 H
RSD8-291	mm. Ø 390 x 140 H
RSD8-292	mm. Ø 390 x 165 H
RSD8-293	mm. Ø 390 x 190 H
RSD8-294	mm. Ø 390 x 240 H
RSD8-295	mm. Ø 390 x 290 H
RSD8-296	mm. Ø 390 x 340 H
RSD8-297	mm. Ø 390 x 490 H



Round Bowls With Cover

RSD8-298	Ø 40 mm ml.30
RSD8-299	Ø 60 mm ml.50
RSD8-300	Ø 80 mm ml.150
RSD8-301	Ø 100 mm ml.250
RSD8-302	Ø 120 mm ml.450
RSD8-303	Ø 150 mm ml.750
RSD8-304	Ø 180 mm ml.1500
RSD8-305	Ø 200 mm ml.2100
RSD8-306	Ø 250 mm ml.3500

Round Bowls With Out Cover

RSD8-307	Ø 40 mm ml.30
RSD8-308	Ø 60 mm ml.50
RSD8-309	Ø 80 mm ml.150
RSD8-310	Ø 100 mm ml.250
RSD8-311	Ø 120 mm ml.450
RSD8-312	Ø 150 mm ml.750
RSD8-313	Ø 180 mm ml.1500
RSD8-314	Ø 200 mm ml.2100
RSD8-315	Ø 250 mm ml.3500

INSTRUMENTS MANAGEMENT SYSTEM



Cassette For 4 Instruments
RSD8-316



Cassette For 5 Instruments
RSD8-317



Cassette For 6 Instruments
RSD8-318



Cassette For 7 Instruments
RSD8-319



Cassette For 10 Instruments
RSD8-320



Cassette For 20 Instruments
RSD8-321



Cassette For Extraction Forceps
RSD8-322



Cassette For Root Elevator
RSD8-323



RSD8-324
Cassette Tray
15 Pcs



RSD8-325
Cassette Tray
15 Pcs

STERILIZING MESH TRAY

Different sizes & dimensions
are available at buyer's demand



MESH TRAYS

RSD8-326

Different sizes & dimensions
are available at buyer's demand



MESH TRAYS

RSD8-327

05 Pcs



MESH TRAYS

RSD8-328

07 Pcs



MESH TRAYS

RSD8-329

10 Pcs



MESH TRAYS

RSD8-330


20 Pcs



RSD8-331


Mesh Cassettes

5 Pieces, 7 Pieces, 10 Pieces,
15 Pieces & 20 Pieces



MESH TRAY

RSD8-332 240x130x50mm
RSD8-333 240x130x60mm
RSD8-334 240x250x50mm
RSD8-335 240x250x60mm
RSD8-336 480x250x50mm
RSD8-337 480x250x60mm
RSD8-338 480x250x70mm



MESH TRAY

RSD8-339 240x250x 50mm
RSD8-340 240x250x 60mm
RSD8-341 240x250x 70mm
RSD8-342 480x250x 50mm
RSD8-343 480x250x 60mm
RSD8-344 480x250x 70mm
RSD8-345 480x250x100mm



MESH TRAY

RSD8-346 240x250x50mm
RSD8-347 240x250x60mm
RSD8-348 240x250x70mm
RSD8-349 280x250x50mm
RSD8-350 280x250x60mm



MESH TRAY

RSD8-351 280x250x 70mm
RSD8-352 280x250x100mm
RSD8-353 280x250x 50mm
RSD8-354 280x280x 60mm
RSD8-355 280x280x 70mm
RSD8-356 300x280x 60mm



RSD8-357
18×7×3.5cm
5 Instruments Cassette
Detachable Lids
+
5 Instruments Cassette Hinged
Lids with button Lock



RSD8-358
29×18.5×3.5cm
20 Instruments Cassette Hinged
Lids with surgical kit fitting,
button lock




Instruments Cassette
RSD8-359 5 Pieces Cassette
RSD8-360 7 Pieces Cassette
RSD8-361 10 Pieces Cassette



RSD8-362
18×11×3.5cm
Extra large cassette for 6
Instruments, with compartment
for small parts hinged lids



RSD8-363
29×18.5×3.5cm
20 Instruments Cassette
Detachable Lids
+
20 Instruments Cassette Hinged
Lids with button Lock



RSD8-364
10×13×2.2cm
10 Instruments Cassette Hinged
Lids, Click Lock Available in 2
designs



RSD8-365
29×18.5×3.5cm
20 Instruments Cassette Detachable
Lids with elevator fitting
+
20 Instruments Cassette Hinged
Lids with elevator fitting , button lock



RSD8-366
18×9×2.2cm
7 Instruments Cassette Hinged
Lids, Click Lock, Available in 2
designs



RSD8-367
21×14×3.5cm
15 Instruments cassette,
Detachable lids, button lock



RSD8-368
21×14×3.5cm
15 Instruments cassette,
detachable lids, button
lock, with fitting of micro
surgery kit



RSD8-369
21×14×4.5cm
8 pcs luxator/elevator
cassette, Detach0
able lids, button lock



RSD8-370
15×9×5cm
Box for dental implant tools and
surgical kits, Capacity:8
tools, 18 tools, 24 Tools, 58 Tools
+ removable compartment



Instruments Tray Perforated
RSD8-371 288 x 187 x 40 mm



Instruments Tray Perforated Base
RSD8-372 288 x 187 x 40 mm



Instruments Tray Solid Base
RSD8-373 288 x 187 x 40 mm



Quality Symbol

RAIZ & SAJJAD SURGICAL PVT. LTD.

Manufactures and Exports of Fine Quality Surgical & Dental Instruments

Sadra Badra Industrial Zone Main Daska Road, Sialkot 51310 - Pakistan.

Tel: +92-52-3575217 - 3575218 Res: +92-52-4273118 Fax: +92-52-3575219

E-mail: msmughal@raizent.com, rsent@brain.net.pk

MQS
(cGMP)
Approved

